

SpectralRadar SDK

4.4

Generated by Doxygen 1.8.4

Mon Jun 26 2017 14:01:41

Contents

1	Spectral Radar SDK	1
1.1	SpectralRadar SDK License	1
1.2	Introduction	1
1.2.1	Overview	1
1.2.2	Data Handle (DataHandle, ColoredDataHandle, ComplexDataHandle, RawDataHandle)	1
1.2.3	OCTDeviceHandle	1
1.2.4	ProcessingHandle	2
1.2.5	ProbeHandle	2
1.2.6	ScanPatternHandle	2
1.2.7	Other Handles	2
1.3	First Steps	3
1.3.1	Initializing The Device	3
1.3.2	Creating Processing Routines	3
1.3.3	Creating A Scan Pattern	3
1.3.4	Acquisition	3
1.4	Error Handling	4
2	Module Index	4
2.1	Modules	4
3	Data Structure Index	5
3.1	Data Structures	5
4	File Index	5
4.1	File List	5
5	Module Documentation	5
5.1	Error Handling	5
5.1.1	Detailed Description	5
5.1.2	Enumeration Type Documentation	5
5.1.3	Function Documentation	6
5.2	Data Access	7
5.2.1	Detailed Description	10
5.2.2	Typedef Documentation	10
5.2.3	Enumeration Type Documentation	11
5.2.4	Function Documentation	12
5.3	Data Creation and Clearing	17
5.3.1	Detailed Description	17
5.3.2	Function Documentation	17
5.4	Hardware	19

5.4.1	Detailed Description	20
5.4.2	Typedef Documentation	20
5.4.3	Enumeration Type Documentation	20
5.4.4	Function Documentation	21
5.5	Internal Values	24
5.5.1	Detailed Description	24
5.5.2	Function Documentation	24
5.6	Pattern Factory/Probe	25
5.6.1	Detailed Description	27
5.6.2	Typedef Documentation	27
5.6.3	Enumeration Type Documentation	27
5.6.4	Function Documentation	28
5.7	Scan Pattern	31
5.7.1	Detailed Description	32
5.7.2	Typedef Documentation	32
5.7.3	Function Documentation	32
5.8	Acquisition	35
5.8.1	Detailed Description	35
5.8.2	Enumeration Type Documentation	35
5.8.3	Function Documentation	36
5.9	Processing	37
5.9.1	Detailed Description	40
5.9.2	Typedef Documentation	41
5.9.3	Enumeration Type Documentation	41
5.9.4	Function Documentation	43
5.10	Export and Import	48
5.10.1	Detailed Description	49
5.10.2	Enumeration Type Documentation	49
5.10.3	Function Documentation	51
5.11	Volume	53
5.11.1	Detailed Description	55
5.11.2	Enumeration Type Documentation	55
5.11.3	Function Documentation	55
5.12	ProbeCalibration	59
5.13	Doppler	60
5.13.1	Detailed Description	60
5.13.2	Typedef Documentation	61
5.13.3	Enumeration Type Documentation	61
5.13.4	Function Documentation	61
5.14	Service	63

5.14.1 Detailed Description	63
5.14.2 Function Documentation	63
5.15 Settings	64
5.15.1 Detailed Description	64
5.15.2 Typedef Documentation	64
5.15.3 Function Documentation	64
5.16 Coloring	66
5.16.1 Detailed Description	67
5.16.2 Typedef Documentation	67
5.16.3 Enumeration Type Documentation	67
5.16.4 Function Documentation	68
5.17 Camera	69
5.17.1 Detailed Description	69
5.17.2 Enumeration Type Documentation	69
5.17.3 Function Documentation	70
5.18 Buffer	71
5.18.1 Detailed Description	71
5.18.2 Typedef Documentation	71
5.18.3 Function Documentation	71
5.19 File Handling	73
5.19.1 Detailed Description	75
5.19.2 Enumeration Type Documentation	75
5.20 BETA_API	77
5.20.1 Detailed Description	77
5.20.2 Enumeration Type Documentation	77
5.21 Post Processing	78
5.21.1 Detailed Description	78
5.21.2 Enumeration Type Documentation	78
5.21.3 Function Documentation	78
5.22 Polarization	80
5.22.1 Detailed Description	80
5.22.2 Typedef Documentation	80
5.22.3 Enumeration Type Documentation	80
6 Data Structure Documentation	81
6.1 ComplexFloat Struct Reference	81
6.1.1 Detailed Description	81
6.1.2 Field Documentation	81
7 File Documentation	81
7.1 SpectralRadar.h File Reference	81

7.1.1	Detailed Description	110
7.1.2	Macro Definition Documentation	110
7.1.3	Typedef Documentation	110
7.1.4	Enumeration Type Documentation	110
7.1.5	Function Documentation	111

1 Spectral Radar SDK

1.1 SpectralRadar SDK License

By using the Thorlabs SpectralRadar SDK you agree to the terms and conditions detailed in the license agreement provided here: [THORLABS SpectralRadar SDK License Agreement](#) (PDF reader required). If this link does not work, you will also find this license agreement in Start Menu -> All Programs -> Thorlabs -> Spectral-Radar-SDK.

1.2 Introduction

This document gives an introduction into using the ANSI C Spectral Radar SDK and demonstrates the use of the most important functions.

1.2.1 Overview

The ANSI C Spectral Radar SDK follows an object-oriented approach. All objects are represented by pointers where appropriate typedefs are provided for convenience. The defined types are called Handles and are used as Return values when created and are passed as value when used. All functionality has been created with full LabVIEW compatibility in mind and it should be possible to use the SDK with most other programming languages as well. The most important handles are given in the following sections.

1.2.2 Data Handle (DataHandle, ColoredDataHandle, ComplexDataHandle, RawDataHandle)

Data acquired and used by the SDK is provided via data objects. A data object can contain

- floating point data (via [DataHandle](#))
- complex floating point data (via [ComplexDataHandle](#))
- ARGB32 colored data (via [ColoredDataHandle](#))
- unprocessed RAW data (via [RawDataHandle](#)) The data objects store all information belonging to them, such as pixel data, spacing between pixels, comments attached to their data, etc. Data objects are automatically resized if necessary and can contain 1-, 2- or 3-dimensional data. The dimensionality can be read by [getDataPropertyInt\(\)](#), etc. Direct access to their memory is possible via [getDataPtr\(\)](#), etc. Data properties can be read out via [getDataPropertyInt\(\)](#), [getDataPropertyFloat\(\)](#), etc. These include sizes along their first, second and third axis, physical spacing between pixels, their total range, etc.

1.2.3 OCTDeviceHandle

A handle specifying the OCT device that is used. In most cases the [OCTDeviceHandle](#) is obtained using the [initDevice\(\)](#) function and needs to be closed after using by [closeDevice\(\)](#). The complete device will be initialized, the SLD will be switched on and all start-up dependent calibration will be performed. All hardware and hardware dependend actions require the [OCTDeviceHandle](#) to be passed. These include for example

- starting and stopping a measurement ([startMeasurement\(\)](#) and [stopMeasurement\(\)](#))
- getting properties of the device ([getDevicePropertyInt\(\)](#) and [getDevicePropertyFloat\(\)](#))

1.2.4 ProcessingHandle

The numerics and processing routines required in order to create A-scans, B-scans and volumes out of directly measured spectra can be accessed via the [ProcessingHandle](#). When the [ProcessingHandle](#) is created, all required temporary memory and routines are initialized and prepared and several threads are started. In most cases the ideal way to create a processing handle is to use [createProcessingForDevice\(\)](#) which creates optimized processing algorithms for the [OCTDeviceHandle](#) specified. If no device is available or the processing routines are to be tweaked manually [createProcessing\(\)](#) must be used. When all required processing is done, [closeProcessing](#) must be used to stop all processing threads and free all temporary memory. All functions whose output is dependent on the processing routines used have a [ProcessingHandle](#) parameter. These include for example

- The [setProcessingParameterInt\(\)](#) and [setProcessingFlag\(\)](#) functions for setting parameters that are used for processing
- The [executeProcessing\(\)](#) function for triggering the processing of raw data

1.2.5 ProbeHandle

The probe is the hardware used for scanning the sample, usually with help of galvanometric scanners. The object referenced by [ProbeHandle](#) is responsible for creating scan patterns and holds all information and settings of the probe attached to the device. It needs to be calibrated to map suitable output voltage (for analog galvo drivers) or digital values (for digital galvo drivers) to scanning angles, inches or millimeters. In most cases this calibration data is provided by *.ini files and the probe is initialized by [initProbe\(\)](#) where the probe configuration file name needs to be specified as a string parameter. Probes calibrated at Thorlabs will usually come with a factory-made probe configuration file which follows the nomenclature Probe + Objective Name.ini, e.g. "ProbeLSM03.ini"

If the probe is to be hardcoded into the software one can also provide an empty string as parameter and provide the configuration manually using the [setProbeParameterInt\(\)](#) and [setProbeParameterFloat\(\)](#) functions. When the Probe object is no longer needed, [closeProbe\(\)](#) must be called to free temporary memory. All actions that depend on the probe configuration require a [ProbeHandle](#) to be specified, such as:

- move galvo scanner to a specific position ([moveScanner\(\)](#)).
- create a scan pattern ([createBScanPattern\(\)](#)), see also [ScanPatternHandle](#).
- set calibration parameters for a specific probe ([setProbeParameterFloat\(\)](#) and [setProbeParameterInt\(\)](#))

1.2.6 ScanPatternHandle

A scan pattern is used to specify the points on the probe to scan during data acquisition, and its information is accessible via the [ScanPatternHandle](#). A dedicated function can be used to create a specific scan pattern, such as [createBScanPattern\(\)](#) for a simple B-scan or [createVolumePattern\(\)](#) for a simple volume scan. When the scan pattern is no longer needed its resources can be freed using [clearScanPattern\(\)](#). The [ScanPatternHandle](#) needs to be specified to all functions that need information on the resulting scan. For example:

- creating a pattern ([createBScanPattern\(\)](#), [createVolumePattern\(\)](#), etc.)
- starting a measurement ([startMeasurement\(\)](#))

1.2.7 Other Handles

Other Handles that are used in the Spectral Radar SDK are

- [DopplerProcessingHandle](#): Handle to Doppler processing routines that can be used to transform complex data to Doppler phase and amplitude signals.
- [SettingsHandle](#): Handle to an INI file that can be read and written to without explicitly taking care of parsing the file.
- [Coloring32BitHandle](#): Handle to processing routines that can map floating point data to 32 bit color data.

1.3 First Steps

The following section describes first steps that are needed to acquire data with the Spectral Radar SDK.

1.3.1 Initializing The Device

The easiest way to initialize the device is to use the [initDevice\(\)](#) function. It returns an appropriate [OCTDeviceHandle](#) that can be used to identify the device:

```
OCTDeviceHandle Dev = initDevice();
// Acquire data, processing, direct hardware access...
closeDevice(Dev);
```

1.3.2 Creating Processing Routines

In most cases raw data acquired by the OCT device needs to be transformed using a Fast Fourier transform and other pre- and postprocessing algorithms. To get a [ProcessingHandle](#) on these algorithms the most convenient way is to use the [createProcessingForDevice\(\)](#) functionality which requires a valid [OCTDeviceHandle](#):

```
// ...
ProcessingHandle Proc = createProcessingForDevice(Dev);
// acquire data and perform processing
closeProcessing(Proc);
// ...
```

1.3.3 Creating A Scan Pattern

In order to scan a sample and acquire B-scan OCT data one needs to specify a scan pattern that describes at which point to acquire data. To get the data of a simple B-Scan one can simply use [createBScanPattern\(\)](#):

```
// ...
ProbeHandle Probe = initProbe(Dev, "Probe");
ScanPatternHandle Pattern = createBScanPattern(Probe, 2.0, 512); // get
                        B-scans with 2.0mm scanning angle and 512 A-scans per B-scan
// acquire data, ...
clearScanPattern(Pattern);
closeProbe(Probe);
// ...
```

1.3.4 Acquisition

The most convenient and fast way to acquire data is to acquire data asynchronously. For this one starts a measurement using [startMeasurement\(\)](#) and retrieves the latest available [getRawData\(\)](#). The memory needed to store the data needs to be allocated first:

```
int i;

RawDataHandle Raw = createRawData();
DataHandle BScan = createRealData();
startMeasurement(Dev, Pattern, Acquisition_ASyncContinuous);

for(i=0; i<1000; ++i) // get 1000 B-scans
{
    getRawData(Raw);
}
```

```

    setProcessedDataOutput(Proc, BScan);
    executeProcessing(Proc, Raw);
    // data is now in BScan...
    // do something with the data...
}

stopMeasurement(Dev);
clearData(BScan);
clearRawData(Raw);

```

1.4 Error Handling

Error handling is done by calling the function `getError()`. The function will return an `ErrorCode` and if the result is not `NO_ERROR` an error string will be provided giving details about the problem.

```

#define ERROR_STRLen 1024;
//...
char error[ERROR_STRLen];
OCTDeviceHandle Dev = initDevice();
if(!getError(error, ERROR_STRLen)) // check whether the previous calls to SDK functions caused an
    error
{
    printf("An error occurred: %s", error);
}
// ...

```

2 Module Index

2.1 Modules

Here is a list of all modules:

Error Handling	5
Data Access	7
Data Creation and Clearing	17
Hardware	19
Internal Values	24
Pattern Factory/Probe	25
Scan Pattern	31
Acquisition	35
Processing	37
Export and Import	48
Volume	53
ProbeCalibration	59
Doppler	60
Service	63
Settings	64
Coloring	66

Camera	69
Buffer	71
File Handling	73
BETA_API	77
Post Processing	78
Polarization	80

3 Data Structure Index

3.1 Data Structures

Here are the data structures with brief descriptions:

ComplexFloat	
A standard complex data type that is used to access complex data	81

4 File Index

4.1 File List

Here is a list of all documented files with brief descriptions:

SpectralRadar.h	
Header containing all functions of the Spectral Radar SDK. This SDK can be used for Callisto, Ganymede, Hyperion and Telesto devices	81

5 Module Documentation

5.1 Error Handling

Error handling.

Enumerations

- enum **ErrorCode** {
 NoError = 0x00000000,
 Error = 0xE0000000 }
This enum is used to describe errors that occur when operating an OCT device.

Functions

- SPECTRALRADAR_API ErrorCode getError** (char *Message, int StringSize)
Returns an error code and a message if an error occurred. The error flag will be cleared.

5.1.1 Detailed Description

Error handling.

5.1.2 Enumeration Type Documentation

5.1.2.1 enum **ErrorCode**

This enum is used to describe errors that occur when operating an OCT device.

Warning

Error codes and error description texts are subject to change in future releases.

Enumerator

NoError No error occurred. This entry can be cast to FALSE.

Error Error occurred. This entry can be cast to TRUE.

5.1.3 Function Documentation

5.1.3.1 **ErrorCode** `getError (char * Message, int StringSize)`

Returns an error code and a message if an error occurred. The error flag will be cleared.

Parameters

<i>Message</i>	Error message describing the error.
<i>StringSize</i>	Size of the string that was given to <i>Message</i> .

5.2 Data Access

Functions for accessing the information stored in data objects.

Data Structures

- struct [ComplexFloat](#)
A standard complex data type that is used to access complex data.

Typedefs

- typedef struct C_RawData * [RawDataHandle](#)
Handle to an object holding the unprocessed raw data.
- typedef struct C_Data * [DataHandle](#)
Handle to an object holding 1-, 2- or 3-dimensional floating point data.
- typedef struct C_ColoredData * [ColoredDataHandle](#)
Handle to an object holding 1-, 2- or 3-dimensional colored data.
- typedef struct C_ComplexData * [ComplexDataHandle](#)
Handle to an object holding complex 1-, 2- or 3-dimensional complex floating point data.
- typedef struct
C_ImageFieldCorrection * [ImageFieldHandle](#)
Handle to the image field description.
- typedef struct C_FileHandling * [OCTFileHandle](#)
Handle to the OCT file class.

Enumerations

- enum [DataPropertyInt](#) {
[Data_Dimensions](#),
[Data_Size1](#),
[Data_Size2](#),
[Data_Size3](#),
[Data_NumberOfElements](#),
[Data_SizeInBytes](#),
[Data_BytesPerElement](#) }
Selects integer point data property.
- enum [RawDataPropertyInt](#) {
[RawData_Size1](#),
[RawData_Size2](#),
[RawData_Size3](#),
[RawData_NumberOfElements](#),
[RawData_SizeInBytes](#),
[RawData_BytesPerElement](#),
[RawData_LostFrames](#) }
Specifies properties of RawData.
- enum [DataPropertyFloat](#) {
[Data_Spacing1](#),
[Data_Spacing2](#),
[Data_Spacing3](#),
[Data_Range1](#),
[Data_Range2](#),
[Data_Range3](#) }
Selects floating point data property.

- enum [DataAnalyzation](#) {
[Data_Min](#),
[Data_Mean](#),
[Data_Max](#),
[Data_MaxDepth](#) }
Selects data property to analyze.
- enum [AScanAnalyzation](#) {
[Data_Noise_dB](#),
[Data_Noise_electrons](#),
[Data_PeakPos_Pixel](#),
[Data_PeakPos_PhysUnits](#),
[Data_PeakHeight_dB](#),
[Data_PeakWidth_6dB](#),
[Data_PeakWidth_20dB](#),
[Data_PeakWidth_40dB](#) }
Selects an appropriate A-Scan analyzation.
- enum [DataOrientation](#) {
[DataOrientation_ZXY](#),
[DataOrientation_ZYX](#),
[DataOrientation_XZY](#),
[DataOrientation_XYZ](#),
[DataOrientation_YXZ](#),
[DataOrientation_YZX](#),
[DataOrientation_ZTX](#),
[DataOrientation_ZXT](#) }
Selects the orientation of the data.

Functions

- [SPECTRALRADAR_API](#) int [getDataPropertyInt](#) ([DataHandle](#) Data, [DataPropertyInt](#) Selection)
Returns the selected integer property of the specified data.
- [SPECTRALRADAR_API](#) double [getDataPropertyFloat](#) ([DataHandle](#) Data, [DataPropertyFloat](#) Selection)
Returns the selected floating point property of the specified data.
- [SPECTRALRADAR_API](#) void [copyData](#) ([DataHandle](#) DataSource, [DataHandle](#) DataDestination)
Copies the content of the specified source to the specified destination.
- [SPECTRALRADAR_API](#) void [copyDataContent](#) ([DataHandle](#) DataSource, float *Destination)
Copies the data in the specified data object ([DataHandle](#)) into the specified pointer.
- [SPECTRALRADAR_API](#) float * [getDataPtr](#) ([DataHandle](#) Data)
Returns a pointer to the content of the specified data.
- [SPECTRALRADAR_API](#) void [reserveData](#) ([DataHandle](#) Data, int Size1, int Size2, int Size3)
Reserves the amount of data specified. This might improve performance if appending data to the [DataHandle](#) as no additional memory needs to be reserved then.
- [SPECTRALRADAR_API](#) void [resizeData](#) ([DataHandle](#) Data, int Size1, int Size2, int Size3)
Resizes the respective data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.
- [SPECTRALRADAR_API](#) void [setDataRange](#) ([DataHandle](#) Data, double range1, double range2, double range3)
Sets the range in mm in the 3 axes represented in the [RealData](#) buffer.
- [SPECTRALRADAR_API](#) void [setDataContent](#) ([DataHandle](#) Data, float *NewContent)
*Sets the data content of the data object. The data chunk pointed to by NewContent needs to be of the size expected by the data object, i. e. Size1*Size2*Size3*sizeof(float).*
- [SPECTRALRADAR_API](#) [DataOrientation](#) [getDataOrientation](#) ([DataHandle](#) Data)
Returns the data orientation of the data object.
- [SPECTRALRADAR_API](#) void [setDataOrientation](#) ([DataHandle](#) Data, [DataOrientation](#) Orientation)

- Sets the data orientation of the data object to the given orientation.*
- **SPECTRALRADAR_API** int **getComplexDataPropertyInt** (**ComplexDataHandle** Data, **DataPropertyInt** Selection)
Returns the selected integer property of the specified data.
 - **SPECTRALRADAR_API** void **copyComplexDataContent** (**ComplexDataHandle** DataSource, **ComplexFloat** *Destination)
Copies the content of the complex data to the pointer specified as destination.
 - **SPECTRALRADAR_API** **ComplexFloat** * **getComplexDataPtr** (**ComplexDataHandle** Data)
*Returns a pointer to the data represented by the **ComplexDataHandle**. The data is still managed by the **ComplexDataHandle** object.*
 - **SPECTRALRADAR_API** void **setComplexDataContent** (**ComplexDataHandle** Data, **ComplexFloat** *NewContent)
*Sets the data content of the **ComplexDataHandle** to the content specified by the pointer.*
 - **SPECTRALRADAR_API** void **reserveComplexData** (**ComplexDataHandle** Data, int Size1, int Size2, int Size3)
*Reserves the amount of data specified. This might improve performance if appending data to the **ComplexDataHandle** as no additional memory needs to be reserved then.*
 - **SPECTRALRADAR_API** void **resizeComplexData** (**ComplexDataHandle** Data, int Size1, int Size2, int Size3)
Resizes the respective data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.
 - **SPECTRALRADAR_API** void **setComplexDataRange** (**ComplexDataHandle** Data, double range1, double range2, double range3)
*Sets the range in mm in the 3 axes represented in the **RealData** buffer.*
 - **SPECTRALRADAR_API** int **getColoredDataPropertyInt** (**ColoredDataHandle** ColData, **DataPropertyInt** Selection)
Returns the selected integer property of the specified colored data.
 - **SPECTRALRADAR_API** double **getColoredDataPropertyFloat** (**ColoredDataHandle** ColData, **DataPropertyFloat** Selection)
Returns the selected integer property of the specified colored data.
 - **SPECTRALRADAR_API** void **copyColoredData** (**ColoredDataHandle** ImageSource, **ColoredDataHandle** ImageDestination)
*Copies the contents of the specified **ColoredDataHandle** to the specified destination **ColoredDataHandle**.*
 - **SPECTRALRADAR_API** void **copyColoredDataContent** (**ColoredDataHandle** Source, unsigned long *Destination)
*Copies the data in the specified colored data object (**ColoredDataHandle**) into the specified pointer.*
 - **SPECTRALRADAR_API** void **copyColoredDataContentAligned** (**ColoredDataHandle** ImageSource, unsigned long *Destination, int Alignment1)
*Copies the data in the specified colored data object (**ColoredDataHandle**) into the specified pointer. This function assures the data to be aligned accordingly.*
 - **SPECTRALRADAR_API** unsigned long * **getColoredDataPtr** (**ColoredDataHandle** ColData)
*Returns a pointer to the content of the specified **ColoredDataHandle**.*
 - **SPECTRALRADAR_API** void **resizeColoredData** (**ColoredDataHandle** ColData, int Size1, int Size2, int Size3)
Resizes the respective colored data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.
 - **SPECTRALRADAR_API** void **reserveColoredData** (**ColoredDataHandle** ColData, int Size1, int Size2, int Size3)
*Reserves the amount of colored data specified. This might improve performance if appending data to the **ColoredDataHandle** as no additional memory needs to be reserved then.*
 - **SPECTRALRADAR_API** void **setColoredDataContent** (**ColoredDataHandle** ColData, unsigned long *NewContent)
*Sets the data content of the colored data object. The data chunk pointed to by NewContent needs to be of the size expected by the data object, i. e. Size1*Size2*Size3*sizeof(unsigned long).*
 - **SPECTRALRADAR_API** void **setColoredDataRange** (**ColoredDataHandle** Data, double range1, double range2, double range3)
Sets the range in mm in the 3 axes represented in the data object buffer.

- **SPECTRALRADAR_API** **DataOrientation** **getColorDataOrientation** (**ColoredDataHandle** Data)
Returns the data orientation of the colored data object.
- **SPECTRALRADAR_API** void **setColoredDataOrientation** (**ColoredDataHandle** Data, **DataOrientation**)
Sets the data orientation of the colored data object to the given orientation.
- **SPECTRALRADAR_API** void **getRawDataSize** (**RawDataHandle** Raw, int *SizeX, int *SizeY, int *SizeZ)
*Returns the size of the specified raw data (**RawDataHandle**).*
- **SPECTRALRADAR_API** void **copyRawDataContent** (**RawDataHandle** RawDataSource, void *DataContent)
Copies the content of the raw data into the specified buffer. The user needs to assure that enough memory is allocated.
- **SPECTRALRADAR_API** void * **getRawDataPtr** (**RawDataHandle** RawDataSource)
*Returns the pointer to the raw data content. The pointer might no longer be valid after additional actions using the **RawDataHandle**.*
- **SPECTRALRADAR_API** int **getRawDataPropertyInt** (**RawDataHandle** RawData, **RawDataPropertyInt** Property)
Returns a raw data property.
- **SPECTRALRADAR_API** void **setRawDataBytesPerPixel** (**RawDataHandle** Raw, int BytesPerPixel)
Sets the bytes per pixel for raw data.
- **SPECTRALRADAR_API** void **resizeRawData** (**RawDataHandle** Raw, int Size1, int Size2, int Size3)
Resizes the specified raw data buffer accordingly.
- **SPECTRALRADAR_API** void **setRawDataContent** (**RawDataHandle** RawDataSource, void *NewContent)
*Sets the content of the raw data buffer. The size of the **RawDataHandle** needs to be adjusted first, as otherwise not all data might be copied.*
- **SPECTRALRADAR_API** void **setScanSpectra** (**RawDataHandle** RawData, int NumberOfScanRegions, int *ScanRegions)
Sets the number of the spectra in the raw data that are used for creating A-scan/B-scan data.
- **SPECTRALRADAR_API** void **setApodizationSpectra** (**RawDataHandle** RawData, int NumberOfScanRegions, int *ApodizationRegions)
Sets the number of the spectra in the raw data that contain data useful as apodization spectra.
- **SPECTRALRADAR_API** int **getNumberOfScanRegions** (**RawDataHandle** Raw)
Returns the number of regions that have been acquired that contain scan data, i. e. spectra that are used to compute A-scans.
- **SPECTRALRADAR_API** int **getNumberOfApodizationRegions** (**RawDataHandle** Raw)
Returns the number of regions in the raw data containing spectra that are supposed to be used for apodization.
- **SPECTRALRADAR_API** void **getScanSpectra** (**RawDataHandle** Raw, int *SpectralIndex)
Returns the indices of spectra that contain scan data, i. e. spectra that are supposed to be used to compute A-scans.
- **SPECTRALRADAR_API** void **getApodizationSpectra** (**RawDataHandle** Raw, int *SpectralIndex)
Returns the indices of spectra that contain apodization data, i. e. spectra that are supposed to be used as input for apodization.

5.2.1 Detailed Description

Functions for accessing the information stored in data objects.

5.2.2 Typedef Documentation

5.2.2.1 ColoredDataHandle

Handle to an object holding 1-, 2- or 3-dimensional colored data.

5.2.2.2 ComplexDataHandle

Handle to an object holding complex 1-, 2- or 3-dimensional complex floating point data.

5.2.2.3 DataHandle

Handle to an object holding 1-, 2- or 3-dimensional floating point data.

5.2.2.4 ImageFieldHandle

Handle to the image field description.

5.2.2.5 OCTFileHandle

Handle to the OCT file class.

5.2.2.6 RawDataHandle

Handle to an object holding the unprocessed raw data.

5.2.3 Enumeration Type Documentation

5.2.3.1 enum AScanAnalyzation

Selects an appropriate A-Scan analyzation.

Enumerator

Data_Noise_dB Noise of the A-scan in dB. This assumes that no signal is present in the A-scan. The noise is computed by averaging all fourier channels larger than 50.

Data_Noise_electrons Noise of the A-scan in electrons. This assumes that no signal is present in the A-scan. The noise is computed by averaging all fourier channels larger than 50.

Data_PeakPos_Pixel Peak position of the highest peak in pixels. The peak position is determined by computing a parable going through the maximum value point and its surrounding pixels. The position of the maximum is used.

Data_PeakPos_PhysUnits Peak position of the highest peak in physical units. The peak position is determined by computing a parable going through the maximum value point and its surrounding pixels. The position of the maximum is used. Physical coordinates are computed by using the calibrated zSpacing property of the device. The concrete physical units of the return value depends on the calibration.

Data_PeakHeight_dB Peak height of the highest peak in dB. The peak hieght is determined by computing a parable going through the maximum value point and its surrounding pixels. The height of the resulting parable is returned.

Data_PeakWidth_6dB Signal width at -6dB. This is the FWHM.

Data_PeakWidth_20dB Signal width at -20dB.

Data_PeakWidth_40dB Signal width at -40dB.

5.2.3.2 enum DataAnalyzation

Selects data property to analyze.

Enumerator

Data_Min Minimum of the values in the data.

Data_Mean Arithmetic mean of all values in the data.

Data_Max Maximum of the values in the data.

Data_MaxDepth The depth of the maximum of the values in the data.

5.2.3.3 enum DataOrientation

Selects the orientation of the data.

5.2.3.4 enum DataPropertyFloat

Selects floating point data property.

Enumerator

Data_Spacing1 Spacing between two subsequent data elements in direction of the first axis in physical units.

Data_Spacing2 Spacing between two subsequent data elements in direction of the second axis in physical units.

Data_Spacing3 Spacing between two subsequent data elements in direction of the third axis in physical units.

Data_Range1 Total range of the data in direction of the first axis in physical units.

Data_Range2 Total range of the data in direction of the second axis in physical units.

Data_Range3 Total range of the data in direction of the third axis in physical units.

5.2.3.5 enum DataPropertyInt

Selects integer point data property.

Enumerator

Data_Dimensions Dimension of the data object. Usually 1, 2 or 3. 0 indicates empty data.

Data_Size1 Size of the first dimension. For OCT data this is usually the longitudinal axis (z)

Data_Size2 Size of the first dimension. For OCT data this is usually a transversal axis (x)

Data_Size3 Size of the first dimension. For OCT data this is usually a transversal axis (y)

Data_NumberOfElements The number of elements in the data object.

Data_SizeInBytes The size of the data object in bytes.

Data_BytesPerElement The number of bytes of a single element.

5.2.3.6 enum RawDataPropertyInt

Specifies properties of RawData.

Enumerator

RawData_Size1 Size of the first dimension. This will be the spectral dimension, i. e. z-dimension prior to Fourier transformation.

RawData_Size2 Size of the second dimension. This is a transversal axis (x).

RawData_Size3 Size of the third dimension. This is a transversal axis (y).

RawData_NumberOfElements The number of elements in the raw data object.

RawData_SizeInBytes The size of the data object in bytes.

RawData_BytesPerElement The number of bytes of a single element, i. e. the data type of the raw data.

RawData_LostFrames The number of lost frames during data acquisition.

5.2.4 Function Documentation

5.2.4.1 void copyColoredData (ColoredDataHandle ImageSource, ColoredDataHandle ImageDestination)

Copies the contents of the specified [ColoredDataHandle](#) to the specified destination [ColoredDataHandle](#).

5.2.4.2 void copyColoredDataContent (ColoredDataHandle Source, unsigned long * Destination)

Copies the data in the specified colored data object ([ColoredDataHandle](#)) into the specified pointer.

5.2.4.3 void copyColoredDataContentAligned (ColoredDataHandle ImageSource, unsigned long * Destination, int Alignment1)

Copies the data in the specified colored data object (ColoredDataHandle) into the specified pointer. This function assures the data to be aligned accordingly.

5.2.4.4 void copyComplexDataContent (ComplexDataHandle DataSource, ComplexFloat * Destination)

Copies the content of the complex data to the pointer specified as destination.

5.2.4.5 void copyData (DataHandle DataSource, DataHandle DataDestination)

Copies the content of the specified source to the specified destination.

5.2.4.6 void copyDataContent (DataHandle DataSource, float * Destination)

Copies the data in the specified data object (DataHandle) into the specified pointer.

5.2.4.7 void copyRawDataContent (RawDataHandle RawDataSource, void * DataContent)

Copies the content of the raw data into the specified buffer. The user needs to assure that enough memory is allocated.

5.2.4.8 void getApodizationSpectra (RawDataHandle Raw, int * SpectralIndex)

Returns the indices of spectra that contain apodization data, i. e. spectra that are supposed to be used as input for apodization.

An array needs to be provided that has twice the size of the number of apodization regions which can be obtained by [getNumberOfApodizationRegions\(\)](#)

5.2.4.9 DataOrientation getColoredDataOrientation (ColoredDataHandle Data)

Returns the data orientation of the colored data object.

5.2.4.10 int getColoredDataPropertyFloat (ColoredDataHandle ColData, DataPropertyFloat Selection)

Returns the selected integer property of the specified colored data.

5.2.4.11 int getColoredDataPropertyInt (ColoredDataHandle ColData, DataPropertyInt Selection)

Returns the selected integer property of the specified colored data.

5.2.4.12 unsigned long * getColoredDataPtr (ColoredDataHandle ColData)

Returns a pointer to the content of the specified ColoredDataHandle.

5.2.4.13 int getComplexDataPropertyInt (ComplexDataHandle Data, DataPropertyInt Selection)

Returns the selected integer property of the specified data.

5.2.4.14 ComplexFloat * getComplexDataPtr (ComplexDataHandle Data)

Returns a pointer to the data represented by the ComplexDataHandle. The data is still managed by the ComplexDataHandle object.

5.2.4.15 void DataOrientation getDataOrientation (DataHandle Data)

Returns the data orientation of the data object.

5.2.4.16 double getDataPropertyFloat (*DataHandle Data*, *DataPropertyFloat Selection*)

Returns the selected floating point property of the specified data.

5.2.4.17 int getDataPropertyInt (*DataHandle Data*, *DataPropertyInt Selection*)

Returns the selected integer property of the specified data.

5.2.4.18 float * getDataPtr (*DataHandle Data*)

Returns a pointer to the content of the specified data.

5.2.4.19 int getNumberOfApodizationRegions (*RawDataHandle Raw*)

Returns the number of regions in the raw data containing spectra that are supposed to be used for apodization.

5.2.4.20 int getNumberOfScanRegions (*RawDataHandle Raw*)

Returns the number of regions that have been acquired that contain scan data, i. e. spectra that are used to compute A-scans.

5.2.4.21 int getRawDataPropertyInt (*RawDataHandle RawData*, *RawDataPropertyInt Property*)

Returns a raw data property.

5.2.4.22 void * getRawDataPtr (*RawDataHandle RawDataSource*)

Returns the pointer to the raw data content. The pointer might no longer after additional actions using the *RawDataHandle*.

5.2.4.23 void getRawDataSize (*RawDataHandle Raw*, int * *SizeX*, int * *SizeY*, int * *SizeZ*)

Returns the size of the specified raw data ([RawDataHandle](#)).

5.2.4.24 void getScanSpectra (*RawDataHandle Raw*, int * *SpectralIndex*)

Returns the indices of spectra that contain scan data, i. e. spectra that are supposed to be used to compute A-scans.

An array needs to be provided that has twice the size of the number of scan regions which can be obtained by [getNumberOfScanRegions\(\)](#)

5.2.4.25 void reserveColoredData (*ColoredDataHandle ColData*, int *Size1*, int *Size2*, int *Size3*)

Reserves the amount of colored data specified. This might improve performance if appending data to the [ColoredDataHandle](#) as no additional memory needs to be reserved then.

5.2.4.26 void reserveComplexData (*ComplexDataHandle Data*, int *Size1*, int *Size2*, int *Size3*)

Reserves the amount of data specified. This might improve performance if appending data to the [ComplexDataHandle](#) as no additional memory needs to be reserved then.

5.2.4.27 void reserveData (*DataHandle Data*, int *Size1*, int *Size2*, int *Size3*)

Reserves the amount of data specified. This might improve performance if appending data to the [DataHandle](#) as no additional memory needs to be reserved then.

5.2.4.28 void resizeColoredData (*ColoredDataHandle ColData*, int *Size1*, int *Size2*, int *Size3*)

Resizes the respective colored data object. In general the data will be 1-dimensional if *Size2* and *Size3* are equal to 1, 2-dimensional if *Size3* is equal to 1 and 3-dimensional if all, *Size1*, *Size2*, *Size3*, are unequal to 1.

5.2.4.29 void resizeComplexData (ComplexDataHandle Data, int Size1, int Size2, int Size3)

Resizes the respective data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.

5.2.4.30 void resizeData (DataHandle Data, int Size1, int Size2, int Size3)

Resizes the respective data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.

5.2.4.31 void resizeRawData (RawDataHandle Raw, int Size1, int Size2, int Size3)

Resizes the specified raw data buffer accordingly.

5.2.4.32 void setApodizationSpectra (RawDataHandle RawData, int NumberOfScanRegions, int * ApodizationRegions)

Sets the number of the spectra in the raw data that contain data useful as apodization spectra.

5.2.4.33 void setColoredDataContent (ColoredDataHandle ColData, unsigned long * NewContent)

Sets the data content of the colored data object. The data chunk pointed to by NewContent needs to be of the size expected by the data object, i. e. Size1*Size2*Size*sizeof(unsigned long).

5.2.4.34 void setColoredDataOrientation (ColoredDataHandle Data, DataOrientation)

Sets the data orientation of the colored data object to the given orientation.

5.2.4.35 void setColoredDataRange (ColoredDataHandle Data, double range1, double range2, double range3)

Sets the range in mm in the 3 axes represented in the data object buffer.

5.2.4.36 void setComplexDataContent (ComplexDataHandle Data, ComplexFloat * NewContent)

Sets the data content of the [ComplexDataHandle](#) to the content specified by the pointer.

5.2.4.37 void setComplexDataRange (ComplexDataHandle Data, double range1, double range2, double range3)

Sets the range in mm in the 3 axes represented in the RealData buffer.

5.2.4.38 void setDataContent (DataHandle Data, float * NewContent)

Sets the data content of the data object. The data chunk pointed to by NewContent needs to be of the size expected by the data object, i. e. Size1*Size2*Size*sizeof(float).

5.2.4.39 void setDataOrientation (DataHandle Data, DataOrientation Orientation)

Sets the data orientation of the data object to the given orientation.

5.2.4.40 void setDataRange (DataHandle Data, double range1, double range2, double range3)

Sets the range in mm in the 3 axes represented in the RealData buffer.

5.2.4.41 void setRawDataBytesPerPixel (RawDataHandle Raw, int BytesPerPixel)

Sets the bytes per pixel for raw data.

5.2.4.42 void setRawDataContent (RawDataHandle RawDataSource, void * NewContent)

Sets the content of the raw data buffer. The size of the RawDataHandle needs to be adjusted first, as otherwise not all data might be copied.

5.2.4.43 void setScanSpectra (*RawDataHandle RawData*, int *NumberOfScanRegions*, int * *ScanRegions*)

Sets the number of the spectra in the raw data that are used for creating A-scan/B-scan data.

5.3 Data Creation and Clearing

Functions to create and clear object containing data.

Functions

- **SPECTRALRADAR_API RawDataHandle createRawData** (void)
Creates a raw data object ([RawDataHandle](#)).
- **SPECTRALRADAR_API void clearRawData** ([RawDataHandle](#) Raw)
Clears a raw data object ([RawDataHandle](#)).
- **SPECTRALRADAR_API DataHandle createData** (void)
Creates a 1-dimensional data object, containing floating point data.
- **SPECTRALRADAR_API DataHandle createGradientData** (int Size)
Creates a 1-dimensional data object, containing floating point data with equidistant arranged values between [0, size-1] with distance 1/(size-1).
- **SPECTRALRADAR_API void clearData** ([DataHandle](#) Data)
Clears the specified [DataHandle](#), [DataHandle](#), [DataHandle](#) or [DataHandle](#) objects.
- **SPECTRALRADAR_API ColoredDataHandle createColoredData** (void)
Creates a colored data object ([ColoredDataHandle](#)).
- **SPECTRALRADAR_API void clearColoredData** ([ColoredDataHandle](#) Volume)
Clears a colored volume object.
- **SPECTRALRADAR_API ComplexDataHandle createComplexData** (void)
Creates a data object holding complex data.
- **SPECTRALRADAR_API void clearComplexData** ([ComplexDataHandle](#) Data)
Clears a data object holding complex data.

5.3.1 Detailed Description

Functions to create and clear object containing data.

5.3.2 Function Documentation

5.3.2.1 void clearColoredData ([ColoredDataHandle](#) Volume)

Clears a colored volume object.

5.3.2.2 void clearComplexData ([ComplexDataHandle](#) Data)

Clears a data object holding complex data.

5.3.2.3 void clearData ([DataHandle](#) Data)

Clears the specified [DataHandle](#), [DataHandle](#), [DataHandle](#) or [DataHandle](#) objects.

5.3.2.4 void clearRawData ([RawDataHandle](#) Raw)

Clears a raw data object ([RawDataHandle](#)).

5.3.2.5 [ColoredDataHandle](#) createColoredData (void)

Creates a colored data object ([ColoredDataHandle](#)).

5.3.2.6 [ComplexDataHandle](#) createComplexData (void)

Creates a data object holding complex data.

5.3.2.7 DataHandle createData (void)

Creates a 1-dimensional data object, containing floating point data.

5.3.2.8 DataHandle createGradientData (int *Size*)

Creates a 1-dimensional data object, containing floating point data with equidistant arranged values between [0, size-1] with distance 1/(size-1).

5.3.2.9 RawDataHandle createRawData (void)

Creates a raw data object ([RawDataHandle](#)).

5.4 Hardware

Functions providing direct access to OCT Hardware functionality.

Typedefs

- typedef struct C_OCTDevice * [OCTDeviceHandle](#)

The OCTDeviceHandle type is used as Handle for using the SpectralRadar.

Enumerations

- enum [DevicePropertyFloat](#) {
[Device_FullWellCapacity](#),
[Device_zSpacing](#),
[Device_zRange](#),
[Device_SignalAmplitudeMin_dB](#),
[Device_SignalAmplitudeLow_dB](#),
[Device_SignalAmplitudeHigh_dB](#),
[Device_SignalAmplitudeMax_dB](#),
[Device_BinToElectronScaling](#),
[Device_Temperature](#),
[Device_SLD_OnTime_sec](#),
[Device_CenterWavelength_nm](#),
[Device_SpectralWidth_nm](#),
[Device_MaxTriggerFrequency_Hz](#) }

Properties of the device that can be read or measured.

- enum [DevicePropertyInt](#) {
[Device_SpectrumElements](#),
[Device_BytesPerElement](#),
[Device_MaxLiveVolumeRenderingScans](#),
[Device_BitDepth](#),
[Device_DataIsSigned](#) }

Properties of the device that can be read or measured.

- enum [ScanAxis](#) {
[ScanAxis_X](#) = 0,
[ScanAxis_Y](#) = 1 }

used to select the axis for manual galvo operations.

- enum [Device_CameraPreset](#) {
[Device_CameraPreset_Default](#),
[Device_CameraPreset_1](#),
[Device_CameraPreset_2](#),
[Device_CameraPreset_3](#),
[Device_CameraPreset_4](#) }

Enum identifying sensitivity and acquisition speed of the device.

Functions

- [SPECTRALRADAR_API OCTDeviceHandle initDevice](#) (void)
Initializes the installed device.
- [SPECTRALRADAR_API void getDeviceType](#) ([OCTDeviceHandle](#) Dev, char DevName[], int BufferSize)
Gives the name of the device type that is given by the [OCTDeviceHandle](#).
- [SPECTRALRADAR_API int getDeviceRevision](#) ([OCTDeviceHandle](#) Dev)
Returns the revision of the device given by the [OCTDeviceHandle](#).

- **SPECTRALRADAR_API** void **getDeviceSerialNumber** (**OCTDeviceHandle** Dev, char DevName[], int Buffer-Size)
Returns the serial number of the device given by the *OCTDeviceHandle*.
- **SPECTRALRADAR_API** int **getDevicePropertyInt** (**OCTDeviceHandle** Dev, **DevicePropertyInt** Selection)
Returns properties of the device belonging to the specified *OCTDeviceHandle*.
- **SPECTRALRADAR_API** double **getDevicePropertyFloat** (**OCTDeviceHandle** Dev, **DevicePropertyFloat** Selection)
Returns properties of the device belonging to the specified *OCTDeviceHandle*.
- **SPECTRALRADAR_API** void **closeDevice** (**OCTDeviceHandle** Dev)
Closes the device opened previously with *initDevice*.
- **SPECTRALRADAR_API** **BOOL** **isDeviceOn** (**OCTDeviceHandle** Handle)
Returns if the device is switched on.
- **SPECTRALRADAR_API** **BOOL** **isVideoCameraAvailable** (**OCTDeviceHandle** Dev)
Returns if the video camera is available.
- **SPECTRALRADAR_API** **BOOL** **isSLDAvailable** (**OCTDeviceHandle** Dev)
Returns whether the SLD is available.
- **SPECTRALRADAR_API** void **setSLD** (**OCTDeviceHandle** Dev, **BOOL** OnOff)
switches the SLD of the SpectralRadar device on and off.
- **SPECTRALRADAR_API** void **moveScanner** (**OCTDeviceHandle** Dev, **ProbeHandle** Probe, **ScanAxis** Axis, double Position)
manually moves the scanner to a given position
- **SPECTRALRADAR_API** void **setLaserDiode** (**OCTDeviceHandle** Dev, **BOOL** OnOff)
switches the LaserDiode of the SpectralRadar device on and off.
- **SPECTRALRADAR_API** double **getWavelengthAtPixel** (**OCTDeviceHandle** Dev, int Pixel)
Returns the wavelength at a specified pixel of the spectrometer.
- **SPECTRALRADAR_API** int **getCameraPreset** (**OCTDeviceHandle** Dev)
Gets the currently used device preset.

5.4.1 Detailed Description

Functions providing direct access to OCT Hardware functionality.

5.4.2 Typedef Documentation

5.4.2.1 OCTDeviceHandle

The OCTDeviceHandle type is used as Handle for using the SpectralRadar.

5.4.3 Enumeration Type Documentation

5.4.3.1 enum Device_CameraPreset

Enum identifying sensitivity and acquisition speed of the device.

Enumerator

- Device_CameraPreset_Default** Default device preset. Most common compromise of acquisition speed and sensitivity.
- Device_CameraPreset_1** Device preset 1.
- Device_CameraPreset_2** Device preset 2.
- Device_CameraPreset_3** Device preset 3.
- Device_CameraPreset_4** Device preset 4.

5.4.3.2 enum DevicePropertyFloat

Properties of the device that can be read or measured.

Enumerator

- Device_FullWellCapacity** The full well capacity of the device.
- Device_zSpacing** The spacing between two pixels in an A-scan.
- Device_zRange** The maximum measurement range for an A-scan.
- Device_SignalAmplitudeMin_dB** The minimum expected dB value for final data.
- Device_SignalAmplitudeLow_dB** The typical low dB value for final data.
- Device_SignalAmplitudeHigh_dB** The typical high dB value for final data.
- Device_SignalAmplitudeMax_dB** The maximum expected dB value for final data.
- Device_BinToElectronScaling** Scaling factor between binary raw data and electrons/photons.
- Device_Temperature** Internal device temperature in degrees C.
- Device_SLD_OnTime_sec** Absolute power-on time of the SLD since first start in seconds.
- Device_CenterWavelength_nm** The center wavelength of the device.
- Device_SpectralWidth_nm** The spectral width of the spectrometer.
- Device_MaxTriggerFrequency_Hz** Maximal valid trigger frequency depending on the chosen camera preset.

5.4.3.3 enum DevicePropertyInt

Properties of the device that can be read or measured.

Enumerator

- Device_SpectrumElements** The number of pixels provided by the spectrometer.
- Device_BytesPerElement** The number of bytes one element of the spectrum occupies.
- Device_MaxLiveVolumeRenderingScans** The maximum number of scans per dimension in the live volume rendering mode.
- Device_BitDepth** Bit depth of the DAQ.
- Device_DataIsSigned** Flag indicating if the data is signed.

5.4.3.4 enum ScanAxis

used to select the axis for manual galvo operations.

Enumerator

- ScanAxis_X** X-Axis of the scanner.
- ScanAxis_Y** Y-Axis of the scanner.

5.4.4 Function Documentation

5.4.4.1 void closeDevice (OCTDeviceHandle Dev)

Closes the device opened previously with initDevice.

Parameters

<i>Dev</i>	The OCTDeviceHandle that was initially provided by <code>initDevice</code> .
------------	--

5.4.4.2 void Device_CameraPreset getCameraPreset (OCTDeviceHandle Dev)

Gets the currently used device preset.

5.4.4.3 double getDevicePropertyFloat (OCTDeviceHandle Dev, DevicePropertyFloat Selection)

Returns properties of the device belonging to the specified [OCTDeviceHandle](#).

5.4.4.4 int getDevicePropertyInt (OCTDeviceHandle Dev, DevicePropertyInt Selection)

Returns properties of the device belonging to the specified [OCTDeviceHandle](#).

5.4.4.5 int getDeviceRevision (OCTDeviceHandle Dev)

Returns the revision of the device given by the [OCTDeviceHandle](#).

5.4.4.6 int void getDeviceSerialNumber (OCTDeviceHandle Dev, char DevName[], int BufferSize)

Returns the serial number of the device given by the [OCTDeviceHandle](#).

5.4.4.7 void getDeviceType (OCTDeviceHandle Dev, char DevName[], int BufferSize)

Gives the name of the device type that is given by the [OCTDeviceHandle](#).

5.4.4.8 double getWavelengthAtPixel (OCTDeviceHandle Dev, int Pixel)

Returns the wavelength at a speicified pixel of the spectrometer.

Warning

This function is still experimental and results might be incorrect.

5.4.4.9 OCTDeviceHandle initDevice (void)

Initializes the installed device.

Returns

Handle to the initialized OCT device.

5.4.4.10 BOOL isDeviceOn (OCTDeviceHandle Handle)

Returns if the device is switched on.

Parameters

<i>Dev</i>	The OCTDeviceHandle that was initially provided by <code>initDevice</code> .
------------	--

5.4.4.11 BOOL isSLDAvailable (OCTDeviceHandle Dev)

Returns whethther the SLD is available.

Parameters

<i>Dev</i>	The OCTDeviceHandle that was initially provided by initDevice.
------------	--

5.4.4.12 **BOOL** isVideoCameraAvailable ([OCTDeviceHandle](#) *Dev*)

Returns if the video camera is available.

Parameters

<i>Dev</i>	The OCTDeviceHandle that was initially provided by initDevice.
------------	--

5.4.4.13 **void** moveScanner ([OCTDeviceHandle](#) *Dev*, [ProbeHandle](#) *Probe*, [ScanAxis](#) *Axis*, double *Position*)

manually moves the scanner to a given position

Parameters

<i>Dev</i>	the OCTDeviceHandle that was initially provided by initDevice.
<i>Probe</i>	A handle to the probe (ProbeHandle); whose galvo position is to be set.
<i>Axis</i>	the axis in which you want to set the position manually
<i>Position</i>	the actual position you want to move the galvo to.

5.4.4.14 **void** setLaserDiode ([OCTDeviceHandle](#) *Dev*, **BOOL** *OnOff*)

switches the LaserDiode of the SpectralRadar device on and off.

Warning

Not all devices are equipped

Parameters

<i>Dev</i>	handle to the OCTDeviceHandle that was initially provided by initDevice.
<i>OnOff</i>	TRUE switches the VisLD on, FALSE swichted the VisLD off.

5.4.4.15 **void** setSLD ([OCTDeviceHandle](#) *Dev*, **BOOL** *OnOff*)

switches the SLD of the SpectralRadar device on and off.

Parameters

<i>Dev</i>	handle to the OCTDeviceHandle that was initially provided by initDevice.
<i>OnOff</i>	TRUE switches the SLD on, FALSE swichted the SLD off.

5.5 Internal Values

Functions for access to all kinds of Digital-to-Analog and Analog-to-Digital on the device.

Functions

- **SPECTRALRADAR_API** int [getNumberOfInternalValues](#) (OCTDeviceHandle Dev)
Returns the number of Analog-to-Digital Converter present in the device.
- **SPECTRALRADAR_API** void [getInternalValueName](#) (OCTDeviceHandle Dev, int Index, char *Name, int NameStringSize, char *Unit, int UnitStringSize)
Returns names and unit for the specified Analog-to-Digital Converter.
- **SPECTRALRADAR_API** double [getInternalValueByName](#) (OCTDeviceHandle Dev, const char *Name)
Returns the value of the specified Analog-to-Digital Converter (ADC);.
- **SPECTRALRADAR_API** double [getInternalValueByIndex](#) (OCTDeviceHandle Dev, int Index)
Returns the value of the selected ADC.

5.5.1 Detailed Description

Functions for access to all kinds of Digital-to-Analog and Analog-to-Digital on the device.

5.5.2 Function Documentation

5.5.2.1 double [getInternalValueByIndex](#) (OCTDeviceHandle Dev, int Index)

Returns the value of the selected ADC.

The index is running number, starting with 0, smaller than the number specified by [getNumberOfInternalValues](#).

5.5.2.2 double [getInternalValueByName](#) (OCTDeviceHandle Dev, const char * Name)

Returns the value of the specified Analog-to-Digital Converter (ADC);.

The ADC is specified by the name returned by [getInternalValueName](#).

5.5.2.3 void [getInternalValueName](#) (OCTDeviceHandle Dev, int Index, char * Name, int NameStringSize, char * Unit, int UnitStringSize)

Returns names and unit for the specified Analog-to-Digital Converter.

The index is running number, starting with 0, smaller than the number specified by [getNumberOfInternalValues](#).

5.5.2.4 int [getNumberOfInternalValues](#) (OCTDeviceHandle Dev)

Returns the number of Analog-to-Digital Converter present in the device.

5.6 Pattern Factory/Probe

Functions setting up a probe that can be used to create scan patterns.

Typedefs

- typedef struct C_Probe * [ProbeHandle](#)
Handle for controlling the galvo scanner.

Enumerations

- enum [ProbeParameterFloat](#) {
[Probe_FactorX](#),
[Probe_OffsetX](#),
[Probe_FactorY](#),
[Probe_OffsetY](#),
[Probe_FlybackTime_Sec](#),
[Probe_ExpansionTime_Sec](#),
[Probe_RotationTime_Sec](#),
[Probe_ExpectedScanRate_Hz](#),
[Probe_CameraScalingX](#),
[Probe_CameraOffsetX](#),
[Probe_CameraScalingY](#),
[Probe_CameraOffsetY](#),
[Probe_CameraAngle](#),
[Probe_WhiteBalanceRed](#),
[Probe_WhiteBalanceGreen](#),
[Probe_WhiteBalanceBlue](#),
[Probe_RangeMaxX](#),
[Probe_RangeMaxY](#),
[Probe_MaximumSlope_XY](#),
[Probe_SpeckleSize](#),
[Probe_ApoPosX](#),
[Probe_ApoPosY](#),
[Probe_ReferenceStageOffset](#) }
Parameters describing the behaviour of the Probe, such as calibration factors and scan parameters.
- enum [ProbeParameterInt](#) {
[Probe_ApodizationCycles](#),
[Probe_Oversampling](#),
[Probe_WhiteBalanceAutomatic](#),
[Probe_Oversampling_SlowAxis](#),
[Probe_SpeckleReduction](#),
[Probe_MaxScanRangeShape](#) }
Parameters describing the behaviour of the Probe, such as calibration factors and scan parameters.
- enum [ProbeFlag](#) {
[Probe_CameraInverted_X](#),
[Probe_CameraInverted_Y](#),
[Probe_HasMEMSScanner](#) }
Boolean parameters describing the behaviour of the Probe.
- enum [ProbeType](#) {
[ProbeType_Standard](#),
[ProbeType_Handheld](#),
[ProbeType_Scientific](#) }
Determines the kind of probe types.

Functions

- **SPECTRALRADAR_API ProbeHandle initProbe** (OCTDeviceHandle Dev, const char *ProbeFile)
Initializes a probe specified by ProbeFile.
- **SPECTRALRADAR_API ProbeHandle initStandardProbe** (OCTDeviceHandle Dev)
Creates a standard probe using the Probe.ini file. If this configuration file is not found, standard parameters without valid calibration will be used.
- **SPECTRALRADAR_API ProbeHandle initProbeWithType** (OCTDeviceHandle Dev, ProbeType Type)
Creates a standard probe for the given probe type but without valid calibration data .
- **SPECTRALRADAR_API void saveProbe** (ProbeHandle Probe, const char *ProbeFile)
Saves the current properties of the ProbeHandle to a specified INI file to be reloaded using the initProbe() function.
- **SPECTRALRADAR_API void setProbeParameterInt** (ProbeHandle Probe, ProbeParameterInt Selection, int Value)
Sets.
- **SPECTRALRADAR_API void setProbeParameterFloat** (ProbeHandle Probe, ProbeParameterFloat Selection, double Value)
Sets floating point parameters of the specified probe.
- **SPECTRALRADAR_API int getProbeParameterInt** (ProbeHandle Probe, ProbeParameterInt Selection)
Gets integer parameters of the specified probe.
- **SPECTRALRADAR_API double getProbeParameterFloat** (ProbeHandle Probe, ProbeParameterFloat Selection)
Gets floating point parameters of the specified probe.
- **SPECTRALRADAR_API BOOL getProbeFlag** (ProbeHandle Probe, ProbeFlag Selection)
Returns the selected boolean value of the specified probe.
- **SPECTRALRADAR_API void getProbeName** (ProbeHandle Probe, char ProbeName[], int BufferSize)
Returns the name of the specified probe.
- **SPECTRALRADAR_API void setProbeName** (ProbeHandle Probe, const char *ProbeName)
Sets the given name of the specified probe.
- **SPECTRALRADAR_API void getProbeSerialNo** (ProbeHandle Probe, char SerialNo[], int BufferSize)
Gets the serial number of the specified probe.
- **SPECTRALRADAR_API void setProbeSerialNo** (ProbeHandle Probe, const char *SerialNo)
Gets the serial number of the specified probe.
- **SPECTRALRADAR_API void getProbeType** (ProbeHandle Probe, char Type[], int BufferSize)
Gets the type of the specified probe.
- **SPECTRALRADAR_API void setProbeType** (ProbeHandle Probe, const char *Type)
Sets the type of the specified probe.
- **SPECTRALRADAR_API void getProbeObjective** (ProbeHandle Probe, char Objective[], int BufferSize)
Gets the objective of the specified probe.
- **SPECTRALRADAR_API void setProbeObjective** (ProbeHandle Probe, const char *Objective)
Sets the given objective of the specified probe.
- **SPECTRALRADAR_API void closeProbe** (ProbeHandle Probe)
Closes the probe and frees all memory associated with it.
- **SPECTRALRADAR_API void blendEnFaceInCamera** (ProbeHandle Probe, ScanPatternHandle Pattern, ColoredDataHandle EnFace2D, ColoredDataHandle Image, float Ratio, BOOL DenseView)
Blends the en-face image of a given volume acquisition on top of the video image. Can be used to calibrate the probe manually.
- **SPECTRALRADAR_API void CameraPixelToPosition** (ProbeHandle Probe, ColoredDataHandle Image, int PixelX, int PixelY, double *PosX, double *PosY)
Computes the physical position of a camera pixel of the video camera in the probe. It needs to be assured that the device is properly calibrated.
- **SPECTRALRADAR_API void PositionToCameraPixel** (ProbeHandle Probe, ColoredDataHandle Image, double PosX, double PosY, int *PixelX, int *PixelY)
Computes the pixel of the video camera corresponding to a physical position. It needs to be assured that the device is properly calibrated.

5.6.1 Detailed Description

Functions setting up a probe that can be used to create scan patterns.

5.6.2 Typedef Documentation

5.6.2.1 ProbeHandle

Handle for controlling the galvo scanner.

5.6.3 Enumeration Type Documentation

5.6.3.1 enum ProbeFlag

Boolean parameters describing the behaviour of the Probe.

Enumerator

Probe_CameraInverted_X Bool if the scan pattern in the video camera image is flipped around x-axis or not.

Probe_CameraInverted_Y Bool if the scan pattern in the video camera image is flipped around y-axis or not.

Probe_HasMEMSScanner Boolean if the probe type uses a MEMS mirror or not, e.g. a handheld probe.

5.6.3.2 enum ProbeParameterFloat

Parameters describing the behaviour of the Probe, such as calibration factors and scan parameters.

Computation of physical position and raw values for the scanner is done by $\text{PhysicalPosition} = \text{Factor} * \text{RawValue} + \text{Offset}$

Enumerator

Probe_FactorX Factor for the x axis.

Probe_OffsetX Offset for the x axis.

Probe_FactorY Factor for the y axis.

Probe_OffsetY Offset for the y axis.

Probe_FlybackTime_Sec Flyback time of the system. This time is usually needed to get from an apodization position to scan position and vice versa.

Probe_ExpansionTime_Sec The scanning range is extended by a number of A-scans equivalent to the expansion time.

Probe_RotationTime_Sec The scan pattern is usually shifted by a number of A-scans equivalent to the rotation time.

Probe_ExpectedScanRate_Hz The expected scan rate.

Warning

In general the expected scan rate is set during initialization of the probe with respect to the attached device. In most cases it should not be altered manually.

Probe_CameraScalingX The px/mm ratio in X direction for the BScan overlay on the video image.

Probe_CameraOffsetX The BScan overlay X offset in pixels.

Probe_CameraScalingY The px/mm ratio in Y direction for the BScan overlay on the video image.

Probe_CameraOffsetY The BScan overlay Y offset in pixels.

Probe_CameraAngle Corrective rotation angle for the BScan overlay.

Probe_WhiteBalanceRed White balance settings will only take effect on initialization of the probe. White balance value for red channel of video camera (if -1, adjustment will be automatic).

Probe_WhiteBalanceGreen White balance value for green channel of video camera (if -1, adjustment will be automatic).

Probe_WhiteBalanceBlue White balance value for blue channel of video camera (if -1, adjustment will be automatic).

Probe_RangeMaxX Maximum scan range in X direction.

Probe_RangeMaxY Maximum scan range in Y direction.

Probe_MaximumSlope_XY Maximum galvo slope (accounting for the distortion capabilities of different galvo types)

Probe_SpeckleSize Speckle size to be used for scan pattern computation if speckle reduction is switched on.

Probe_ApoPosX X-position used to acquire the apodization spectrum.

Probe_ApoPosY Y-position used to acquire the apodization spectrum.

Probe_ReferenceStageOffset Offset for reference stage marking the zero delay line.

5.6.3.3 enum ProbeParameterInt

Parameters describing the behaviour of the Probe, such as calibration factors and scan parameters.

Enumerator

Probe_ApodizationCycles The number of cycles used for apodization.

Probe_Oversampling A factor used as oversampling.

Probe_WhiteBalanceAutomatic Automatic white balance for video camera, 0 == off, not 0 == on.

Probe_Oversampling_SlowAxis A factor used as oversampling of the slow scanner axis.

Probe_SpeckleReduction Number of speckles that are scanned over for averaging. Requires Oversampling \geq SpeckleReduction.

Probe_MaxScanRangeShape Shape of the maximum scan range: 0 is a rectangle, 1 is an ellipse.

5.6.3.4 enum ProbeType

Determines the kind of probe types.

Enumerator

ProbeType_Standard Specifies the standard or general probe.

ProbeType_Handheld Specifies the handheld probe.

ProbeType_Scientific Specifies the scientific probe.

5.6.4 Function Documentation

5.6.4.1 void blendEnFaceInCamera (ProbeHandle Probe, ScanPatternHandle Pattern, ColoredDataHandle EnFace2D, ColoredDataHandle Image, float Ratio, BOOL DenseView)

Blends the en-face image of a given volume acquisition on top of the video image. Can be used to calibrate the probe manually.

5.6.4.2 void CameraPixelToPosition (ProbeHandle Probe, ColoredDataHandle Image, int PixelX, int PixelY, double * PosX, double * PosY)

Computes the physical position of a camera pixel of the video camera in the probe. It needs to be assured that the device is properly calibrated.

5.6.4.3 void closeProbe (ProbeHandle Probe)

Closes the probe and frees all memory associated with it.

5.6.4.4 BOOL getProbeFlag (ProbeHandle Probe, ProbeFlag Selection)

Returns the selected boolean value of the specified probe.

5.6.4.5 void getProbeName (ProbeHandle Probe, char ProbeName[], int BufferSize)

Returns the name of the specified probe.

5.6.4.6 void getProbeObjective (ProbeHandle Probe, char Objective[], int BufferSize)

Gets the objective of the specified probe.

5.6.4.7 double getProbeParameterFloat (ProbeHandle Probe, ProbeParameterFloat Selection)

Gets floating point parameters of the specified probe.

5.6.4.8 int getProbeParameterInt (ProbeHandle Probe, ProbeParameterInt Selection)

Gets integer parameters of the specified probe.

5.6.4.9 void getProbeSerialNo (ProbeHandle Probe, char SerialNo[], int BufferSize)

Gets the serial number of the specified probe.

5.6.4.10 void getProbeType (ProbeHandle Probe, char Type[], int BufferSize)

Gets the type of the specified probe.

5.6.4.11 ProbeHandle initProbe (OCTDeviceHandle Dev, const char * ProbeFile)

Initializes a probe specified by ProbeFile.

In older systems up until a manufacturing date of May 2011 either "Handheld" or "Microscope" are used. An according ini-file (i. e. "Handheld.ini" or "Microscope.ini"); will be loaded from the config path of the SpectralRadar installation containing all necessary information. With systems manufactured after May 2011 "Probe" should be used.

It is recommended to use #initStandardProbe for systems manufactured in or after May 2011.

5.6.4.12 ProbeHandle initProbeWithType (OCTDeviceHandle Dev, ProbeType Type)

Creates a standard probe for the given probe type but without valid calibration data .

5.6.4.13 ProbeHandle initStandardProbe (OCTDeviceHandle Dev)

Creates a standard probe using the Probe.ini file. If this configuration file is not found, standard parameters without valid calibration will be used.

5.6.4.14 void PositionToCameraPixel (ProbeHandle Probe, ColoredDataHandle Image, double PosX, double PosY, int * PixelX, int * PixelY)

Computes the pixel of the video camera corresponding to a physical position. It needs to be assured that the device is properly calibrated.

5.6.4.15 void saveProbe (ProbeHandle Probe, const char * ProbeFile)

Saves the current properties of the [ProbeHandle](#) to a specified INI file to be reloaded using the [initProbe\(\)](#) function.

5.6.4.16 void setProbeName (ProbeHandle *Probe*, const char * *ProbeName*)

Sets the given name of the specified probe.

5.6.4.17 void setProbeObjective (ProbeHandle *Probe*, const char * *Objective*)

Sets the given objective of the specified probe.

5.6.4.18 void setProbeParameterFloat (ProbeHandle *Probe*, ProbeParameterFloat *Selection*, double *Value*)

Sets floating point parameters of the specified probe.

5.6.4.19 void setProbeParameterInt (ProbeHandle *Probe*, ProbeParameterInt *Selection*, int *Value*)

Sets.

5.6.4.20 void setProbeSerialNo (ProbeHandle *Probe*, const char * *SerialNo*)

Gets the serial number of the specified probe.

5.6.4.21 void setProbeType (ProbeHandle *Probe*, const char * *Type*)

Sets the type of the specified probe.

5.7 Scan Pattern

Functions that describe the movement of the Scanner during measurement.

Typedefs

- typedef C_ScanPattern * [ScanPatternHandle](#)
Handle for controlling the scan pattern.

Functions

- [SPECTRALRADAR_API ScanPatternHandle createPointScanPattern](#) ([ProbeHandle](#) Probe, int Size, double PosX, double PosY)
Creates a scan pattern that consists of a single point (PosX, PosY). The galvo doesn't move from there. Use this pattern for point scans and/or non-scanning probes.
- [SPECTRALRADAR_API ScanPatternHandle createNoScanPattern](#) ([ProbeHandle](#) Probe, int Scans, int NumberOfScans)
Creates a simple scan pattern that does not move the galvo. Use this pattern for point scans and/or non-scanning probes.
- [SPECTRALRADAR_API ScanPatternHandle createTriggerPattern](#) ([ProbeHandle](#) Probe, int Scans)
Creates a pattern only consisting of a specified amount of trigger signals.
- [SPECTRALRADAR_API ScanPatternHandle createBScanPattern](#) ([ProbeHandle](#) Probe, double Range, int AScans, [BOOL](#) apodization)
Creates a simple B-scan pattern that moves the galvo over a specified range.
- [SPECTRALRADAR_API ScanPatternHandle createBilateralBScanPattern](#) ([ProbeHandle](#) Probe, double Range, int AScans, double Shift)
Creates a bilateral scan pattern. The contouring error can be influenced using the Shift parameter.
- [SPECTRALRADAR_API ScanPatternHandle createBScanPatternManual](#) ([ProbeHandle](#) Probe, double StartX, double StartY, double StopX, double StopY, int AScans, [BOOL](#) apodization)
Creates a B-scan pattern specified by start and end points.
- [SPECTRALRADAR_API ScanPatternHandle createIdealBScanPattern](#) ([ProbeHandle](#) Probe, double Range, int AScans)
Creates an ideal B-scan pattern assuming scanners with infinite speed. No correction factors are taken into account. This is only used for internal purposes and not as a scan pattern designed to be output to the galvo drivers.
- [SPECTRALRADAR_API ScanPatternHandle createCirclePattern](#) ([ProbeHandle](#) Probe, double Radius, int AScans)
Creates a circle scan pattern.
- [SPECTRALRADAR_API ScanPatternHandle createVolumePattern](#) ([ProbeHandle](#) Probe, double RangeX, int SizeX, double RangeY, int SizeY)
Creates a simple volume pattern.
- [SPECTRALRADAR_API ScanPatternHandle createBScanStackPattern](#) ([ProbeHandle](#) Probe, double RangeX, int SizeX, double RangeY, int SizeY)
Creates a simple stack pattern.
- [SPECTRALRADAR_API ScanPatternHandle createFreeformScanPattern](#) ([ProbeHandle](#) Probe, float *positions, int size_x, int size_y, [BOOL](#) apodization)
Creates a freeform scan pattern based on an array of positions.
- [SPECTRALRADAR_API ScanPatternHandle createFragmentedScanPattern](#) ([ProbeHandle](#) Probe, int ChunkSize, int NumberOfChunks)
Creates a scan pattern which can be used to acquire a dataset of <NumberOfChunks> times <ChunkSize> A-scans at position 0/0. The Fragmented scan pattern can be compared in structure to a B-scan stack pattern with x and y ranges of 0; however the fragmented scan pattern behaves like a volume pattern in that it shows no delay between the respective chunks.
- [SPECTRALRADAR_API void updateScanPattern](#) ([ScanPatternHandle](#) Pattern)

- Updates the specified pattern (ScanPatternHandle);.*
- **SPECTRALRADAR_API** void **rotateScanPattern** (ScanPatternHandle Pattern, double Angle)
- Rotates the specified pattern (ScanPatternHandle);.*
- **SPECTRALRADAR_API** void **rotateScanPatternExt** (ScanPatternHandle Pattern, double Angle, int index)
- Rotates the scan #index (0-based) of the specified pattern (ScanPatternHandle).*
- **SPECTRALRADAR_API** void **shiftScanPattern** (ScanPatternHandle Pattern, double ShiftX, double ShiftY)
- Shifts the specified pattern (ScanPatternHandle).*
- **SPECTRALRADAR_API** void **shiftScanPatternExt** (ScanPatternHandle Pattern, double ShiftX, double ShiftY, BOOL ShiftApo, int Index)
- Shifts the scan #index (0-based) of the specified pattern (ScanPatternHandle).*
- **SPECTRALRADAR_API** void **zoomScanPattern** (ScanPatternHandle Pattern, double Factor)
- Zooms the specified pattern (ScanPatternHandle).*
- **SPECTRALRADAR_API** int **getScanPatternLUTSize** (ScanPatternHandle Pattern)
- Returns the number of data points the specified pattern (ScanPatternHandle) used.*
- **SPECTRALRADAR_API** void **getScanPatternLUT** (ScanPatternHandle Pattern, double *PosX, double *PosY)
- Returns the actual positions to be scanned with the specified pattern (ScanPatternHandle).*
- **SPECTRALRADAR_API** void **clearScanPattern** (ScanPatternHandle Pattern)
- Clears the specified scan pattern (ScanPatternHandle).*

5.7.1 Detailed Description

Functions that describe the movement of the Scanner during measurement.

5.7.2 Typedef Documentation

5.7.2.1 ScanPatternHandle

Handle for controlling the scan pattern.

5.7.3 Function Documentation

5.7.3.1 void clearScanPattern (ScanPatternHandle Pattern)

Clears the specified scan pattern (ScanPatternHandle).

5.7.3.2 ScanPatternHandle createBilateralBScanPattern (ProbeHandle Probe, double Range, int AScans, double Shift)

Creates a bilateral scan pattern. The contouring error can be influenced using the Shift parameter.

5.7.3.3 ScanPatternHandle createBScanPattern (ProbeHandle Probe, double Range, int AScans, BOOL apodization)

Creates a simple B-scan pattern that moves the galvo over a specified range.

5.7.3.4 ScanPatternHandle createBScanPatternManual (ProbeHandle Probe, double StartX, double StartY, double StopX, double StopY, int AScans, BOOL apodization)

Creates a B-scan pattern specified by start and end points.

5.7.3.5 ScanPatternHandle createBScanStackPattern (ProbeHandle Probe, double RangeX, int SizeX, double RangeY, int SizeY)

Creates a simple stack pattern.

The BScan stack pattern is a volume measurement which consists of several shifted B-Scan measurements. The resulting data will be identical to a volume (see [createVolumePattern\(\)](#)) but an apodization is performed for each slice (B-scan). The volume will be returned slice-by-slice by calling [getRawData\(\)](#).

5.7.3.6 **ScanPatternHandle** createCirclePattern (**ProbeHandle** *Probe*, double *Radius*, int *AScans*)

Creates a circle scan pattern.

Warning

Circle patterns cannot be rotated properly.

5.7.3.7 **ScanPatternHandle** createFragmentedScanPattern (**ProbeHandle** *Probe*, int *ChunkSize*, int *NumberOfChunks*)

Creates a scan pattern which can be used to acquire a dataset of <NumberOfChunks> times <ChunkSize> A-scans at position 0/0. The Fragmented scan pattern can be compared in structure to a B-scan stack pattern with x and y ranges of 0; however the fragmented scan pattern behaves like a volume pattern in that it shows no delay between the respective chunks.

5.7.3.8 **ScanPatternHandle** createFreeformScanPattern (**ProbeHandle** *Probe*, float * *positions*, int *size_x*, int *size_y*, **BOOL** *apodization*)

Creates a freeform scan pattern based on an array of positions.

The positions array must consist of pairs of x/y coordinates inside the valid scanning limits of the probe and contain all points of the scan pattern. *size_x/size_y* can be used analogue to the B-scan stack pattern to discern the length (*size_x*) and number (*size_y*) of single lines or sections inside the freeform scan pattern. The position array is taken as-is, so care must be taken to use sensible values with regard to the capabilities of the utilized scanner system and to the resolution of the system resp. the desired resolution of your scan pattern.

5.7.3.9 **ScanPatternHandle** createIdealBScanPattern (**ProbeHandle** *Probe*, double *Range*, int *AScans*)

Creates an ideal B-scan pattern assuming scanners with infinite speed. No correction factors are taken into account. This is only used for internal purposes and not as a scan pattern designed to be output to the galvo drivers.

5.7.3.10 **ScanPatternHandle** createNoScanPattern (**ProbeHandle** *Probe*, int *Scans*, int *NumberOfScans*)

Creates a simple scan pattern that does not move the galvo. Use this pattern for point scans and/or non-scanning probes.

5.7.3.11 **ScanPatternHandle** createPointScanPattern (**ProbeHandle** *Probe*, int *Size*, double *PosX*, double *PosY*)

Creates a scan pattern that consists of a single point (PosX, PosY). The galvo doesn't move from there. Use this pattern for point scans and/or non-scanning probes.

5.7.3.12 **ScanPatternHandle** createTriggerPattern (**ProbeHandle** *Probe*, int *Scans*)

Creates a pattern only consisting of a specified amount of trigger signals.

5.7.3.13 **ScanPatternHandle** createVolumePattern (**ProbeHandle** *Probe*, double *RangeX*, int *SizeX*, double *RangeY*, int *SizeY*)

Creates a simple volume pattern.

The volume pattern consists of a single uninterrupted scan and all data is acquired in a single measurement. In contrast to a B-scan stack pattern (see [createBScanStackPattern\(\)](#)) only one apodization is performed for the complete volume. The complete volume will be returned in one raw data ([RawDataHandle](#)).

5.7.3.14 void getScanPatternLUT (**ScanPatternHandle** *Pattern*, double * *PosX*, double * *PosY*)

Returns the actual positions to be scanned with the specified pattern ([ScanPatternHandle](#)).

5.7.3.15 int getScanPatternLUTSize (**ScanPatternHandle** *Pattern*)

Returns the number of data points the specified pattern ([ScanPatternHandle](#)) used.

5.7.3.16 void rotateScanPattern (ScanPatternHandle *Pattern*, double *Angle*)

Rotates the specified pattern ([ScanPatternHandle](#));.

5.7.3.17 void rotateScanPatternExt (ScanPatternHandle *Pattern*, double *Angle*, int *index*)

Rotates the scan #index (0-based) of the specified pattern ([ScanPatternHandle](#)).

5.7.3.18 void shiftScanPattern (ScanPatternHandle *Pattern*, double *ShiftX*, double *ShiftY*)

Shifts the specified pattern ([ScanPatternHandle](#)).

5.7.3.19 void shiftScanPatternExt (ScanPatternHandle *Pattern*, double *ShiftX*, double *ShiftY*, BOOL *ShiftApo*, int *Index*)

Shifts the scan #index (0-based) of the specified pattern ([ScanPatternHandle](#)).

5.7.3.20 void updateScanPattern (ScanPatternHandle *Pattern*)

Updates the specified pattern ([ScanPatternHandle](#));.

5.7.3.21 void zoomScanPattern (ScanPatternHandle *Pattern*, double *Factor*)

Zooms the specified pattern ([ScanPatternHandle](#)).

5.8 Acquisition

Functions for acquisition.

Enumerations

- enum `AcquisitionType` {
`Acquisition_AsyncContinuous`,
`Acquisition_AsyncFinite`,
`Acquisition_Sync` }

Determines the kind of acquisition process.

Functions

- `SPECTRALRADAR_API` void `startMeasurement` (`OCTDeviceHandle` Dev, `ScanPatternHandle` Pattern, `AcquisitionType` type)
starts a continuous measurement BScans.
- `SPECTRALRADAR_API` void `getRawData` (`OCTDeviceHandle` Dev, `RawDataHandle` RawData)
Acquires data and stores the data unprocessed.
- `SPECTRALRADAR_API` void `getRawDataEx` (`OCTDeviceHandle` Dev, `RawDataHandle` RawData, int CameraIdx)
Acquires data with the specific camera given with camera index and stores the data unprocessed.
- `SPECTRALRADAR_API` void `stopMeasurement` (`OCTDeviceHandle` Dev)
stops the current measurement.
- `SPECTRALRADAR_API` void `measureSpectra` (`OCTDeviceHandle` Dev, int NumberOfSpectra, `RawDataHandle` Raw)
Acquires N spectra of raw data without moving galvo scanners.
- `SPECTRALRADAR_API` void `measureSpectraEx` (`OCTDeviceHandle` Dev, int N, `RawDataHandle` Raw, int CameraIndex)
Acquires N spectra of raw data without moving galvo scanners. Supports multiple cameras (e.g. PS-OCT).

5.8.1 Detailed Description

Functions for acquisition.

5.8.2 Enumeration Type Documentation

5.8.2.1 enum `AcquisitionType`

Determines the kind of acquisition process.

Enumerator

`Acquisition_AsyncContinuous` Specifies an asynchronous infinite/continuous measurement. The internal memory management is constructed with a loop of buffers. Note that you may lose data if the acquisition is faster than the collection with `getRawData()`.

`Acquisition_AsyncFinite` Specifies an asynchronous finite measurement. The reserved memory for the data is as large as required to make sure that no data will be lost.

`Acquisition_Sync` Specifies a synchronous measurement. Each measurement of e.g. a single B-scan will be started with `getRawData()`. Since the acquisition is synchronized with the software it is not a continuous measurement over the time of e.g. several B-scans.

5.8.3 Function Documentation

5.8.3.1 void getRawData (OCTDeviceHandle Dev, RawDataHandle RawData)

Acquires data and stores the data unprocessed.

5.8.3.2 void void getRawDataEx (OCTDeviceHandle Dev, RawDataHandle RawData, int CameraIdx)

Acquires data with the specific camera given with camera index and stores the data unprocessed.

5.8.3.3 void measureSpectra (OCTDeviceHandle Dev, int NumberOfSpectra, RawDataHandle Raw)

Acquires N spectra of raw data without moving galvo scanners.

5.8.3.4 void measureSpectraEx (OCTDeviceHandle Dev, int N, RawDataHandle Raw, int CameraIndex)

Acquires N spectra of raw data without moving galvo scanners. Supports multiple cameras (e.g. PS-OCT).

5.8.3.5 void startMeasurement (OCTDeviceHandle Dev, ScanPatternHandle Pattern, AcquisitionType type)

starts a continuous measurement BScans.

Scanning takes place according to the specified scan pattern handle. Data can be recorded using the [getRawData\(\)](#) function. If you are done, call [stopMeasurement\(\)](#).

Parameters

<i>Dev</i>	The OCTDeviceHandle that was initially provided by <code>initDevice</code> .
<i>Pattern</i>	the ScanPatternHandle

5.8.3.6 void stopMeasurement (OCTDeviceHandle Dev)

stops the current measurement.

Parameters

<i>Dev</i>	The OCTDeviceHandle that was initially provided by <code>initDevice</code> .
------------	--

5.9 Processing

Standard Processing Routines.

Typedefs

- typedef struct C_Processing * [ProcessingHandle](#)
Handle for a processing routine.

Enumerations

- enum [ProcessingType](#) {
 [Processing_StandardFFT](#),
 [Processing_StandardNDFT](#),
 [Processing_iFFT1](#),
 [Processing_iFFT2](#),
 [Processing_iFFT3](#),
 [Processing_iFFT4](#),
 [Processing_NFFT1](#),
 [Processing_NFFT2](#),
 [Processing_NFFT3](#),
 [Processing_NFFT4](#) }
defines the algorithm used for dechirping the input signal and Fourier transformation
- enum [ApodizationWindow](#) {
 [Apodization_Hann](#) = 0,
 [Apodization_Hamming](#) = 1,
 [Apodization_Gauss](#) = 2,
 [Apodization_TaperedCosine](#) = 3,
 [Apodization_Blackman](#) = 4,
 [Apodization_BlackmanHarris](#) = 5,
 [Apodization_LightSourceBased](#) = 6,
 [Apodization_Unknown](#) = 999 }
To select the apodization window function.
- enum [ProcessingParameterInt](#) {
 [Processing_SpectrumAveraging](#),
 [Processing_AScanAveraging](#),
 [Processing_BScanAveraging](#),
 [Processing_ZeroPadding](#),
 [Processing_NumberOfThreads](#),
 [Processing_FourierAveraging](#) }
Parameters that set the behaviour of the processing algorithms.
- enum [ProcessingParameterFloat](#) {
 [Processing_ApodizationDamping](#),
 [Processing_MinElectrons](#) }
Parameters that set the behaviour of the processing algorithms.
- enum [CalibrationData](#) {
 [Calibration_OffsetErrors](#),
 [Calibration_ApodizationSpectrum](#),
 [Calibration_ApodizationVector](#),
 [Calibration_Dispersion](#),
 [Calibration_Chirp](#),
 [Calibration_ExtendedAdjust](#),
 [Calibration_FixedPattern](#) }
Data describing the calibration of the processing routines.

- enum `ProcessingFlag` {
`Processing_UseOffsetErrors`,
`Processing_RemoveDCSpectrum`,
`Processing_RemoveAdvancedDCSpectrum`,
`Processing_UseApodization`,
`Processing_UseScanForApodization`,
`Processing_UseUndersamplingFilter`,
`Processing_UseDispersionCompensation`,
`Processing_UseDechirp`,
`Processing_UseExtendedAdjust`,
`Processing_FullRangeOutput`,
`Processing_FilterDC`,
`Processing_UseAutocorrCompensation`,
`Processing_UseDEFR`,
`Processing_OnlyWindowing`,
`Processing_RemoveFixedPattern` }

Flags that set the behaviour of the processing algorithms.

- enum `ProcessingAveragingAlgorithm` {
`Processing_Averaging_Min`,
`Processing_Averaging_Mean`,
`Processing_Averaging_Median`,
`Processing_Averaging_Norm2`,
`Processing_Averaging_Max`,
`Processing_Averaging_Fourier_Min`,
`Processing_Averaging_Fourier_Norm4`,
`Processing_Averaging_Fourier_Max`,
`Processing_Averaging_StandardDeviationAbs` }

This sets the averaging algorithm to be used for processing.

- enum `ApodizationWindowParameter` {
`ApodizationWindowParameter_Sigma`,
`ApodizationWindowParameter_Ratio`,
`ApodizationWindowParameter_Frequency` }

Sets certain parameters that are used by the window functions to be applied during apodization.

Functions

- `SPECTRALRADAR_API ProcessingHandle createProcessingForDevice (OCTDeviceHandle Dev)`
Creates suitable standard processing routines for the specified device (`OCTDeviceHandle`).
- `SPECTRALRADAR_API ProcessingHandle createProcessingForDeviceEx (OCTDeviceHandle Dev, int CameraIndex)`
Creates suitable standard processing routines for the specified device (`OCTDeviceHandle`) with camera index.
- `SPECTRALRADAR_API int getInputSize (ProcessingHandle Proc)`
Returns the expected input size (pixels per spectrum); of the processing algorithms.
- `SPECTRALRADAR_API int getAScanSize (ProcessingHandle Handle)`
gives the number of pixels in an A-Scan of the SpectralRadar device. This number is identical to the number of rows in a finished B-Scan.
- `SPECTRALRADAR_API void setApodizationWindow (ProcessingHandle Proc, ApodizationWindow Window)`
Sets the window function that is to be used for apodization. The selected function will be used in all subsequent processings.
- `SPECTRALRADAR_API int getApodizationWindow (ProcessingHandle Proc)`
Gets the window function that is being used for apodization.
- `SPECTRALRADAR_API void setApodizationWindowParameter (ProcessingHandle Proc, ApodizationWindowParameter Selection, double Value)`
Sets the apodization window parameter, such as window width or ratio between constant and cosine part.

- **SPECTRALRADAR_API** double `getApodizationWindowParameter` (`ProcessingHandle` Proc, `ApodizationWindowParameter` Selection)
Gets the apodization window parameter, such as window width or ratio between constant and cosine part.
- **SPECTRALRADAR_API** void `setDechirpAlgorithm` (`ProcessingHandle` Proc, `ProcessingType` Type)
Sets the algorithm that is to be sued for dechirping the input spectra.
- **SPECTRALRADAR_API** void `setProcessingParameterInt` (`ProcessingHandle` Proc, `ProcessingParameterInt` Selection, int Value)
Sets the specified integer value processing parameter.
- **SPECTRALRADAR_API** int `getProcessingParameterInt` (`ProcessingHandle` Proc, `ProcessingParameterInt` Selection)
Returns the specified integer value processing parameter.
- **SPECTRALRADAR_API** double `getProcessingParameterFloat` (`ProcessingHandle` Proc, `ProcessingParameterFloat` Selection)
Gets the specified processing paramter.
- **SPECTRALRADAR_API** void `setProcessingFlag` (`ProcessingHandle` Proc, `ProcessingFlag` Flag, `BOOL` Value)
Sets the specified processing flag.
- **SPECTRALRADAR_API** `BOOL` `getProcessingFlag` (`ProcessingHandle` Proc, `ProcessingFlag` Flag)
Returns TRUE if the specified processing flag is set, FALSE otherwise.
- **SPECTRALRADAR_API** void `setProcessingAveragingAlgorithm` (`ProcessingHandle` Proc, `ProcessingAveragingAlgorithm` Algorithm)
Sets the algorithm that is used for averaing by the processing.
- **SPECTRALRADAR_API** void `setCalibration` (`ProcessingHandle` Proc, `CalibrationData` Selection, `DataHandle` Data)
Sets the current active calibration data.
- **SPECTRALRADAR_API** void `getCalibration` (`ProcessingHandle` Proc, `CalibrationData` Selection, `DataHandle` Data)
Returns the currently active calibration parameter.
- **SPECTRALRADAR_API** void `measureCalibration` (`OCTDeviceHandle` Dev, `ProcessingHandle` Proc, `CalibrationData` Selection)
Measures the specified calibration parameters and uses them in subsequent processing.
- **SPECTRALRADAR_API** void `measureCalibrationEx` (`OCTDeviceHandle` Dev, `ProcessingHandle` Proc, `CalibrationData` Selection, int CameraIndex)
Measures the specified calibration parameters and uses them in subsequent processing with specified camera index.
- **SPECTRALRADAR_API** void `measureSpectrum` (`OCTDeviceHandle` Dev, `ProbeHandle` Probe, `ProcessingHandle` Proc, `BOOL` moveToApoPos)
Measures apodization spectrum and uses them in subsequent processing.
- **SPECTRALRADAR_API** void `saveCalibrationAuto` (`ProcessingHandle` Proc, `CalibrationData` Selection)
Saves the selected calibration in its default path.
- **SPECTRALRADAR_API** void `saveCalibration` (`ProcessingHandle` Proc, `CalibrationData` Selection, const char Path[])
Saves the selected calibration in the specified path.
- **SPECTRALRADAR_API** void `loadCalibration` (`ProcessingHandle` Proc, `CalibrationData` Selection, const char Path[])
Will load a specified calibration file and use for subsequent processing.
- **SPECTRALRADAR_API** void `setSpectrumOutput` (`ProcessingHandle` Proc, `DataHandle` Spectrum)
Sets the location for the resulting spectral data.
- **SPECTRALRADAR_API** void `setOffsetCorrectedSpectrumOutput` (`ProcessingHandle` Proc, `DataHandle` OffsetCorrectedSpectrum)
Sets the location for the resulting offset corrected spectral data.
- **SPECTRALRADAR_API** void `setDCCorrectedSpectrumOutput` (`ProcessingHandle` Proc, `DataHandle` ProcessingCorrectedSpectrum)
Sets the location for the resulting DC removed spectral data.

- **SPECTRALRADAR_API** void **setApodizedSpectrumOutput** (**ProcessingHandle** Proc, **DataHandle** Apodized-Spectrum)
Sets the location for the resulting apodized spectral data.
- **SPECTRALRADAR_API** void **setComplexDataOutput** (**ProcessingHandle** Proc, **ComplexDataHandle** ComplexBScan)
Sets the pointer the resulting complex B-Scan of the next processing is written to.
- **SPECTRALRADAR_API** void **setProcessedDataOutput** (**ProcessingHandle** Proc, **DataHandle** Scan)
Sets the pointer the resulting B-Scan of the next processing is written to.
- **SPECTRALRADAR_API** void **setHorMirroredDataOutput** (**ProcessingHandle** Proc, **DataHandle** Scan)
Sets the pointer the resulting B-Scan of the next processing is written to. The result will be written mirrored at the horizontal axis.
- **SPECTRALRADAR_API** void **setColoredDataOutput** (**ProcessingHandle** Proc, **ColoredDataHandle** BScan, **Coloring32BitHandle** Color)
Sets the pointer the resulting colored B-Scan of the next processing is written to.
- **SPECTRALRADAR_API** void **setTransposedColoredDataOutput** (**ProcessingHandle** Proc, **ColoredDataHandle** BScan, **Coloring32BitHandle** Color)
Sets the pointer the resulting colored B-Scan of the next processing is written to. The data will be transposed so that the first axis is the x-axis.
- **SPECTRALRADAR_API** void **executeProcessing** (**ProcessingHandle** Proc, **RawDataHandle** RawData)
Execute the processing.
- **SPECTRALRADAR_API** void **closeProcessing** (**ProcessingHandle** Proc)
Closes the processing and frees all temporary memory that was associated with it. Processing threads will be stopped.
- **SPECTRALRADAR_API** void **computeDispersion** (**DataHandle** Spectrum1In, **DataHandle** Spectrum2In, **DataHandle** ChirpOut, **DataHandle** DispOut)
Computes the dispersion and chirp of the two provided spectra, where both spectra need to have been subjected to same dispersion mismatch. Both spectra need to have been acquired for different path length differences.
- **SPECTRALRADAR_API** void **computeDispersionByCoeff** (double QuadraticIn, **DataHandle** ChirpIn, **DataHandle** DispOut)
Computes dispersion by a quadratic approximation specified by the quadratic factor.
- **SPECTRALRADAR_API** void **computeDispersionByImage** (**DataHandle** LinearKSpectralIn, **DataHandle** ChirpIn, **DataHandle** DispOut)
Guesses the dispersion based on the raw data specified. The raw data needs to be linearized in k before applying to this function.
- **SPECTRALRADAR_API** int **getNumberOfDispersionPresets** (**ProcessingHandle** Proc)
Gets the number of dispersion presets.
- **SPECTRALRADAR_API** const char * **getDispersionPresetName** (**ProcessingHandle** Proc, int Index)
Gets the name of the dispersion preset specified with index.
- **SPECTRALRADAR_API** void **setDispersionPresetByName** (**ProcessingHandle** Proc, const char *Name)
Sets the dispersion preset specified with name.
- **SPECTRALRADAR_API** void **setDispersionPresetByIndex** (**ProcessingHandle** Proc, int Index)
Sets the dispersion preset specified with index.
- **SPECTRALRADAR_API** void **setDispersionPresets** (**ProcessingHandle** Proc, **ProbeHandle** Probe)
Sets the dispersion presets for the probe.
- **SPECTRALRADAR_API** void **computeLinearKRawData** (**ComplexDataHandle** ComplexDataAfterFFT, **DataHandle** LinearKData)
Computes the linear k raw data of the complex data after FFT by an inverse Fourier transform.
- **SPECTRALRADAR_API** void **linearizeSpectralData** (**DataHandle** SpectralIn, **DataHandle** SpectraOut, **DataHandle** Chirp)
Linearizes the spectral data using the given chirp vector.

5.9.1 Detailed Description

Standard Processing Routines.

5.9.2 Typedef Documentation

5.9.2.1 ProcessingHandle

Handle for a processing routine.

5.9.3 Enumeration Type Documentation

5.9.3.1 enum ApodizationWindow

To select the apodization window function.

Enumerator

Apodization_Hann Hann window function.

Apodization_Hamming Hamming window function.

Apodization_Gauss Gaussian window function.

Apodization_TaperedCosine Tapered cosine window function.

Apodization_Blackman Blackman window function.

Apodization_BlackmanHarris 4-Term Blackman-Harris window function

Apodization_LightSourceBased The apodization function is determined, based on the shape of the light source at hand.

Warning

{This feature is still experimental.}

Apodization_Unknown Unknown apodization window.

5.9.3.2 enum ApodizationWindowParameter

Sets certain parameters that are used by the window functions to be applied during apodization.

Enumerator

ApodizationWindowParameter_Sigma Sets the width of a Gaussian apodization window.

ApodizationWindowParameter_Ratio Sets the ratio of the constant to the cosine part when using a tapered cosine window.

ApodizationWindowParameter_Frequency Sets the corner frequency of the filter applied when using a light-source based apodization.

Warning

{Light source based apodization is still experimental and might contain bugs or decrease performance of the OCT system.}

5.9.3.3 enum CalibrationData

Data describing the calibration of the processing routines.

Enumerator

Calibration_OffsetErrors Calibration vector used as offset.

Calibration_ApodizationSpectrum Calibration data used as reference spectrum.

Calibration_ApodizationVector Calibration data used as apodization multipliers.

Calibration_Dispersion Calibration data used to compensate for dispersion.

Calibration_Chirp Calibration data used for dechirping spectral data.

Calibration_ExtendedAdjust Calibration data used as extended adjust.

Calibration_FixedPattern Calibration data used as fixed scan pattern data.

5.9.3.4 enum ProcessingAveragingAlgorithm

This sets the averaging algorithm to be used for processing.

Warning

{This features is still experimental and might contain bugs.}

5.9.3.5 enum ProcessingFlag

Flags that set the behaviour of the processing algorithms.

Enumerator

Processing_UseOffsetErrors Flag identifying whether to apply offset error removal. This flag is activated by default.

Processing_RemoveDCSpectrum Flag sets whether the DC spectrum as measured is to be removed from the spectral data. This flag is activated by default.

Processing_RemoveAdvancedDCSpectrum Flag sets whether the DC spectrum to be removed is rescaled by the respective spectrum intensity it is applied to. This flag is activated by default.

Processing_UseApodization Flag identifying whether to apply apodization. This flag is activated by default.

Processing_UseScanForApodization Flag to determine whether the acquired data is to be averaged in order to compute an apodization spectrum. This flag is deactivated by default.

Processing_UseUndersamplingFilter Flag to activate or deactivate a filter removing undersampled signals from the A-scan. This flag is deactivated by default.

Processing_UseDispersionCompensation Flag activating or deactivating dispersion compensation. This flag is deactivated by default.

Processing_UseDechirp Flag identifying whether to apply dechirp. This flag is activated by default.

Processing_UseExtendedAdjust Flag identifying whether to use extended adjust. This flag is deactivated by default.

Processing_FullRangeOutput Flag identifying whether to use full range output. This flag is deactivated by default.

Processing_FilterDC Experimental: Flag for an experimental lateral DC filtering algorithm. This flag is deactivated by default.

Processing_UseAutocorrCompensation Flag activating or deactivating autocorrelation compensation. This flag is deactivated by default.

Processing_UseDEFR Exprtimental: Toggles dispersion encoded full range processing mode, eliminating folding of the signal at the top. This flag is deactivated by default.

Processing_OnlyWindowing Flag deactivating deconvolution in apodization processing, using windowing only. This flag is deactivated by default.

Processing_RemoveFixedPattern Flag for removal of fixed pattern noise, used for swept source OCT systems. This flag is deactivated by default.

5.9.3.6 enum ProcessingParameterFloat

Parameters that set the behaviour of the processing algorithms.

Enumerator

Processing_ApodizationDamping Sets how much influence newly acquired apodizations have compared to older ones.

Processing_MinElectrons Determines the minimum signal intensity on the edge channels of the spectra.

Warning

{Setting this value may seriously reduce performance of the system.}

5.9.3.7 enum **ProcessingParameterInt**

Parameters that set the behaviour of the processing algorithms.

Enumerator

Processing_SpectrumAveraging Identifier for averaging of several subsequent spectra prior to Fourier transform.

Processing_AScanAveraging Identifier for averaging the absolute values of several subsequent A-scan after Fourier transform.

Processing_BScanAveraging Averaging of subsequent B-scans.

Processing_ZeroPadding Identifier for zero padding prior to Fourier transformation.

Processing_NumberOfThreads The maximum number of threads to be used by processing. A value of 0 indicates automatic selection, equal to the number of cores in the host PC.

Processing_FourierAveraging Averaging of fourier spectra.

5.9.3.8 enum **ProcessingType**

defines the algorithm used for dechirping the input signal and Fourier transformation

Enumerator

Processing_StandardFFT FFT with no dechirp algorithm applied.

Processing_StandardNDFT Full matrix multiplication ("filter bank"). Mathematical precise dechirp, but rather slow.

Processing_iFFT1 Linear interpolation prior to FFT.

Processing_iFFT2 Linear interpolation with 2x oversampling prior to FFT.

Processing_iFFT3 Linear interpolation with 3x oversampling prior to FFT.

Processing_iFFT4 Linear interpolation with 4x oversampling prior to FFT.

Processing_NFFT1 NFFT algorithm with parameter m=1.

Processing_NFFT2 NFFT algorithm with parameter m=2.

Processing_NFFT3 NFFT algorithm with parameter m=3.

Processing_NFFT4 NFFT algorithm with parameter m=4.

5.9.4 Function Documentation

5.9.4.1 void **closeProcessing** (**ProcessingHandle Proc**)

Closes the processing and frees all temporary memory that was associated with it. Processing threads will be stopped.

5.9.4.2 void **computeDispersion** (**DataHandle Spectrum1In**, **DataHandle Spectrum2In**, **DataHandle ChirpOut**, **DataHandle DispOut**)

Computes the dispersion and chirp of the two provided spectra, where both spectra need to have been subjected to same dispersion mismatch. Both spectra need to have been acquired for different path length differences.

5.9.4.3 void **computeDispersionByCoeff** (**double QuadraticIn**, **DataHandle ChirpIn**, **DataHandle DispOut**)

Computes dispersion by a quadratic approximation specified by the quadratic factor.

5.9.4.4 void **computeDispersionByImage** (**DataHandle LinearKSpectralIn**, **DataHandle ChirpIn**, **DataHandle DispOut**)

Guesses the dispersion based on the raw data specified. The raw data needs to be linearized in k before applying to this function.

5.9.4.5 void computeLinearKRawData (ComplexDataHandle *ComplexDataAfterFFT*, DataHandle *LinearKData*)

Computes the linear k raw data of the complex data after FFT by an inverse Fourier transform.

5.9.4.6 ProcessingHandle createProcessingForDevice (OCTDeviceHandle *Dev*)

Creates suitable standard processing routines for the specified device ([OCTDeviceHandle](#)).

5.9.4.7 ProcessingHandle createProcessingForDeviceEx (OCTDeviceHandle *Dev*, int *CameraIndex*)

Creates suitable standard processing routines for the specified device ([OCTDeviceHandle](#)) with camera index.

5.9.4.8 void executeProcessing (ProcessingHandle *Proc*, RawDataHandle *RawData*)

Execute the processing.

The specified raw data will be transformed. Results will be written to data objects specified by [setProcessedDataOutput\(\)](#), [setComplexDataOutput\(\)](#), [setColoredDataOutput\(\)](#), etc.

5.9.4.9 int getApodizationWindow (ProcessingHandle *Proc*)

Gets the window function that is being used for apodization.

Parameters

<i>Proc</i>	handle to the OCTDeviceHandle that was initially provided by <code>initDevice</code> .
-------------	--

5.9.4.10 double getApodizationWindowParameter (ProcessingHandle *Proc*, ApodizationWindowParameter *Selection*)

Gets the apodization window parameter, such as window width or ratio between constant and cosine part.

5.9.4.11 int getAScanSize (ProcessingHandle *Proc*)

gives the number of pixels in an A-Scan of the SpectralRadar device. This number is identical to the number of rows in a finished B-Scan.

Parameters

<i>Proc</i>	Processing that is used to get the A-Scan.
-------------	--

Returns

The number of pixels in an A-Scan of the SpectralRadar device.

5.9.4.12 void getCalibration (ProcessingHandle *Proc*, CalibrationData *Selection*, DataHandle *Data*)

Returns the currently active calibration parameter.

5.9.4.13 const char * getDispersionPresetName (ProcessingHandle *Proc*, int *Index*)

Gets the name of the dispersion preset specified with index.

5.9.4.14 int getInputSize (ProcessingHandle *Proc*)

Returns the expected input size (pixels per spectrum); of the processing algorithms.

This function is provided for convenience as processing routines can be used independently of the device.

5.9.4.15 int getNumberOfDispersionPresets (ProcessingHandle *Proc*)

Gets the number of dispersion presets.

5.9.4.16 `BOOL getProcessingFlag (ProcessingHandle Proc, ProcessingFlag Flag)`

Returns TRUE if the specified processing flag is set, FALSE otherwise.

5.9.4.17 `double getProcessingParameterFloat (ProcessingHandle Proc, ProcessingParameterFloat Selection)`

Gets the specified processing parameter.

5.9.4.18 `int getProcessingParameterInt (ProcessingHandle Proc, ProcessingParameterInt Selection)`

Returns the specified integer value processing parameter.

5.9.4.19 `void linearizeSpectralData (DataHandle SpectralIn, DataHandle SpectraOut, DataHandle Chirp)`

Linearizes the spectral data using the given chirp vector.

5.9.4.20 `void loadCalibration (ProcessingHandle Proc, CalibrationData Selection, const char Path[])`

Will load a specified calibration file and use for subsequent processing.

5.9.4.21 `void measureCalibration (OCTDeviceHandle Dev, ProcessingHandle Proc, CalibrationData Selection)`

Measures the specified calibration parameters and uses them in subsequent processing.

5.9.4.22 `void void measureCalibrationEx (OCTDeviceHandle Dev, ProcessingHandle Proc, CalibrationData Selection, int CameraIndex)`

Measures the specified calibration parameters and uses them in subsequent processing with specified camera index.

5.9.4.23 `void measureSpectrum (OCTDeviceHandle Dev, ProbeHandle Probe, ProcessingHandle Proc, BOOL moveToApoPos)`

Measures apodization spectrum and uses them in subsequent processing.

5.9.4.24 `void saveCalibration (ProcessingHandle Proc, CalibrationData Selection, const char Path[])`

Saves the selected calibration in the specified path.

Warning

This will override your default calibration of the device if you specify the default path.

5.9.4.25 `void saveCalibrationAuto (ProcessingHandle Proc, CalibrationData Selection)`

Saves the selected calibration in its default path.

Warning

This will override your default calibration of the device.

5.9.4.26 `void setApodizationWindow (ProcessingHandle Proc, ApodizationWindow Window)`

Sets the window function that is to be used for apodization. The selected function will be used in all subsequent processings.

If this function is not explicitly called a Hann window will be used.

Parameters

<i>Handle</i>	Processing handle.
<i>Window</i>	The apodization window that is used for data processing.

5.9.4.27 void setApodizationWindowParameter (**ProcessingHandle Proc**, **ApodizationWindowParameter Selection**, double *Value*)

Sets the apodization window parameter, such as window width or ratio between constant and cosine part.

5.9.4.28 void setApodizedSpectrumOutput (**ProcessingHandle Proc**, **DataHandle ApodizedSpectrum**)

Sets the location for the resulting apodized spectral data.

5.9.4.29 void setCalibration (**ProcessingHandle Proc**, **CalibrationData Selection**, **DataHandle Data**)

Sets the current active calibration data.

5.9.4.30 void setColoredDataOutput (**ProcessingHandle Proc**, **ColoredDataHandle BScan**, **Coloring32BitHandle Color**)

Sets the pointer the resulting colored B-Scan of the next processing is written to.

5.9.4.31 void setComplexDataOutput (**ProcessingHandle Proc**, **ComplexDataHandle ComplexBScan**)

Sets the pointer the resulting complex B-Scan of the next processing is written to.

If set to 0 no complex data result will be created in the next processing.

5.9.4.32 void setDCCorrectedSpectrumOutput (**ProcessingHandle Proc**, **DataHandle ProcessingCorrectedSpectrum**)

Sets the location for the resulting DC removed spectral data.

5.9.4.33 void setDechirpAlgorithm (**ProcessingHandle Proc**, **ProcessingType Type**)

Sets the algorithm that is to be used for dechirping the input spectra.

5.9.4.34 void setDispersionPresetByIndex (**ProcessingHandle Proc**, int *Index*)

Sets the dispersion preset specified with index.

5.9.4.35 void setDispersionPresetByName (**ProcessingHandle Proc**, const char * *Name*)

Sets the dispersion preset specified with name.

5.9.4.36 void setDispersionPresets (**ProcessingHandle Proc**, **ProbeHandle Probe**)

Sets the dispersion presets for the probe.

5.9.4.37 void setHorMirroredDataOutput (**ProcessingHandle Proc**, **DataHandle Scan**)

Sets the pointer the resulting B-Scan of the next processing is written to. The result will be written mirrored at the horizontal axis.

If set to 0 no floating point processed data in dB will be created in the next processing.

5.9.4.38 void setOffsetCorrectedSpectrumOutput (**ProcessingHandle Proc**, **DataHandle OffsetCorrectedSpectrum**)

Sets the location for the resulting offset corrected spectral data.

5.9.4.39 void setProcessedDataOutput (ProcessingHandle *Proc*, DataHandle *Scan*)

Sets the pointer the resulting B-Scan of the next processing is written to.

If set to 0 no processed floating point data in dB will be created in the next processing.

5.9.4.40 void setProcessingAveragingAlgorithm (ProcessingHandle *Proc*, ProcessingAveragingAlgorithm *Algorithm*)

Sets the algorithm that is used for averaging by the processing.

5.9.4.41 void setProcessingFlag (ProcessingHandle *Proc*, ProcessingFlag *Flag*, BOOL *Value*)

Sets the specified processing flag.

5.9.4.42 setProcessingParameterInt (ProcessingHandle *Proc*, ProcessingParameterInt *Selection*, int *Value*)

Sets the specified integer value processing parameter.

5.9.4.43 void setSpectrumOutput (ProcessingHandle *Proc*, DataHandle *Spectrum*)

Sets the location for the resulting spectral data.

5.9.4.44 void setTransposedColoredDataOutput (ProcessingHandle *Proc*, ColoredDataHandle *BScan*, Coloring32BitHandle *Color*)

Sets the pointer the resulting colored B-Scan of the next processing is written to. The data will be transposed so that the first axis is the x-axis.

5.10 Export and Import

Export functionality to store data to disk and load it from there.

Enumerations

- enum [Data1DExportFormat](#) {
[Data1DExport_RAW](#),
[Data1DExport_TXT](#),
[Data1DExport_CSV](#),
[Data1DExport_TableTXT](#),
[Data1DExport_Fits](#) }
Export format for 1-dimensional data ([DataHandle](#)).
- enum [Data2DExportFormat](#) {
[Data2DExport_SRM](#),
[Data2DExport_RAW](#),
[Data2DExport_TXT](#),
[Data2DExport_CSV](#),
[Data2DExport_TableTXT](#),
[Data2DExport_Fits](#) }
Export format for 2-dimensional data ([DataHandle](#)).
- enum [Data3DExportFormat](#) {
[Data3DExport_SRM](#),
[Data3DExport_RAW](#),
[Data3DExport_TXT](#),
[Data3DExport_CSV](#),
[Data3DExport_VFF](#),
[Data3DExport_VTK](#),
[Data3DExport_Fits](#),
[Data3DExport_TIFF](#) }
Export format for 3-dimensional data ([DataHandle](#)).
- enum [ComplexDataExportFormat](#) { [ComplexDataExport_RAW](#) }
Export format for complex data.
- enum [ColoredDataExportFormat](#) {
[ColoredDataExport_SRM](#),
[ColoredDataExport_RAW](#),
[ColoredDataExport_BMP](#),
[ColoredDataExport_PNG](#),
[ColoredDataExport_JPG](#),
[ColoredDataExport_PDF](#),
[ColoredDataExport_TIFF](#) }
Export format for images ([ColoredDataHandle](#)).
- enum [DataImportFormat](#) { [DataImport_SRM](#) }
Supported import format to load data from disk.
- enum [RawDataExportFormat](#) {
[RawDataExport_RAW](#),
[RawDataExport_SRR](#) }
Supported raw data export formats to store data to disk.
- enum [RawDataImportFormat](#) { [RawDataImport_SRR](#) }
Supported raw data import formats to load data from disk.

Functions

- [SPECTRALRADAR_API](#) void [exportData1D](#) ([DataHandle](#) Data, [Data1DExportFormat](#) Format, const char *Path)

Exports 1-dimensional data ([DataHandle](#)).

- [SPECTRALRADAR_API](#) void [exportData2D](#) ([DataHandle](#) Data, [Data2DExportFormat](#) Format, const char *Path)

Exports 2-dimensional data ([DataHandle](#)).

- [SPECTRALRADAR_API](#) void [exportData3D](#) ([DataHandle](#) Volume, [Data3DExportFormat](#) Format, const char *Path)

Exports 3-dimensional data ([DataHandle](#)).

- [SPECTRALRADAR_API](#) void [exportComplexData](#) ([ComplexDataHandle](#), [ComplexDataExportFormat](#), const char *)

Exports 1-, 2- and 3-dimensional complex data ([ComplexDataHandle](#)).

- [SPECTRALRADAR_API](#) void [exportColoredData](#) ([ColoredDataHandle](#) Image, [ColoredDataExportFormat](#) Format, const char *fileName)

Exports colored data ([ColoredDataHandle](#)).

- [SPECTRALRADAR_API](#) void [importColoredData](#) ([ColoredDataHandle](#) ColoredData, [DataImportFormat](#) Format, const char *Path)

Imports colored data ([ColoredDataHandle](#)) with the specified format and copied it into a data object ([ColoredDataHandle](#)).

- [SPECTRALRADAR_API](#) void [importData](#) ([DataHandle](#) Data, [DataImportFormat](#) Format, const char *Path)

Imports data with the specified format and copies it into a data object ([DataHandle](#)).

- [SPECTRALRADAR_API](#) void [exportRawData](#) ([RawDataHandle](#) Raw, [RawDataExportFormat](#) Format, const char *Path)

Exports the specified data to disk.

- [SPECTRALRADAR_API](#) void [importRawData](#) ([RawDataHandle](#) Raw, [RawDataImportFormat](#) Format, const char *Path)

Imports the specified data from disk.

5.10.1 Detailed Description

Export functionality to store data to disk and load it from there.

5.10.2 Enumeration Type Documentation

5.10.2.1 enum [ColoredDataExportFormat](#)

Export format for images ([ColoredDataHandle](#)).

Enumerator

[ColoredDataExport_SRM](#) Spectral Radar Metaformat, containing no data but all additional parameters, such as spacing, size, etc.

[ColoredDataExport_RAW](#) RAW data format containing the data of the object as binary, 32-bit unsigned integer values, little endian. The concrete format of the data depends on the colored data object ([ColoredDataHandle](#)). In most cases it will be RGB32 or RGBA32.

[ColoredDataExport_BMP](#) BMP - Bitmap image format.

[ColoredDataExport_PNG](#) PNG image format.

[ColoredDataExport_JPG](#) JPG/JPEG image format.

[ColoredDataExport_PDF](#) PDF image format.

[ColoredDataExport_TIFF](#) TIFF image format.

5.10.2.2 enum **ComplexDataExportFormat**

Export format for complex data.

Enumerator

ComplexDataExport_RAW RAW data format containing binary data.

5.10.2.3 enum **Data1DExportFormat**

Export format for 1-dimensional data ([DataHandle](#)).

Enumerator

Data1DExport_RAW RAW data format containing the data of the object as binary, single precision floating point values, little endian.

Data1DExport_TXT TXT is a text file having all values stored space separated and human readable.

Data1DExport_CSV CSV (Comma Separated Values) is a text file having all values stored, comma separated and human readable.

Data1DExport_TableTXT TableTXT is a human readable text-file in a table like format, having the physical 1- and 2-axis as first two columns and the data value as third.

Data1DExport_Fits FITS Data format.

5.10.2.4 enum **Data2DExportFormat**

Export format for 2-dimensional data ([DataHandle](#)).

Enumerator

Data2DExport_SRM Spectral Radar Metaformat, containing no data but all additional parameters, such as spacing, size, etc.

Data2DExport_RAW RAW data format containing the data of the object as binary, single precision floating point values, little endian.

Data2DExport_TXT TXT is a text file having all values stored space separated and human readable.

Data2DExport_CSV CSV (Comma Separated Values) is a text file having all values stored, comma separated and human readable.

Data2DExport_TableTXT TableTXT is a human readable text-file in a table like format, having the physical 1- and 2-axis as first two columns and the data value as third.

Data2DExport_Fits FITS Data format.

5.10.2.5 enum **Data3DExportFormat**

Export format for 3-dimensional data ([DataHandle](#)).

Enumerator

Data3DExport_SRM Spectral Radar Metaformat, containing no data but all additional parameters, such as spacing, size, etc.

Data3DExport_RAW RAW data format containing the data of the object as binary, single precision floating point values, little endian.

Data3DExport_TXT TXT is a text file having all values stored space separated and human readable.

Data3DExport_CSV CSV (Comma Separated Values) is a text file having all values stored, comma separated and human readable.

Data3DExport_VFF VFF data format.

Data3DExport_VTK VTK data format.

Data3DExport_Fits FITS Data format.

Data3DExport_TIFF TIFF Data format.

5.10.2.6 enum **DataImportFormat**

Supported import format to load data from disk.

Enumerator

DataImport_SRM Spectral Radar Metaformat, containing no data but all additional parameters, such as spacing, size, etc. It is searched for an appropriate file with same name but different extension containing the according data.

5.10.2.7 enum **RawDataExportFormat**

Supported raw data export formats to store data to disk.

Enumerator

RawDataExport_RAW Single precision floating point raw data.

RawDataExport_SRR Spectral Radar raw data format, specified additional information such as apodization scans, scan range, etc.

5.10.2.8 enum **RawDataImportFormat**

Supported raw data import formats to load data from disk.

Enumerator

RawDataImport_SRR Spectral Radar raw data format, specified additional information such as apodization scans, scan range, etc.

5.10.3 Function Documentation

5.10.3.1 void exportColoredData (ColoredDataHandle *Image*, ColoredDataExportFormat *Format*, const char * *fileName*)

Exports colored data ([ColoredDataHandle](#))

5.10.3.2 void exportComplexData (ComplexDataHandle , ComplexDataExportFormat , const char *)

Exports 1-, 2- and 3-dimensional complex data ([ComplexDataHandle](#))

5.10.3.3 void exportData1D (DataHandle *Data*, Data1DExportFormat *Format*, const char * *Path*)

Exports 1-dimensional data ([DataHandle](#)).

5.10.3.4 void exportData2D (DataHandle *Data*, Data2DExportFormat *Format*, const char * *Path*)

Exports 2-dimensional data ([DataHandle](#)).

5.10.3.5 void exportData3D (DataHandle *Volume*, Data3DExportFormat *Format*, const char * *Path*)

Exports 3-dimensional data ([DataHandle](#)).

5.10.3.6 void exportRawData (RawDataHandle *Raw*, RawDataExportFormat *Format*, const char * *Path*)

Exports the specified data to disk.

5.10.3.7 void importColoredData (ColoredDataHandle *ColoredData*, DataImportFormat *Format*, const char * *Path*)

Imports colored data ([ColoredDataHandle](#)) with the specified format and copied it into a data object ([ColoredDataHandle](#))

5.10.3.8 `void importData (DataHandle Data, DataImportFormat Format, const char * Path)`

Imports data with the specified format and copies it into a data object ([DataHandle](#)).

5.10.3.9 `void importRawData (RawDataHandle Raw, RawDataImportFormat Format, const char * Path)`

Imports the specified data from disk.

5.11 Volume

Functionality to store and access volume data.

Enumerations

- enum [Direction](#) {
[Direction_1](#),
[Direction_2](#),
[Direction_3](#) }
Specifies a direction.
- enum [Plane2D](#) {
[Plane2D_12](#),
[Plane2D_23](#),
[Plane2D_13](#) }
Planes for slices of the volume data.

Functions

- [SPECTRALRADAR_API](#) void [appendRawData](#) ([RawDataHandle](#) Data, [RawDataHandle](#) DataToAppend, [Direction](#) direction)
Appends the new raw data to the old raw data in the specified direction.
- [SPECTRALRADAR_API](#) void [getRawDataSliceIndex](#) ([RawDataHandle](#) Data, [RawDataHandle](#) Slice, [Direction](#) SliceNormalDirection, int Index)
Returns a slice of raw data in the specified direction at the specified index.
- [SPECTRALRADAR_API](#) double [analyzeData](#) ([DataHandle](#) Data, [DataAnalyzation](#) Selection)
Performs the selected analyzation of the specified data and returns the resulting value.
- [SPECTRALRADAR_API](#) double [analyzeAScan](#) ([DataHandle](#) Data, [AScanAnalyzation](#) Selection)
Performs the selected analyzation of the specified A-scan and returns the resulting value.
- [SPECTRALRADAR_API](#) void [determineDynamicRange](#) ([DataHandle](#) Data, float *MinRange_dB, float *MaxRange_dB)
Gives a rough estimation of the dynamic range of the specified data object.
- [SPECTRALRADAR_API](#) void [transpose](#) ([DataHandle](#) DataIn, [DataHandle](#) DataOut)
Transposes the given data and writes the result to DataOut.
- [SPECTRALRADAR_API](#) void [transposeInplace](#) ([DataHandle](#) Data)
Transposes the given Data.
- [SPECTRALRADAR_API](#) void [transposeAndScaleData](#) ([DataHandle](#) DataIn, [DataHandle](#) DataOut, float Min, float Max)
Transposes the given data and scales it to the range [Min, Max].
- [SPECTRALRADAR_API](#) void [normalizeData](#) ([DataHandle](#) Data, float Min, float Max)
Scales the given data to the range [Min, Max].
- [SPECTRALRADAR_API](#) void [lockData](#) ([DataHandle](#) Data)
Locks the given data.
- [SPECTRALRADAR_API](#) void [unlockData](#) ([DataHandle](#) Data)
Unlocks the given data.
- [SPECTRALRADAR_API](#) void [getDataSlicePos](#) ([DataHandle](#) Data, [DataHandle](#) Slice, [Direction](#) SliceNormalDirection, double Pos)
Returns a slice of data in the specified direction at the specified position.
- [SPECTRALRADAR_API](#) void [getDataSliceIndex](#) ([DataHandle](#) Data, [DataHandle](#) Slice, [Direction](#) SliceNormalDirection, int Index)
Returns a slice of data in the specified direction at the specified index.
- [SPECTRALRADAR_API](#) void [getDataSliceAnalyzed](#) ([DataHandle](#) Data, [DataHandle](#) Slice, [Direction](#) SliceNormalDirection, [DataAnalyzation](#) Selection)

- Returns a slice of data that has been computed of all slice using the specified analyzation method.*
- **SPECTRALRADAR_API** void **appendData** (**DataHandle** Data, **DataHandle** NewData, **Direction** direction)
Appends the new data to the old data in the specified direction.
 - **SPECTRALRADAR_API** void **cropData** (**DataHandle** Data, **Direction** direction, int Index)
Crops the data at the specific direction at the given index. The result will contain the data with range [0, index] at the cropping direction.
 - **SPECTRALRADAR_API** void **cropDataEx** (**DataHandle** Data, **Direction** direction, int IndexMax, int IndexMin)
Crops the data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping direction.
 - **SPECTRALRADAR_API** void **separateData** (**DataHandle** Data1, **DataHandle** Data2, int SeparationIndex, **Direction** Dir)
Separates the data at the given index at specific separation direction. The first part of the separated data will be in Data1, the second separated in Data2.
 - **SPECTRALRADAR_API** void **flipData** (**DataHandle** Data, **Direction** FlippingDirection)
Flips the data around the specific direction.
 - **SPECTRALRADAR_API** void **getComplexDataSlicePos** (**ComplexDataHandle** Data, **ComplexDataHandle** Slice, **Direction** SliceNormalDirection, double Pos)
Returns a slice of data in the specified direction at the specified position.
 - **SPECTRALRADAR_API** void **getComplexDataSliceIndex** (**ComplexDataHandle** Data, **ComplexDataHandle** Slice, **Direction** SliceNormalDirection, int Index)
Returns a slice of data in the specified direction at the specified index.
 - **SPECTRALRADAR_API** void **appendComplexData** (**ComplexDataHandle** Data, **ComplexDataHandle** DataToAppend, **Direction** direction)
Appends the new data to the old data in the specified direction.
 - **SPECTRALRADAR_API** void **cropComplexData** (**ComplexDataHandle** Data, **Direction** CroppingDirection, int IndexMax, int IndexMin)
Crops the complex data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping direction.
 - **SPECTRALRADAR_API** void **cropColoredData** (**ColoredDataHandle** Data, **Direction** CroppingDirection, int IndexMax, int IndexMin)
Crops the colored data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping direction.
 - **SPECTRALRADAR_API** void **appendColoredData** (**ColoredDataHandle** Data, **ColoredDataHandle** DataToAppend, **Direction** AppendingDirection)
Appends the new colored data to the old colored data in the specified direction.
 - **SPECTRALRADAR_API** void **getColoredDataSlicePos** (**ColoredDataHandle** Data, **ColoredDataHandle** Slice, **Direction** SliceNormalDirection, double Pos)
Get a slice of the colored data with specific slicing direction at given index position.
 - **SPECTRALRADAR_API** void **getColoredDataSliceIndex** (**ColoredDataHandle** Data, **ColoredDataHandle** Slice, **Direction** SliceNormalDirection, int Index)
Get a slice of the colored data with specific slicing direction at given index.
 - **SPECTRALRADAR_API** **ImageFieldHandle** **createImageField** (void)
Creates an object holding image field data.
 - **SPECTRALRADAR_API** void **clearImageField** (**ImageFieldHandle** ImageField)
Frees an object holding image field data.
 - **SPECTRALRADAR_API** void **saveImageField** (**ImageFieldHandle** ImageField, const char *Path)
Saves data containing image field data.
 - **SPECTRALRADAR_API** void **loadImageField** (**ImageFieldHandle** ImageField, const char *Path)
Loads data containing image field data.
 - **SPECTRALRADAR_API** void **determineImageField** (**ImageFieldHandle** ImageField, **DataHandle** Surface)
Determines the image field correction of the surface.
 - **SPECTRALRADAR_API** void **determineImageFieldForProbe** (**ProbeHandle** Probe, **DataHandle** Surface)
Determines the image field correction of the surface for the specified probe handle.

- **SPECTRALRADAR_API** void **determineImageFieldForProbeWithMap** (ProbeHandle Probe, DataHandle Surface, DataHandle Map)
Determines the image field correction of the surface for the specified probe handle using the given map. Values != 0 in the map specifies to use the data in the surface handle otherwise they will be interpolated.
- **SPECTRALRADAR_API** void **correctImageField** (ImageFieldHandle ImageField, ScanPatternHandle Pattern, DataHandle Data)
Applies the image field correction to the B-Scan or volume data .
- **SPECTRALRADAR_API** void **correctSurface** (ImageFieldHandle ImageField, DataHandle Surface)
Applies the image field correction to the given Surface.
- **SPECTRALRADAR_API** void **correctImageFieldFromProbe** (ProbeHandle Probe, ScanPatternHandle Pattern, DataHandle Data)
Applies the image field correction saved in the probe handle to the B-Scan or volume data .

5.11.1 Detailed Description

Functionality to store and access volume data.

5.11.2 Enumeration Type Documentation

5.11.2.1 enum Direction

Specifies a direction.

Enumerator

Direction_1 The 1-axis direction.

Direction_2 The 2-axis direction.

Direction_3 The 3-axis direction.

5.11.2.2 enum Plane2D

Planes for slices of the volume data.

Enumerator

Plane2D_12 The 12 (XZ) plane, orthogonal to the 3 (Y) axis.

Plane2D_23 The 23 (XY) plane, orthogonal to the 3 (Z) axis.

Plane2D_13 The 13 (ZY) plane, orthogonal to the 2 (X) axis.

5.11.3 Function Documentation

5.11.3.1 double analyzeAScan (DataHandle Data, AScanAnalysis Selection)

Performs the selected analyzation of the specified A-scan and returns the resulting value.

5.11.3.2 double analyzeData (DataHandle Data, DataAnalysis Selection)

Performs the selected analyzation of the specified data and returns the resulting value.

5.11.3.3 void appendColoredData (ColoredDataHandle Data, ColoredDataHandle DataToAppend, Direction AppendingDirection)

Appends the new colored data to the old colored data in the specified direction.

5.11.3.4 void appendComplexData (ComplexDataHandle Data, ComplexDataHandle DataToAppend, Direction direction)

Appends the new data to the old data in the specified direction.

5.11.3.5 void appendData (DataHandle Data, DataHandle NewData, Direction direction)

Appends the new data to the old data in the specified direction.

5.11.3.6 void appendRawData (RawDataHandle Data, RawDataHandle DataToAppend, Direction direction)

Appends the new raw data to the old raw data in the specified direction.

5.11.3.7 void clearImageField (ImageFieldHandle ImageField)

Frees an object holding image field data.

5.11.3.8 void correctImageField (ImageFieldHandle ImageField, ScanPatternHandle Pattern, DataHandle Data)

Applies the image field correction to the B-Scan or volume data .

5.11.3.9 void correctImageFieldFromProbe (ProbeHandle Probe, ScanPatternHandle Pattern, DataHandle Data)

Applies the image field correction saved in the probe handle to the B-Scan or volume data .

5.11.3.10 void correctSurface (ImageFieldHandle ImageField, DataHandle Surface)

Applies the image field correction to the given Surface.

5.11.3.11 ImageFieldHandle creatImageField (void)

Creates an object holding image field data.

5.11.3.12 void cropColoredData (ColoredDataHandle Data, Direction CroppingDirection, int IndexMax, int IndexMin)

Crops the colored data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping dircition.

5.11.3.13 void cropComplexData (ComplexDataHandle Data, Direction CroppingDirection, int IndexMax, int IndexMin)

Crops the complex data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping dircition.

5.11.3.14 void cropData (DataHandle Data, Direction direction, int Index)

Crops the data at the specific direction at the given index. The result will contain the data with range [0, index] at the cropping dircition.

5.11.3.15 void cropDataEx (DataHandle Data, Direction direction, int IndexMax, int IndexMin)

Crops the data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping dircition.

5.11.3.16 void determineDynamicRange (DataHandle Data, float * MinRange_dB, float * MaxRange_dB)

Gives a rough estimation of the dynamic range of the specified data object.

This functions assumes that the data contains an A-scan and performs A-scan specific analysis on it.

5.11.3.17 void determinImageField (ImageFieldHandle ImageField, DataHandle Surface)

Determines the image field correction of the surface.

5.11.3.18 **void determinelImageFieldForProbe (*ProbeHandle Probe*, *DataHandle Surface*)**

Determines the image field correction of the surface for the specified probe handle.

5.11.3.19 **void determinelImageFieldForProbeWithMap (*ProbeHandle Probe*, *DataHandle Surface*, *DataHandle Map*)**

Determines the image field correction of the surface for the specified probe handle using the given map. Values != 0 in the map specifies to use the data in the surface handle otherwise they will be interpolated.

5.11.3.20 **void flipData (*DataHandle Data*, *Direction FlippingDirection*)**

Flips the data around the specific direction.

5.11.3.21 **void getColoredDataSliceIndex (*ColoredDataHandle Data*, *ColoredDataHandle Slice*, *Direction SliceNormalDirection*, int *Index*)**

Get a slice of the colored data with specific slicing direction at given index.

5.11.3.22 **void getColoredDataSlicePos (*ColoredDataHandle Data*, *ColoredDataHandle Slice*, *Direction SliceNormalDirection*, double *Pos*)**

Get a slice of the colored data with specific slicing direction at given index position.

5.11.3.23 **void getComplexDataSliceIndex (*ComplexDataHandle Data*, *ComplexDataHandle Slice*, *Direction SliceNormalDirection*, int *Index*)**

Returns a slice of data in the specified direction at the specified index.

5.11.3.24 **void getComplexDataSlicePos (*ComplexDataHandle Data*, *ComplexDataHandle Slice*, *Direction SliceNormalDirection*, double *Pos*)**

Returns a slice of data in the specified direction at the specified position.

5.11.3.25 **void getDataSliceAnalyzed (*DataHandle Data*, *DataHandle Slice*, *Direction SliceNormalDirection*, *DataAnalyzation Selection*)**

Returns a slice of data that has been computed of all slice using the specified analyzation method.

5.11.3.26 **void getDataSliceIndex (*DataHandle Data*, *DataHandle Slice*, *Direction SliceNormalDirection*, int *Index*)**

Returns a slice of data in the specified direction at the specified index.

5.11.3.27 **void getDataSlicePos (*DataHandle Data*, *DataHandle Slice*, *Direction SliceNormalDirection*, double *Pos*)**

Returns a slice of data in the specified direction at the specified position.

5.11.3.28 **void getRawDataSliceIndex (*RawDataHandle Data*, *RawDataHandle Slice*, *Direction SliceNormalDirection*, int *Index*)**

Returns a slice of raw data in the specified direction at the specified index.

5.11.3.29 **void loadImageField (*ImageFieldHandle ImageField*, const char * *Path*)**

Loads data containing image field data.

5.11.3.30 **void lockData (*DataHandle Data*)**

Locks the given data.

5.11.3.31 **void normalizeData (*DataHandle Data*, float *Min*, float *Max*)**

Scales the given data to the range [Min, Max].

5.11.3.32 **void saveImageField (ImageFieldHandle *ImageField*, const char * *Path*)**

Saves data containing image field data.

5.11.3.33 **void separateData (DataHandle *Data1*, DataHandle *Data2*, int *SeparationIndex*, Direction *Dir*)**

Separates the data at the given index at specific separation direction. The first part of the separated data will be in Data1, the second separated in Data2.

5.11.3.34 **void transpose (DataHandle *DataIn*, DataHandle *DataOut*)**

Transposes the given data and writes the result to DataOut.

5.11.3.35 **void transposeAndScaleData (DataHandle *DataIn*, DataHandle *DataOut*, float *Min*, float *Max*)**

Transposes the given data and scales it to the range [Min, Max].

5.11.3.36 **void transposeInplace (DataHandle *Data*)**

Transposes the given Data.

5.11.3.37 **void unlockData (DataHandle *Data*)**

Unlocks the given data.

5.12 ProbeCalibration

Functionality to perform the probe calibration. Please use the ThorImageOCT software to perform the probe calibration if necessary.

Functionality to perform the probe calibration. Please use the ThorImageOCT software to perform the probe calibration if necessary.

5.13 Doppler

Doppler Processing Routines.

Typedefs

- typedef struct
C_DopplerProcessing * [DopplerProcessingHandle](#)
Handle used for Doppler processing.

Enumerations

- enum [DopplerPropertyInt](#) {
[DopplerAveraging_1](#),
[DopplerAveraging_2](#),
[DopplerStride_1](#),
[DopplerStride_2](#) }
Values that determine the behaviour of the Doppler processing routines.
- enum [DopplerPropertyFloat](#) { [DopplerRefractiveIndex](#) }
Values that determine the behaviour of the Doppler processing routines.
- enum [DopplerFlag](#) { [DopplerVelocityScaling](#) }
Flags that determine the behaviour of the Doppler processing routines.

Functions

- [SPECTRALRADAR_API](#)
[DopplerProcessingHandle](#) [createDopplerProcessing](#) (void)
- [SPECTRALRADAR_API](#) void [createDopplerProcessingForProcessing](#) ([DopplerProcessingHandle](#) *Doppler, [ProcessingHandle](#) Proc)
- [SPECTRALRADAR_API](#) void [setDopplerPropertyInt](#) ([DopplerProcessingHandle](#) Doppler, [DopplerPropertyInt](#) Property, int Value)
Sets Doppler processing properties.
- [SPECTRALRADAR_API](#) void [setDopplerPropertyFloat](#) ([DopplerProcessingHandle](#) Doppler, [DopplerPropertyFloat](#) Property, float Value)
Sets Doppler processing properties.
- [SPECTRALRADAR_API](#) void [setDopplerFlag](#) ([DopplerProcessingHandle](#) Doppler, [DopplerFlag](#) Flag, [BOOL](#) OnOff)
Sets the Doppler processing flags.
- [SPECTRALRADAR_API](#) void [setDopplerAmplitudeOutput](#) ([DopplerProcessingHandle](#) Doppler, [DataHandle](#) AmpOut)
Sets the location of the resulting doppler amplitude output.
- [SPECTRALRADAR_API](#) void [setDopplerPhaseOutput](#) ([DopplerProcessingHandle](#) Doppler, [DataHandle](#) PhasesOut)
Sets the location of the resulting doppler phase output.
- [SPECTRALRADAR_API](#) void [executeDopplerProcessing](#) ([DopplerProcessingHandle](#) Doppler, [ComplexDataHandle](#) Input)
Executes the Doppler processing of the input data and returns phases and amplitudes.
- [SPECTRALRADAR_API](#) void [closeDopplerProcessing](#) ([DopplerProcessingHandle](#) Doppler)
Closes the Doppler processing routines and frees the memory that has been allocated for these to work properly.

5.13.1 Detailed Description

Doppler Processing Routines.

5.13.2 Typedef Documentation

5.13.2.1 DopplerProcessingHandle

Handle used for Doppler processing.

5.13.3 Enumeration Type Documentation

5.13.3.1 enum DopplerFlag

Flags that determine the behaviour of the Doppler processing routines.

Enumerator

DopplerVelocityScaling Averaging along the first axis, usually the longitudinal axis (z)

5.13.3.2 enum DopplerPropertyFloat

Values that determine the behaviour of the Doppler processing routines.

Enumerator

DopplerRefractiveIndex Averaging along the first axis, usually the longitudinal axis (z)

5.13.3.3 enum DopplerPropertyInt

Values that determine the behaviour of the Doppler processing routines.

Enumerator

DopplerAveraging_1 Averaging along the first axis, usually the longitudinal axis (z)

DopplerAveraging_2 Averaging along the first axis, usually the first transversal axis (x)

DopplerStride_1 Step size for calculating the doppler processing in the longitudinal axis (z). Stride needs to be smaller or equal to DopplerAveraging_1 and larger or equal to 1.

DopplerStride_2 Step size for calculating the doppler processing in the transversal axis (x). Stride needs to be smaller or equal to DopplerAveraging_2 and larger or equal to 1.

5.13.4 Function Documentation

5.13.4.1 void closeDopplerProcessing (DopplerProcessingHandle Doppler)

Closes the Doppler processing routines and frees the memory that has been allocated for these to work properly.

5.13.4.2 DopplerProcessingHandle createDopplerProcessing (void)

Returns

[DopplerProcessingHandle](#) to the created Doppler routines.

5.13.4.3 void createDopplerProcessingForProcessing (DopplerProcessingHandle * Doppler, ProcessingHandle Proc)

Parameters

<i>Doppler</i>	Handle to the Doppler processing.
<i>Proc</i>	Handle to the Processing routines whose results are used as input for Doppler routines.

5.13.4.4 void executeDopplerProcessing (DopplerProcessingHandle *Doppler*, ComplexDataHandle *Input*)

Executes the Doppler processing of the input data and returns phases and amplitudes.

5.13.4.5 void setDopplerAmplitudeOutput (DopplerProcessingHandle *Doppler*, DataHandle *AmpOut*)

Sets the location of the resulting doppler amplitude output.

5.13.4.6 void setDopplerFlag (DopplerProcessingHandle *Doppler*, DopplerFlag *Flag*, BOOL *OnOff*)

Sets the Doppler processing flags.

5.13.4.7 void setDopplerPhaseOutput (DopplerProcessingHandle *Doppler*, DataHandle *PhasesOut*)

Sets the location of the resulting doppler phase output.

5.13.4.8 void setDopplerPropertyFloat (DopplerProcessingHandle *Doppler*, DopplerPropertyFloat *Property*, float *Value*)

Sets Doppler processing properties.

5.13.4.9 void setDopplerPropertyInt (DopplerProcessingHandle *Doppler*, DopplerPropertyInt *Property*, int *Value*)

Sets Doppler processing properties.

5.14 Service

Service functions for additional analyzing of OCT functionality.

Functions

- [SPECTRALRADAR_API](#) void `calcContrast` ([DataHandle](#) ApodizedSpectrum, [DataHandle](#) Contrast)
Computes the contrast for the specified (apodized); spectrum.

5.14.1 Detailed Description

Service functions for additional analyzing of OCT functionality.

5.14.2 Function Documentation

5.14.2.1 void `calcContrast` ([DataHandle](#) ApodizedSpectrum, [DataHandle](#) Contrast)

Computes the contrast for the specified (apodized); spectrum.

5.15 Settings

Direct access to INI files and settings.

Typedefs

- typedef struct C_Settings * [SettingsHandle](#)

Handle for saving settings on disk.

Functions

- [SPECTRALRADAR_API SettingsHandle loadSettingsFile](#) (const char *Path)
*Loads a settings file (usually *.ini); and prepares its properties to be read.*
- [SPECTRALRADAR_API int getSettingsEntryInt](#) ([SettingsHandle](#) SettingsFile, const char *Node, int DefaultValue)
Gets an integer number from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.
- [SPECTRALRADAR_API double getSettingsEntryDouble](#) ([SettingsHandle](#) SettingsFile, const char *Node, double DefaultValue)
Gets an floating point number from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.
- [SPECTRALRADAR_API void getSettingsEntryString](#) ([SettingsHandle](#) SettingsFile, const char *Node, const char *Default, char *Data, int MaxDataSize)
Gets a string from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.
- [SPECTRALRADAR_API void setSettingsEntryInt](#) ([SettingsHandle](#) SettingsFile, const char *Node, int Value)
Sets an integer entry in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.
- [SPECTRALRADAR_API void setSettingsEntryFloat](#) ([SettingsHandle](#) SettingsFile, const char *Node, double Value)
Sets a floating point entry in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.
- [SPECTRALRADAR_API void setSettingsEntryString](#) ([SettingsHandle](#) SettingsFile, const char *Node, const char *Value)
Sets a string in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.
- [SPECTRALRADAR_API void saveSettings](#) ([SettingsHandle](#) SettingsFile)
Saves the changes to the specified Settings file.
- [SPECTRALRADAR_API void closeSettingsFile](#) ([SettingsHandle](#) SettingsFile)
Closes the specified ini file and stores the set entries (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.1 Detailed Description

Direct access to INI files and settings.

5.15.2 Typedef Documentation

5.15.2.1 SettingsHandle

Handle for saving settings on disk.

5.15.3 Function Documentation

5.15.3.1 void closeSettingsFile ([SettingsHandle](#) SettingsFile)

Closes the specified ini file and stores the set entries (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.3.2 `double getSettingsEntryDouble (SettingsHandle SettingsFile, const char * Node, double DefaultValue)`

Gets an floating point number from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.3.3 `int getSettingsEntryInt (SettingsHandle SettingsFile, const char * Node, int DefaultValue)`

Gets an integer number from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.3.4 `void getSettingsEntryString (SettingsHandle SettingsFile, const char * Node, const char * Default, char * Data, int MaxDataSize)`

Gets a string from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.3.5 `SettingsHandle loadSettingsFile (const char * Path)`

Loads a settings file (usually *.ini); and prepares its properties to be read.

5.15.3.6 `void saveSettings (SettingsHandle SettingsFile)`

Saves the changes to the specified Settings file.

5.15.3.7 `void setSettingsEntryFloat (SettingsHandle SettingsFile, const char * Node, double Value)`

Sets a floating point entry in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.3.8 `void setSettingsEntryInt (SettingsHandle SettingsFile, const char * Node, int Value)`

Sets an integer entry in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.15.3.9 `void setSettingsEntryString (SettingsHandle SettingsFile, const char * Node, const char * Value)`

Sets a string in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.

5.16 Coloring

Functions used for coloring of floating point data.

Typedefs

- typedef struct C_Coloring32Bit * [Coloring32BitHandle](#)
Handle for routines that color available scans for displaying.

Enumerations

- enum [ColorScheme](#) {
[ColorScheme_BlackAndWhite](#) = 0,
[ColorScheme_Inverted](#) = 1,
[ColorScheme_Color](#) = 2,
[ColorScheme_BlackAndOrange](#) = 3,
[ColorScheme_BlackAndRed](#) = 4,
[ColorScheme_BlackRedAndYellow](#) = 5,
[ColorScheme_DopplerPhase](#) = 6,
[ColorScheme_BlueAndBlack](#) = 7,
[ColorScheme_PolarizationRetardation](#) = 8,
[ColorScheme_GreenBlueAndBlack](#) = 9,
[ColorScheme_BlackAndRedYellow](#) = 10 }
selects the ColorScheme of the data to transform real data to colored data.
- enum [ColoringByteOrder](#) {
[Coloring_RGBA](#) = 0,
[Coloring_BGRA](#) = 1,
[Coloring_ARGB](#) = 2 }
Selects the byte order of the coloring to be applied.
- enum [ColorEnhancement](#) {
[ColorEnhancement_None](#) = 0,
[ColorEnhancement_Sine](#) = 1,
[ColorEnhancement_Parable](#) = 2,
[ColorEnhancement_Cubic](#) = 3,
[ColorEnhancement_Sqrt](#) = 4 }
Selects the byte order of the coloring to be applied.

Functions

- [SPECTRALRADAR_API](#)
[Coloring32BitHandle](#) [createColoring32Bit](#) ([ColorScheme](#) Color, [ColoringByteOrder](#) ByteOrder)
Creates processing that can be used to color given floating point B-scans to 32 bit colored images.
- [SPECTRALRADAR_API](#)
[Coloring32BitHandle](#) [createCustomColoring32Bit](#) (int LUTSize, unsigned long LUT[])
Create custom coloring using the specified color look-up-table.
- [SPECTRALRADAR_API](#) void [setColoringBoundaries](#) ([Coloring32BitHandle](#) Coloring, float Min_dB, float Max_dB)
Sets the boundaries in dB which are used by the coloring algorithm to map colors to floating point values in dB.
- [SPECTRALRADAR_API](#) void [setColoringEnhancement](#) ([Coloring32BitHandle](#) Coloring, [ColorEnhancement](#) Enhancement)
Selects a function for non-linear coloring to enhance (subjective) image impression.
- [SPECTRALRADAR_API](#) void [colorizeData](#) ([Coloring32BitHandle](#) Coloring, [DataHandle](#) Data, [ColoredDataHandle](#) ColoredData, [BOOL](#) Transpose)
Colors a given data object ([DataHandle](#)) into a given colored object ([ColoredDataHandle](#)).

- **SPECTRALRADAR_API** void **colorizeDopplerData** (**Coloring32BitHandle** AmpColoring, **Coloring32BitHandle** PhaseColoring, **DataHandle** AmpData, **DataHandle** PhaseData, **ColoredDataHandle** Output, double MinSignal_dB, **BOOL** Transpose)
*Colors a two given data object (**DataHandle**) using overlay and intensity to represent phase and amplitude data. Used for Doppler imaging.*
- **SPECTRALRADAR_API** void **clearColoring32Bit** (**Coloring32BitHandle** Coloring)
*Clears the coloring previously created by **createColoring32Bit()**.*

5.16.1 Detailed Description

Functions used for coloring of floating point data.

5.16.2 Typedef Documentation

5.16.2.1 **Coloring32BitHandle**

Handle for routines that color available scans for displaying.

5.16.3 Enumeration Type Documentation

5.16.3.1 **enum ColorEnhancement**

Selects the byte order of the coloring to be applied.

Enumerator

- ColorEnhancement_None** Use no color enhancement.
- ColorEnhancement_Sine** Apply a sine function as enhancement.
- ColorEnhancement_Parable** Apply a parable as enhancement.
- ColorEnhancement_Cubic** Apply a cubic function as enhancement.
- ColorEnhancement_Sqrt** Apply a sqrt function as enhancement.

5.16.3.2 **enum ColoringByteOrder**

Selects the byte order of the coloring to be applied.

Enumerator

- Coloring_RGBA** Byte order RGBA.
- Coloring_BGRA** Byte order BGRA.
- Coloring_ARGB** Byte order ARGB.

5.16.3.3 **enum ColorScheme**

selects the ColorScheme of the data to transform real data to colored data.

Enumerator

- ColorScheme_BlackAndWhite** Black and white (monochrome) coloring.
- ColorScheme_Inverted** Black and white inverted (monochrome inverted) coloring.
- ColorScheme_Color** colored
- ColorScheme_BlackAndOrange** orange and black coloring
- ColorScheme_BlackAndRed** red and black coloring

ColorScheme_BlackRedAndYellow black, red and yellow coloring

ColorScheme_DopplerPhase Doppler phase data coloring. Red and blue allways colored in a range from -pi to +pi. Setting the boundaries for this color scheme is only allowed inbetween +pi and -pi

ColorScheme_BlueAndBlack blue and black coloring

ColorScheme_PolarizationRetardation colorful colorscheme

5.16.4 Function Documentation

5.16.4.1 void clearColoring32Bit (Coloring32BitHandle Coloring)

Clears the coloring previously created by [createColoring32Bit\(\)](#).

5.16.4.2 void colorizeData (Coloring32BitHandle Coloring, DataHandle Data, ColoredDataHandle ColoredData, BOOL Transpose)

Colors a given data object ([DataHandle](#)) into a given colored object ([ColoredDataHandle](#)).

5.16.4.3 oid colorizeDopplerData (Coloring32BitHandle AmpColoring, Coloring32BitHandle PhaseColoring, DataHandle AmpData, DataHandle PhaseData, ColoredDataHandle Output, double MinSignal_dB, BOOL Transpose)

Colors a two given data object ([DataHandle](#)) using overlay and intensity to represent phase and amplitude data. Used for Doppler imaging.

5.16.4.4 Coloring32BitHandle createColoring32Bit (ColorScheme Color, ColoringByteOrder ByteOrder)

Creates processing that can be used to color given floating point B-scans to 32 bit colored images.

Parameters

<i>Color</i>	The color-table to be used
<i>ByteOrder</i>	The byte order the coloring is supposed to use.

Returns

The handle ([Coloring32BitHandle](#)) to the coloring algorithm.

5.16.4.5 Coloring32BitHandle createCustomColoring32Bit (int LUTSize, unsigned long LUT[])

Create custom coloring using the specified color look-up-table.

5.16.4.6 void setColoringBoundaries (Coloring32BitHandle Colorng, float Min_dB, float Max_dB)

Sets the boundaries in dB which are used by the coloring algorithm to map colors to floating point values in dB.

5.16.4.7 void setColoringEnhancement (Coloring32BitHandle Coloring, ColorEnhancement Enhancement)

Selects a function for non-linear coloring to enhance (subjective) image impression.

5.17 Camera

Functions for acquiring camera video images.

Enumerations

- enum [CameraPropertyFloat](#) {
[Camera_Saturation](#),
[Camera_Brightness](#),
[Camera_Contrast](#),
[Camera_WB_Red](#),
[Camera_WB_Green](#),
[Camera_WB_Blue](#),
[Camera_WB_Auto](#) }

Enum identifying properties of the camera.

Functions

- [SPECTRALRADAR_API](#) void [getMaxCameraImageSize](#) ([OCTDeviceHandle](#) Dev, int *SizeX, int *SizeY)
Returns the maximum possible camera image size for the current device.
- [SPECTRALRADAR_API](#) void [getCameraImage](#) ([OCTDeviceHandle](#) Dev, int SizeX, int SizeY, [ColoredDataHandle](#) Image)
Gets a camera image.
- [SPECTRALRADAR_API](#) void [getMirroredCameraImage](#) ([OCTDeviceHandle](#) Dev, int SizeX, int SizeY, [ColoredDataHandle](#) Image)
Gets a camera image.
- [SPECTRALRADAR_API](#) void [setCameraPropertyFloat](#) ([OCTDeviceHandle](#) Dev, [CameraPropertyFloat](#) Selection, double Value)
Sets saturation, brightness and contrast for the camera images if this option is available for the current device.
- [SPECTRALRADAR_API](#) void [setCameraShowScanPattern](#) ([OCTDeviceHandle](#) Dev, [BOOL](#) Value)
Enables to turn on/off the scan pattern overlay in the video camera image.
- [SPECTRALRADAR_API](#) void [visualizeScanPattern](#) ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe, [ScanPatternHandle](#) Pattern, [BOOL](#) showRawPattern)
Visualizes the scan pattern in top of the camera image.

5.17.1 Detailed Description

Functions for acquiring camera video images.

5.17.2 Enumeration Type Documentation

5.17.2.1 enum [CameraPropertyFloat](#)

Enum identifying properties of the camera.

Enumerator

- [Camera_Saturation](#)** Saturation of the video camera.
- [Camera_Brightness](#)** Brightness of the video camera.
- [Camera_Contrast](#)** Contrast of the video camera.
- [Camera_WB_Red](#)** Red white-balance value.
- [Camera_WB_Green](#)** Red white-balance value.
- [Camera_WB_Blue](#)** Red white-balance value.
- [Camera_WB_Auto](#)** Automatic setting of white balance values.

5.17.3 Function Documentation

5.17.3.1 **void getCameraImage (OCTDeviceHandle *Dev*, int *SizeX*, int *SizeY*, ColoredDataHandle *Image*)**

Gets a camera image.

5.17.3.2 **void getMaxCameraImageSize (OCTDeviceHandle *Dev*, int * *SizeX*, int * *SizeY*)**

Returns the maximum possible camera image size for the current device.

5.17.3.3 **void getMirroredCameraImage (OCTDeviceHandle *Dev*, int *SizeX*, int *SizeY*, ColoredDataHandle *Image*)**

Gets a camera image.

The returned camera image is mirrored in the X axis.

5.17.3.4 **void setCameraPropertyFloat (OCTDeviceHandle *Dev*, CameraPropertyFloat *Selection*, double *Value*)**

Sets saturation, brightness and contrast for the camera images if this option is available for the current device.

5.17.3.5 **void setCameraShowScanPattern (OCTDeviceHandle *Dev*, BOOL *Value*)**

Enables to turn on/off the scan pattern overlay in the video camera image.

5.17.3.6 **void void visualizeScanPattern (OCTDeviceHandle *Dev*, ProbeHandle *Probe*, ScanPatternHandle *Pattern*,
BOOL *showRawPattern*)**

Visualizes the scan pattern in top of the camera image.

5.18 Buffer

Functions for acquiring camera video images.

Typedefs

- typedef struct C_Buffer * [BufferHandle](#)
The BufferHandle identifies a data buffer.

Functions

- [SPECTRALRADAR_API BufferHandle createMemoryBuffer](#) (void)
Creates a buffer holding data and colored data.
- [SPECTRALRADAR_API void appendToBuffer](#) (BufferHandle, DataHandle, ColoredDataHandle)
Appends specified data and colored data to the requested buffer.
- [SPECTRALRADAR_API int getBufferSize](#) (BufferHandle)
Returns the currently available data sets in the buffer.
- [SPECTRALRADAR_API DataHandle getBufferData](#) (BufferHandle, int Index)
Returns the data in the buffer.
- [SPECTRALRADAR_API ColoredDataHandle getColoredBufferData](#) (BufferHandle, int Index)
Returns the colored data in the buffer.
- [SPECTRALRADAR_API void clearBuffer](#) (BufferHandle)
Clears the buffer and frees all data and colored data objects in it.

5.18.1 Detailed Description

Functions for acquiring camera video images.

5.18.2 Typedef Documentation

5.18.2.1 BufferHandle

The BufferHandle identifies a data buffer.

5.18.3 Function Documentation

5.18.3.1 void appendToBuffer (BufferHandle , DataHandle , ColoredDataHandle)

Appends specified data and colored data to the requested buffer.

If insufficient memory is available the oldest items in the buffer will be freed automatically.

5.18.3.2 void clearBuffer (BufferHandle)

Clears the buffer and frees all data and colored data objects in it.

5.18.3.3 BufferHandle createMemoryBuffer (void)

Creates a buffer holding data and colored data.

5.18.3.4 DataHandle getBufferData (BufferHandle , int Index)

Returns the data in the buffer.

5.18.3.5 int getBufferSize (BufferHandle)

Returns the currently available data sets in the buffer.

5.18.3.6 ColoredDataHandle getColoredBufferData (BufferHandle , int Index)

Returns the colored data in the buffer.

5.19 File Handling

Enumerations

- enum `OCTFileFormat` {
FF_Unknown,
FF_SRM,
FF_RAW,
FF_SDR,
FF_PHS,
FF_IMG,
FF_CSV,
FF_TXT,
FF_TABLETXT,
FF_OCITY,
FF_FITS,
FF_VTK,
FF_VFF,
FF_TIFF,
FF_JPG,
FF_BMP,
FF_PNG }

Enum identifying possible file formats.

- enum `DataKind` {
dkReal,
dkColored,
dkComplex,
dkRaw,
dkBinary,
dkText,
dkUnknown = 999 }

Enum identifying.

- enum `FileMetadataFloatField` {
FMD_RefractiveIndex,
FMD_RangeX,
FMD_RangeY,
FMD_RangeZ,
FMD_CenterX,
FMD_CenterY,
FMD_Angle,
FMD_BinToElectronScaling,
FMD_CentralWavelength,
FMD_SourceBandwidth,
FMD_ElectronFloor,
FMD_DynamicRange_Lower,
FMD_DynamicRange_Upper,
FMD_Rotation3D_X,
FMD_Rotation3D_Y,
FMD_ClipPlaneDepth3D,
FMD_QuadraticDispersionCorrectionFactor,
FMD_SpeckleVarianceThreshold,
FMD_ScanTime,
FMD_ReferenceIntensity,
FMD_ScanPause,
FMD_Zoom,
FMD_MinPointDistance,
FMD_MaxPointDistance,
FMD_MaxExternalTriggerFrequency }

Enum identifying file metadata fields of floating point type.

- enum `FileMetadataIntField` {
`FMD_ProcessState`,
`FMD_SizeX`,
`FMD_SizeY`,
`FMD_SizeZ`,
`FMD_Oversampling`,
`FMD_IntensityAveragedSpectra`,
`FMD_IntensityAveragedAScans`,
`FMD_IntensityAveragedBScans`,
`FMD_DopplerAverageX`,
`FMD_DopplerAverageZ`,
`FMD_ApoWindow`,
`FMD_DeviceBitDepth`,
`FMD_SpectrometerElements`,
`FMD_Colormap`,
`FMD_Aspect`,
`FMD_ExperimentNumber`,
`FMD_DevicePreset`,
`FMD_Timestamp`,
`FMD_CompressionLevel`,
`FMD_DeviceBytesPerPixel`,
`FMD_SpeckleAveragingFastAxis`,
`FMD_SpeckleAveragingSlowAxis` }

Enum identifying file metadata fields of integral type.

- enum `FileMetadataStringField` {
`FMD_AxisOrder`,
`FMD_DeviceName`,
`FMD_Serial`,
`FMD_Comment`,
`FMD_CustomInfo`,
`FMD_AcquisitionMode`,
`FMD_Study`,
`FMD_DispersionPreset`,
`FMD_ProbeName`,
`FMD_DevicePresetDescription` }

Enum identifying file metadata fields of character string type.

- enum `FileMetadataBoolField` {
`FMD_OffsetApplied`,
`FMD_DCSubtracted`,
`FMD_ApoApplied`,
`FMD_De chirpApplied`,
`FMD_UndersamplingFilterApplied`,
`FMD_DispersionCompensationApplied`,
`FMD_QuadraticDispersionCorrectionUsed`,
`FMD_ImageFieldCorrectionApplied`,
`FMD_ScanLineShown`,
`FMD_AutoCorrCompensationUsed`,
`FMD_BScanCrossCorrelation`,
`FMD_DCSubtractedAdvanced`,
`FMD_OnlyWindowing`,
`FMD_RawDataIsSigned`,
`FMD_FreeformScanPatternIsActive`,
`FMD_FreeformScanPatternCloseLoop`,
`FMD_FreeformScanPatternSplineInterpolation`,
`FMD_ExternalTriggerActive` }

Enum identifying file metadata fields of bool type.

5.19.1 Detailed Description

5.19.2 Enumeration Type Documentation

5.19.2.1 enum DataKind

Enum identifying.

5.19.2.2 enum FileMetadataBoolField

Enum identifying file metadata fields of bool type.

5.19.2.3 enum FileMetadataFloatField

Enum identifying file metadata fields of floating point type.

Enumerator

FMD_RefractiveIndex The refractive index applied to the whole image.

FMD_RangeX The FOV in longitudinal axis (z) in mm.

FMD_RangeY The FOV in axial direction (x) in mm.

FMD_RangeZ The FOV in axial direction (y) in mm.

FMD_CenterX The center of the scan pattern in axial direction (x) in mm.

FMD_CenterY The center of the scan pattern in axial direction (y) in mm.

FMD_Angle The angle between the scanner and the video camera image.

FMD_BinToElectronScaling Ratio between the binary value from the camera to the count of electrons.

FMD_CentralWavelength Central wavelength of the device.

FMD_SourceBandwidth Bandwidth of the light source.

FMD_DynamicRange_Lower Lower border of the dynamic range.

FMD_DynamicRange_Upper Upper border of the dynamic range.

FMD_SpeckleVarianceThreshold Threshold for speckle variance mode.

FMD_ScanTime Time needed for data acquisition. The processing and saving time is not included.

FMD_ReferenceIntensity Value for the reference intensity.

FMD_Zoom Zooms the scan pattern.

FMD_MinPointDistance Minimum distance between two points of the scan pattern used for freeform scan patterns.

FMD_MaxPointDistance Maximum distance between two points of the scan pattern used for freeform scan patterns.

FMD_MaxExternalTriggerFrequency maximal external trigger frequency for the chosen preset

5.19.2.4 enum FileMetadataIntField

Enum identifying file metadata fields of integral type.

Enumerator

FMD_ProcessState Contains the specific data format.

FMD_SizeX Number of pixels in x.

FMD_SizeY Number of pixels in y.

FMD_SizeZ Number of pixels in z.

FMD_Oversampling Oversampling parameter.

FMD_IntensityAveragedAScans A-scan averaging.

FMD_IntensityAveragedBScans B-scan averaging.

FMD_DopplerAverageX Averaging for doppler processing in x-direction.

FMD_DopplerAverageZ Averaging for doppler processing in z-direction.

FMD_ApoWindow Type of window used for apodization.

FMD_DeviceBitDepth Bits per pixel of the camera.

FMD_SpectrometerElements Number of elements of the spectrometer.

FMD_Colormap Actual colormap.

FMD_Aspect Ratio between pixel and real size.

FMD_ExperimentNumber Serial number of the dataset.

FMD_DevicePreset Preset for scan speed.

FMD_Timestamp Timestamp of acquired dataset.

FMD_CompressionLevel Compression used to create the zipped container (.oct-file). No compression used right now.

FMD_DeviceBytesPerPixel Bytes per pixel of the camera.

FMD_SpeckleAveragingFastAxis Averaging parameter of the fast scan axis in speckle variance mode.

FMD_SpeckleAveragingSlowAxis Averaging parameter of the slow scan axis in speckle variance mode.

5.19.2.5 enum FileMetadataStringField

Enum identifying file metadata fields of character string type.

Enumerator

FMD_AxisOrder Order of the axis, e.g. ZXY.

FMD_DeviceName Name of the OCT device.

FMD_Serial Serial number of the OCT device.

FMD_Comment Comment of the OCT data file.

FMD_AcquisitionMode Acquisition mode of the OCT data file.

FMD_Study Study of the OCT data file.

FMD_DispersionPreset Dispersion Preset of the OCT data file.

FMD_ProbeName Name of the probe.

FMD_DevicePresetDescription Description of the scan speed, e.g. "Default (48kHz)".

5.19.2.6 enum OCTFileFormat

Enum identifying possible file formats.

5.20 BETA_API

BETA_API.

Enumerations

- enum `InterpolationMethod` {
 `Interpolation_Linear`,
 `Interpolation_Spline` }
 Selects the interpolation method.
- enum `BoundaryCondition` {
 `BoundaryCondition_standard`,
 `BoundaryCondition_natural`,
 `BoundaryCondition_periodic` }
 Selects the boundary conditions for the interpolation.

5.20.1 Detailed Description

BETA_API.

5.20.2 Enumeration Type Documentation

5.20.2.1 enum `BoundaryCondition`

Selects the boundary conditions for the interpolation.

5.20.2.2 enum `InterpolationMethod`

Selects the interpolation method.

Enumerator

- `Interpolation_Linear`** Linear interpolation.
- `Interpolation_Spline`** Cubic B-Spline interpolation.

5.21 Post Processing

Algorithms and functions used for post processing of floating point data.

Enumerations

- enum `PepperFilterType` {
`PepperFilter_Horizontal`,
`PepperFilter_Vertical`,
`PepperFilter_Star`,
`PepperFilter_Block` }

Specifies the type of pepper filter to be applied.

Functions

- `SPECTRALRADAR_API` void `medianFilter` (`DataHandle` Data, int Rank)
Computes a median filter on the specified 2D data.
- `SPECTRALRADAR_API` void `levelData` (`DataHandle` Data)
Levels the specified data and removes tilt.
- `SPECTRALRADAR_API` void `pepperFilter` (`DataHandle` Data, `PepperFilterType` Type, float Threshold)
Removes pepper-noise (very low values, i. e. dark spots in the data). This enhances the visual (colored) representation of the data.

5.21.1 Detailed Description

Algorithms and functions used for post processing of floating point data.

5.21.2 Enumeration Type Documentation

5.21.2.1 enum `PepperFilterType`

Specifies the type of pepper filter to be applied.

Enumerator

`PepperFilter_Horizontal` Values along the horizontal axis are taken into account for the pepper filter.

`PepperFilter_Vertical` Values along the vertical axis are taken into account for the pepper filter.

`PepperFilter_Star` Values along the vertical and horizontal axis (star shape) are taken into account for the pepper filter.

`PepperFilter_Block` Values in a block surrounding the destination pixel are taken into account.

5.21.3 Function Documentation

5.21.3.1 void `levelData` (`DataHandle` Data)

Levels the specified data and removes tilt.

5.21.3.2 void `medianFilter` (`DataHandle` Data, int Rank)

Computes a median filter on the specified 2D data.

5.21.3.3 void pepperFilter (DataHandle *Data*, PepperFilterType *Type*, float *Threshold*)

Removes pepper-noise (very low values, i. e. dark spots in the data). This enhances the visual (colored) representation of the data.

The pepper filter compares all pixels to a mean of surrounding pixels. The surrounding pixels taking into account are specified by [PepperFilterType](#). If the pixels is lower than specified by the Threshold the pixel will be replaced by the mean.

5.22 Polarization

Polarization code.

Typedefs

- typedef struct
C_PolarizationProcessing * [PolarizationProcessingHandle](#)
Handle used for Polarization processing.

Enumerations

- enum [PolarizationPropertyFloat](#) { **PolarizationProcessing_SurfaceThreshold** }
Values that determine the behaviour of the Polarization processing routines.
- enum [PolarizationFlag](#) { [PolarizationVelocityScaling](#) }
Flags that determine the behaviour of the Polarization processing routines.

5.22.1 Detailed Description

Polarization code. Polarization Sensitive OCT Processing Routines.

5.22.2 Typedef Documentation

5.22.2.1 PolarizationProcessingHandle

Handle used for Polarization processing.

5.22.3 Enumeration Type Documentation

5.22.3.1 enum PolarizationFlag

Flags that determine the behaviour of the Polarization processing routines.

Enumerator

PolarizationVelocityScaling Averaging along the first axis, usually the longitudinal axis (z)

5.22.3.2 enum PolarizationPropertyFloat

Values that determine the behaviour of the Polarization processing routines.

6 Data Structure Documentation

6.1 ComplexFloat Struct Reference

A standard complex data type that is used to access complex data.

Data Fields

- float [data](#) [2]
data[0] is the real part and data[1] is the imaginary part.

6.1.1 Detailed Description

A standard complex data type that is used to access complex data.

6.1.2 Field Documentation

6.1.2.1 float data[2]

data[0] is the real part and data[1] is the imaginary part.

7 File Documentation

7.1 SpectralRadar.h File Reference

Header containing all functions of the Spectral Radar SDK. This SDK can be used for Callisto, Ganymede, Hyperion and Telesio devices.

Data Structures

- struct [ComplexFloat](#)
A standard complex data type that is used to access complex data.

Macros

- #define [SPECTRALRADAR_API](#) __declspec(dllimport)
Export/Import of define of DLL members.
- #define [TRUE](#) 1
TRUE for use with data type [BOOL](#).
- #define [FALSE](#) 0
FALSE for use with data type [BOOL](#).

Typedefs

- typedef int [BOOL](#)
A standard boolean data type used in the API.
- typedef struct C_RawData * [RawDataHandle](#)
Handle to an object holding the unprocessed raw data.
- typedef struct C_Data * [DataHandle](#)

- Handle to an object holding 1-, 2- or 3-dimensional floating point data.*

 - typedef struct C_ColoredData * [ColoredDataHandle](#)

Handle to an object holding 1-, 2- or 3-dimensional colored data.
- typedef struct C_ComplexData * [ComplexDataHandle](#)

Handle to an object holding complex 1-, 2- or 3-dimensional complex floating point data.
- typedef struct C_Buffer * [BufferHandle](#)

The BufferHandle identifies a data buffer.
- typedef struct C_OCTDevice * [OCTDeviceHandle](#)

The OCTDeviceHandle type is used as Handle for using the SpectralRadar.
- typedef struct C_Probe * [ProbeHandle](#)

Handle for controlling the galvo scanner.
- typedef C_ScanPattern * [ScanPatternHandle](#)

Handle for controlling the scan pattern.
- typedef struct C_Processing * [ProcessingHandle](#)

Handle for a processing routine.
- typedef struct C_DopplerProcessing * [DopplerProcessingHandle](#)

Handle used for Doppler processing.
- typedef struct C_PolarizationProcessing * [PolarizationProcessingHandle](#)

Handle used for Polarization processing.
- typedef struct C_Coloring32Bit * [Coloring32BitHandle](#)

Handle for routines that color available scans for displaying.
- typedef struct C_ImageFieldCorrection * [ImageFieldHandle](#)

Handle to the image field description.
- typedef struct C_VisualCalibration * [VisualCalibrationHandle](#)

Handle to the visual galvo calibration class.
- typedef struct C_MarkerList * [MarkerListHandle](#)

Handle to the marker list class.
- typedef struct C_FileHandling * [OCTFileHandle](#)

Handle to the OCT file class.
- typedef struct C_Settings * [SettingsHandle](#)

Handle for saving settings on disk.
- typedef struct C_SpeckleVariance * **SpeckleVarianceHandle**
- typedef struct C_FullRange * **FullRangeHandle**
- typedef void(__stdcall * **cbRefstageStatusChanged**)(RefstageStatus)
- typedef void(__stdcall * **cbRefstagePositionChanged**)(double)

Enumerations

- enum [ErrorCode](#) {
[NoError](#) = 0x00000000,
[Error](#) = 0xE0000000 }
- This enum is used to describe errors that occur when operating an OCT device.*
- enum [RawDataPropertyInt](#) {
[RawData_Size1](#),
[RawData_Size2](#),
[RawData_Size3](#),
[RawData_NumberOfElements](#),
[RawData_SizeInBytes](#),
[RawData_BytesPerElement](#),
[RawData_LostFrames](#) }

Specifies properties of RawData.

- enum [DataPropertyInt](#) {
[Data_Dimensions](#),
[Data_Size1](#),
[Data_Size2](#),
[Data_Size3](#),
[Data_NumberOfElements](#),
[Data_SizeInBytes](#),
[Data_BytesPerElement](#) }

Selects integer point data property.

- enum [DataPropertyFloat](#) {
[Data_Spacing1](#),
[Data_Spacing2](#),
[Data_Spacing3](#),
[Data_Range1](#),
[Data_Range2](#),
[Data_Range3](#) }

Selects floating point data property.

- enum [DataAnalyzation](#) {
[Data_Min](#),
[Data_Mean](#),
[Data_Max](#),
[Data_MaxDepth](#) }

Selects data property to analyze.

- enum [AScanAnalyzation](#) {
[Data_Noise_dB](#),
[Data_Noise_electrons](#),
[Data_PeakPos_Pixel](#),
[Data_PeakPos_PhysUnits](#),
[Data_PeakHeight_dB](#),
[Data_PeakWidth_6dB](#),
[Data_PeakWidth_20dB](#),
[Data_PeakWidth_40dB](#) }

Selects an appropriate A-Scan analyzation.

- enum [DataOrientation](#) {
[DataOrientation_ZXY](#),
[DataOrientation_ZYX](#),
[DataOrientation_XZY](#),
[DataOrientation_XYZ](#),
[DataOrientation_YXZ](#),
[DataOrientation_YZX](#),
[DataOrientation_ZTX](#),
[DataOrientation_ZXT](#) }

Selects the orientation of the data.

- enum [DevicePropertyFloat](#) {
[Device_FullWellCapacity](#),
[Device_zSpacing](#),
[Device_zRange](#),
[Device_SignalAmplitudeMin_dB](#),
[Device_SignalAmplitudeLow_dB](#),
[Device_SignalAmplitudeHigh_dB](#),
[Device_SignalAmplitudeMax_dB](#),
[Device_BinToElectronScaling](#),
[Device_Temperature](#),
[Device_SLD_OnTime_sec](#),
[Device_CenterWavelength_nm](#),
[Device_SpectralWidth_nm](#),

`Device_MaxTriggerFrequency_Hz` }

Properties of the device that can be read or measured.

- enum `DevicePropertyInt` {
`Device_SpectrumElements`,
`Device_BytesPerElement`,
`Device_MaxLiveVolumeRenderingScans`,
`Device_BitDepth`,
`Device_DatalsSigned` }

Properties of the device that can be read or measured.

- enum `ScanAxis` {
`ScanAxis_X` = 0,
`ScanAxis_Y` = 1 }

used to select the axis for manual galvo operations.

- enum `Device_CameraPreset` {
`Device_CameraPreset_Default`,
`Device_CameraPreset_1`,
`Device_CameraPreset_2`,
`Device_CameraPreset_3`,
`Device_CameraPreset_4` }

Enum identifying sensitivity and acquisition speed of the device.

- enum `ProbeParameterFloat` {
`Probe_FactorX`,
`Probe_OffsetX`,
`Probe_FactorY`,
`Probe_OffsetY`,
`Probe_FlybackTime_Sec`,
`Probe_ExpansionTime_Sec`,
`Probe_RotationTime_Sec`,
`Probe_ExpectedScanRate_Hz`,
`Probe_CameraScalingX`,
`Probe_CameraOffsetX`,
`Probe_CameraScalingY`,
`Probe_CameraOffsetY`,
`Probe_CameraAngle`,
`Probe_WhiteBalanceRed`,
`Probe_WhiteBalanceGreen`,
`Probe_WhiteBalanceBlue`,
`Probe_RangeMaxX`,
`Probe_RangeMaxY`,
`Probe_MaximumSlope_XY`,
`Probe_SpeckleSize`,
`Probe_ApoPosX`,
`Probe_ApoPosY`,
`Probe_ReferenceStageOffset` }

Parameters describing the behaviour of the Probe, such as calibration factors and scan parameters.

- enum `ProbeParameterInt` {
`Probe_ApodizationCycles`,
`Probe_Oversampling`,
`Probe_WhiteBalanceAutomatic`,
`Probe_Oversampling_SlowAxis`,
`Probe_SpeckleReduction`,
`Probe_MaxScanRangeShape` }

Parameters describing the behaviour of the Probe, such as calibration factors and scan parameters.

- enum `ProbeFlag` {
`Probe_CameraInverted_X`,
`Probe_CameraInverted_Y`,
`Probe_HasMEMSScanner` }

Boolean parameters describing the behaviour of the Probe.

- enum `ProbeType` {
`ProbeType_Standard`,
`ProbeType_Handheld`,
`ProbeType_Scientific` }

Determines the kind of probe types.

- enum `AcquisitionType` {
`Acquisition_AsyncContinuous`,
`Acquisition_AsyncFinite`,
`Acquisition_Sync` }

Determines the kind of acquisition process.

- enum `ProcessingType` {
`Processing_StandardFFT`,
`Processing_StandardNDFT`,
`Processing_iFFT1`,
`Processing_iFFT2`,
`Processing_iFFT3`,
`Processing_iFFT4`,
`Processing_NFFT1`,
`Processing_NFFT2`,
`Processing_NFFT3`,
`Processing_NFFT4` }

defines the algorithm used for dechirping the input signal and Fourier transformation

- enum `ApodizationWindow` {
`Apodization_Hann` = 0,
`Apodization_Hamming` = 1,
`Apodization_Gauss` = 2,
`Apodization_TaperedCosine` = 3,
`Apodization_Blackman` = 4,
`Apodization_BlackmanHarris` = 5,
`Apodization_LightSourceBased` = 6,
`Apodization_Unknown` = 999 }

To select the apodization window function.

- enum `ProcessingParameterInt` {
`Processing_SpectrumAveraging`,
`Processing_AScanAveraging`,
`Processing_BScanAveraging`,
`Processing_ZeroPadding`,
`Processing_NumberOfThreads`,
`Processing_FourierAveraging` }

Parameters that set the behaviour of the processing algorithms.

- enum `ProcessingParameterFloat` {
`Processing_ApodizationDamping`,
`Processing_MinElectrons` }

Parameters that set the behaviour of the processing algorithms.

- enum `CalibrationData` {
`Calibration_OffsetErrors`,
`Calibration_ApodizationSpectrum`,
`Calibration_ApodizationVector`,
`Calibration_Dispersion`,
`Calibration_Chirp`,
`Calibration_ExtendedAdjust`,
`Calibration_FixedPattern` }

Data describing the calibration of the processing routines.

- enum `ProcessingFlag` {
`Processing_UseOffsetErrors`,
`Processing_RemoveDCSpectrum`,
`Processing_RemoveAdvancedDCSpectrum`,
`Processing_UseApodization`,
`Processing_UseScanForApodization`,
`Processing_UseUndersamplingFilter`,
`Processing_UseDispersionCompensation`,
`Processing_UseDechirp`,
`Processing_UseExtendedAdjust`,
`Processing_FullRangeOutput`,
`Processing_FilterDC`,
`Processing_UseAutocorrCompensation`,
`Processing_UseDEFR`,
`Processing_OnlyWindowing`,
`Processing_RemoveFixedPattern` }

Flags that set the behaviour of the processing algorithms.

- enum `ProcessingAveragingAlgorithm` {
`Processing_Averaging_Min`,
`Processing_Averaging_Mean`,
`Processing_Averaging_Median`,
`Processing_Averaging_Norm2`,
`Processing_Averaging_Max`,
`Processing_Averaging_Fourier_Min`,
`Processing_Averaging_Fourier_Norm4`,
`Processing_Averaging_Fourier_Max`,
`Processing_Averaging_StandardDeviationAbs` }

This sets the averaging algorithm to be used for processing.

- enum `ApodizationWindowParameter` {
`ApodizationWindowParameter_Sigma`,
`ApodizationWindowParameter_Ratio`,
`ApodizationWindowParameter_Frequency` }

Sets certain parameters that are used by the window functions to be applied during apodization.

- enum `Data1DExportFormat` {
`Data1DExport_RAW`,
`Data1DExport_TXT`,
`Data1DExport_CSV`,
`Data1DExport_TableTXT`,
`Data1DExport_Fits` }

Export format for 1-dimensional data ([DataHandle](#)).

- enum `Data2DExportFormat` {
`Data2DExport_SRM`,
`Data2DExport_RAW`,
`Data2DExport_TXT`,
`Data2DExport_CSV`,
`Data2DExport_TableTXT`,
`Data2DExport_Fits` }

Export format for 2-dimensional data ([DataHandle](#)).

- enum `Data3DExportFormat` {
`Data3DExport_SRM`,
`Data3DExport_RAW`,
`Data3DExport_TXT`,
`Data3DExport_CSV`,
`Data3DExport_VFF`,
`Data3DExport_VTK`,
`Data3DExport_Fits`,
`Data3DExport_TIFF` }

- Export format for 3-dimensional data ([DataHandle](#)).*

 - enum [ComplexDataExportFormat](#) { [ComplexDataExport_RAW](#) }

Export format for complex data.
- enum [ColoredDataExportFormat](#) {
[ColoredDataExport_SRM](#),
[ColoredDataExport_RAW](#),
[ColoredDataExport_BMP](#),
[ColoredDataExport_PNG](#),
[ColoredDataExport_JPG](#),
[ColoredDataExport_PDF](#),
[ColoredDataExport_TIFF](#) }

Export format for images ([ColoredDataHandle](#)).
- enum [Direction](#) {
[Direction_1](#),
[Direction_2](#),
[Direction_3](#) }

Specifies a direction.
- enum [DataImportFormat](#) { [DataImport_SRM](#) }

Supported import format to load data from disk.
- enum [RawDataExportFormat](#) {
[RawDataExport_RAW](#),
[RawDataExport_SRR](#) }

Supported raw data export formats to store data to disk.
- enum [RawDataImportFormat](#) { [RawDataImport_SRR](#) }

Supported raw data import formats to load data from disk.
- enum [Plane2D](#) {
[Plane2D_12](#),
[Plane2D_23](#),
[Plane2D_13](#) }

Planes for slices of the volume data.
- enum **CornerPos** {
NotRecognized,
BottomLeft,
BottomRight,
TopRight,
TopLeft }
- enum [DopplerPropertyInt](#) {
[DopplerAveraging_1](#),
[DopplerAveraging_2](#),
[DopplerStride_1](#),
[DopplerStride_2](#) }

Values that determine the behaviour of the Doppler processing routines.
- enum [DopplerPropertyFloat](#) { [DopplerRefractiveIndex](#) }

Values that determine the behaviour of the Doppler processing routines.
- enum [DopplerFlag](#) { [DopplerVelocityScaling](#) }

Flags that determine the behaviour of the Doppler processing routines.
- enum [ColorScheme](#) {
[ColorScheme_BlackAndWhite](#) = 0,
[ColorScheme_Inverted](#) = 1,
[ColorScheme_Color](#) = 2,
[ColorScheme_BlackAndOrange](#) = 3,
[ColorScheme_BlackAndRed](#) = 4,
[ColorScheme_BlackRedAndYellow](#) = 5,
[ColorScheme_DopplerPhase](#) = 6,
[ColorScheme_BlueAndBlack](#) = 7,
[ColorScheme_PolarizationRetardation](#) = 8,
ColorScheme_GreenBlueAndBlack = 9,

ColorScheme_BlackAndRedYellow = 10 }

selects the ColorScheme of the data to transform real data to colored data.

- enum **ColoringByteOrder** {
Coloring_RGBA = 0,
Coloring_BGRA = 1,
Coloring_ARGB = 2 }

Selects the byte order of the coloring to be applied.

- enum **ColorEnhancement** {
ColorEnhancement_None = 0,
ColorEnhancement_Sine = 1,
ColorEnhancement_Parable = 2,
ColorEnhancement_Cubic = 3,
ColorEnhancement_Sqrt = 4 }

Selects the byte order of the coloring to be applied.

- enum **CameraPropertyFloat** {
Camera_Saturation,
Camera_Brightness,
Camera_Contrast,
Camera_WB_Red,
Camera_WB_Green,
Camera_WB_Blue,
Camera_WB_Auto }

Enum identifying properties of the camera.

- enum **OCTFileFormat** {
FF_Unknown,
FF_SRM,
FF_RAW,
FF_SDR,
FF_PHS,
FF_IMG,
FF_CSV,
FF_TXT,
FF_TABLETXT,
FF_OCITY,
FF_FITS,
FF_VTK,
FF_VFF,
FF_TIFF,
FF_JPG,
FF_BMP,
FF_PNG }

Enum identifying possible file formats.

- enum **DataKind** {
dkReal,
dkColored,
dkComplex,
dkRaw,
dkBinary,
dkText,
dkUnknown = 999 }

Enum identifying.

- enum **FileMetadataFloatField** {

```

FMD_RefractiveIndex,
FMD_RangeX,
FMD_RangeY,
FMD_RangeZ,
FMD_CenterX,
FMD_CenterY,
FMD_Angle,
FMD_BinToElectronScaling,
FMD_CentralWavelength,
FMD_SourceBandwidth,
FMD_ElectronFloor,
FMD_DynamicRange_Lower,
FMD_DynamicRange_Upper,
FMD_Rotation3D_X,
FMD_Rotation3D_Y,
FMD_ClipPlaneDepth3D,
FMD_QuadraticDispersionCorrectionFactor,
FMD_SpeckleVarianceThreshold,
FMD_ScanTime,
FMD_ReferenceIntensity,
FMD_ScanPause,
FMD_Zoom,
FMD_MinPointDistance,
FMD_MaxPointDistance,
FMD_MaxExternalTriggerFrequency }

```

Enum identifying file metadata fields of floating point type.

- enum `FileMetadataIntField` {


```

FMD_ProcessState,
FMD_SizeX,
FMD_SizeY,
FMD_SizeZ,
FMD_Oversampling,
FMD_IntensityAveragedSpectra,
FMD_IntensityAveragedAScans,
FMD_IntensityAveragedBScans,
FMD_DopplerAverageX,
FMD_DopplerAverageZ,
FMD_ApoWindow,
FMD_DeviceBitDepth,
FMD_SpectrometerElements,
FMD_Colormap,
FMD_Aspect,
FMD_ExperimentNumber,
FMD_DevicePreset,
FMD_Timestamp,
FMD_CompressionLevel,
FMD_DeviceBytesPerPixel,
FMD_SpeckleAveragingFastAxis,
FMD_SpeckleAveragingSlowAxis }

```

Enum identifying file metadata fields of integral type.

- enum `FileMetadataStringField` {

```

FMD_AxisOrder,
FMD_DeviceName,
FMD_Serial,
FMD_Comment,
FMD_CustomInfo,
FMD_AcquisitionMode,
FMD_Study,
FMD_DispersionPreset,
FMD_ProbeName,
FMD_DevicePresetDescription }

```

Enum identifying file metadata fields of character string type.

- enum **FileMetadataBoolField** {


```

FMD_OffsetApplied,
FMD_DCSubtracted,
FMD_ApoApplied,
FMD_De chirpApplied,
FMD_UndersamplingFilterApplied,
FMD_DispersionCompensationApplied,
FMD_QuadraticDispersionCorrectionUsed,
FMD_ImageFieldCorrectionApplied,
FMD_ScanLineShown,
FMD_AutoCorrCompensationUsed,
FMD_BScanCrossCorrelation,
FMD_DCSubtractedAdvanced,
FMD_OnlyWindowing,
FMD_RawDataIsSigned,
FMD_FreeformScanPatternsActive,
FMD_FreeformScanPatternCloseLoop,
FMD_FreeformScanPatternSplineInterpolation,
FMD_ExternalTriggerActive }

```

Enum identifying file metadata fields of bool type.

- enum **SpeckleVarianceType** {


```

SpeckleVariance_LogscaleVariance_Linear,
SpeckleVariance_LogscaleVariance_Logscale,
SpeckleVariance_LinearVariance_Linear,
SpeckleVariance_LinearVariance_Logscale,
SpeckleVariance_ComplexVariance_Linear,
SpeckleVariance_ComplexVariance_Logscale }

```
- enum **Device_TriggerType** {


```

Trigger_FreeRunning,
Trigger_TrigBoard_ExternalStart,
Trigger_External_AScan }

```
- enum **ScanPatternPropertyInt** {


```

ScanPattern_SizeTotal,
ScanPattern_Cycles,
ScanPattern_SizeCycle,
ScanPattern_SizePreparationCycle,
ScanPattern_SizeImagingCycle }

```
- enum **ScanPattern_AcquisitionOrder** {


```

ScanPattern_AcqOrderFrameByFrame,
ScanPattern_AcqOrderAll }

```
- enum **ScanPatternPropertyFloat** {


```

ScanPattern_RangeX,
ScanPattern_RangeY,
ScanPattern_CenterX,
ScanPattern_CenterY,
ScanPattern_Angle }

```

- enum [InterpolationMethod](#) {
[Interpolation_Linear](#),
[Interpolation_Spline](#) }
Selects the interpolation method.
- enum [BoundaryCondition](#) {
BoundaryCondition_standard,
BoundaryCondition_natural,
BoundaryCondition_periodic }
Selects the boundary conditions for the interpolation.
- enum [PepperFilterType](#) {
[PepperFilter_Horizontal](#),
[PepperFilter_Vertical](#),
[PepperFilter_Star](#),
[PepperFilter_Block](#) }
Specifies the type of pepper filter to be applied.
- enum **GaussianFilterType** {
GaussianFilter_3x3,
GaussianFilter_5x5 }
- enum **LaplacianFilterType** { **LaplacianFilter_3x3** }
- enum **PrewittFilterType** {
PrewittFilter_Horizontal_3x3,
PrewittFilter_Vertical_3x3 }
- enum **SobelFilterType** {
SobelFilter_Horizontal_3x3,
SobelFilter_Vertical_3x3 }
- enum **FilterType** {
Filter_DarkField,
Filter_BrightField,
Filter_PhaseContrast }
- enum **PolarizationPropertyInt**
- enum [PolarizationPropertyFloat](#) { **PolarizationProcessing_SurfaceThreshold** }
Values that determine the behaviour of the Polarization processing routines.
- enum [PolarizationFlag](#) { [PolarizationVelocityScaling](#) }
Flags that determine the behaviour of the Polarization processing routines.
- enum **USBProbeButtonID** {
USB_BTN1 = 0xA0u,
USB_BTN2 = 0xA1u,
USB_BTN3 = 0xA2u,
USB_BTN4 = 0xA3u,
USB_BTN5 = 0xA4u }
- enum **USBProbeMessage** {
USB_INT_ERR_INVALID_CFG = 0x10,
USB_INT_BTN = 0xA0u,
USB_INT_BTN1 = 0xA0u,
USB_INT_BTN2 = 0xA1u,
USB_INT_BTN3 = 0xA2u,
USB_INT_BTN4 = 0xA3u,
USB_INT_BTN5 = 0xA4u,
USB_INT_MSG_START = 0xB0u,
USB_INT_MSG_STOP = 0xB1u,
USB_INT_MSG_SNAPSHOT = 0xB2u,
USB_INT_MSG_LED1_ON = 0xB3u,
USB_INT_MSG_LED1_OFF = 0xB4u,
USB_INT_MSG_LED2_ON = 0xB5u,
USB_INT_MSG_LED2_OFF = 0xB6u,
USB_INT_MSG_FILESIZE = 0xB7u,
USB_INT_MSG_FIRMWARE_VER = 0xB8u,

- ```

USB_INT_MSG_BITMAPSIZE = 0xB9u }
• enum USBProbeCommand {
 USB_CMD_CONFIGURE = 0x80u,
 USB_CMD_RINGLIGHT_SET = 0x81u,
 USB_CMD_LED1_TOGGLE = 0x82u,
 USB_CMD_LED2_TOGGLE = 0x83u,
 USB_CMD_GET_FIRMWARE_VER = 0x84u,
 USB_CMD_READ_USER_SIG = 0x87u,
 USB_CMD_READ_FILE = 0x88u,
 USB_CMD_READ_KEY_BITMAP = 0x89u }
• enum USBProbeConfiguration {
 USB_CFG_RINGLIGHT_UP = 0x90u,
 USB_CFG_RINGLIGHT_DOWN = 0x91u,
 USB_CFG_ACQ_STARTSTOP = 0x92u,
 USB_CFG_ACQ_SNAPSHOT = 0x93u,
 USB_CFG_NOFUNCTION = 0x9fu }
• enum USBProbeMisc {
 USB_INT_NOFUNCTION = 0x9Fu,
 USB_INT_ACK = 0x99u }
• enum RefstageStatus {
 REFSTAGE_STATUS_IDLE = 0,
 REFSTAGE_STATUS_HOMING = 1,
 REFSTAGE_STATUS_MOVING = 2,
 REFSTAGE_STATUS_MOVING_TO = 3,
 REFSTAGE_STATUS_UNDEFINED = -1 }
• enum RefstageSpeed {
 REFSTAGE_SPEED_SLOW = 0,
 REFSTAGE_SPEED_FAST = 1,
 REFSTAGE_SPEED_SLOWER = 2,
 REFSTAGE_SPEED_FASTER = 3 }

```

## Functions

- **SPECTRALRADAR\_API** unsigned long **InterpretReferenceIntensity** (float intensity)  
*interprets the reference intensity and gives a color code that reflects its state.*
- **SPECTRALRADAR\_API** **ErrorCode** **getError** (char \*Message, int StringSize)  
*Returns an error code and a message if an error occurred. The error flag will be cleared.*
- **SPECTRALRADAR\_API** int **getDataPropertyInt** (**DataHandle** Data, **DataPropertyInt** Selection)  
*Returns the selected integer property of the specified data.*
- **SPECTRALRADAR\_API** double **getDataPropertyFloat** (**DataHandle** Data, **DataPropertyFloat** Selection)  
*Returns the selected floating point property of the specified data.*
- **SPECTRALRADAR\_API** void **copyData** (**DataHandle** DataSource, **DataHandle** DataDestination)  
*Copies the content of the specified source to the specified destination.*
- **SPECTRALRADAR\_API** void **copyDataContent** (**DataHandle** DataSource, float \*Destination)  
*Copies the data in the specified data object (**DataHandle**) into the specified pointer.*
- **SPECTRALRADAR\_API** float \* **getDataPtr** (**DataHandle** Data)  
*Returns a pointer to the content of the specified data.*
- **SPECTRALRADAR\_API** void **reserveData** (**DataHandle** Data, int Size1, int Size2, int Size3)  
*Reserves the amount of data specified. This might improve performance if appending data to the **DataHandle** as no additional memory needs to be reserved then.*
- **SPECTRALRADAR\_API** void **resizeData** (**DataHandle** Data, int Size1, int Size2, int Size3)  
*Resizes the respective data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.*
- **SPECTRALRADAR\_API** void **setDataRange** (**DataHandle** Data, double range1, double range2, double range3)



- Sets the range in mm in the 3 axes represented in the RealData buffer.*
- **SPECTRALRADAR\_API** void **setDataContent** ([DataHandle](#) Data, float \*NewContent)  
*Sets the data content of the data object. The data chunk pointed to by NewContent needs to be of the size expected by the data object, i. e. Size1\*Size2\*Size\*sizeof(float).*
  - **SPECTRALRADAR\_API** [DataOrientation](#) **getDataOrientation** ([DataHandle](#) Data)  
*Returns the data orientation of the data object.*
  - **SPECTRALRADAR\_API** void **setDataOrientation** ([DataHandle](#) Data, [DataOrientation](#) Orientation)  
*Sets the data orientation of the data object to the given orientation.*
  - **SPECTRALRADAR\_API** int **getComplexDataPropertyInt** ([ComplexDataHandle](#) Data, [DataPropertyInt](#) Selection)  
*Returns the selected integer property of the specified data.*
  - **SPECTRALRADAR\_API** void **copyComplexDataContent** ([ComplexDataHandle](#) DataSource, [ComplexFloat](#) \*Destination)  
*Copies the content of the complex data to the pointer specified as destination.*
  - **SPECTRALRADAR\_API** [ComplexFloat](#) \* **getComplexDataPtr** ([ComplexDataHandle](#) Data)  
*Returns a pointer to the data represented by the [ComplexDataHandle](#). The data is still managed by the [ComplexDataHandle](#) object.*
  - **SPECTRALRADAR\_API** void **setComplexDataContent** ([ComplexDataHandle](#) Data, [ComplexFloat](#) \*NewContent)  
*Sets the data content of the [ComplexDataHandle](#) to the content specified by the pointer.*
  - **SPECTRALRADAR\_API** void **reserveComplexData** ([ComplexDataHandle](#) Data, int Size1, int Size2, int Size3)  
*Reserves the amount of data specified. This might improve performance if appending data to the [ComplexDataHandle](#) as no additional memory needs to be reserved then.*
  - **SPECTRALRADAR\_API** void **resizeComplexData** ([ComplexDataHandle](#) Data, int Size1, int Size2, int Size3)  
*Resizes the respective data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.*
  - **SPECTRALRADAR\_API** void **setComplexDataRange** ([ComplexDataHandle](#) Data, double range1, double range2, double range3)  
*Sets the range in mm in the 3 axes represented in the RealData buffer.*
  - **SPECTRALRADAR\_API** int **getColoredDataPropertyInt** ([ColoredDataHandle](#) ColData, [DataPropertyInt](#) Selection)  
*Returns the selected integer property of the specified colored data.*
  - **SPECTRALRADAR\_API** double **getColoredDataPropertyFloat** ([ColoredDataHandle](#) ColData, [DataPropertyFloat](#) Selection)  
*Returns the selected integer property of the specified colored data.*
  - **SPECTRALRADAR\_API** void **copyColoredData** ([ColoredDataHandle](#) ImageSource, [ColoredDataHandle](#) ImageDestination)  
*Copies the contents of the specified [ColoredDataHandle](#) to the specified destination [ColoredDataHandle](#).*
  - **SPECTRALRADAR\_API** void **copyColoredDataContent** ([ColoredDataHandle](#) Source, unsigned long \*Destination)  
*Copies the data in the specified colored data object ([ColoredDataHandle](#)) into the specified pointer.*
  - **SPECTRALRADAR\_API** void **copyColoredDataContentAligned** ([ColoredDataHandle](#) ImageSource, unsigned long \*Destination, int Alignment1)  
*Copies the data in the specified colored data object ([ColoredDataHandle](#)) into the specified pointer. This function assures the data to be aligned accordingly.*
  - **SPECTRALRADAR\_API** unsigned long \* **getColoredDataPtr** ([ColoredDataHandle](#) ColData)  
*Returns a pointer to the content of the specified [ColoredDataHandle](#).*
  - **SPECTRALRADAR\_API** void **resizeColoredData** ([ColoredDataHandle](#) ColData, int Size1, int Size2, int Size3)  
*Resizes the respective colored data object. In general the data will be 1-dimensional if Size2 and Size3 are equal to 1, 2-dimensional if Size3 is equal to 1 and 3-dimensional if all, Size1, Size2, Size3, are unequal to 1.*
  - **SPECTRALRADAR\_API** void **reserveColoredData** ([ColoredDataHandle](#) ColData, int Size1, int Size2, int Size3)  
*Reserves the amount of colored data specified. This might improve performance if appending data to the [ColoredDataHandle](#) as no additional memory needs to be reserved then.*

- **SPECTRALRADAR\_API** void **setColoredDataContent** (**ColoredDataHandle** ColData, unsigned long \*NewContent)  
*Sets the data content of the colored data object. The data chung pointed to by NewContent needs to be of the size expected by the data object, i. e. Size1\*Size2\*Size\*sizeof(unsigned long).*
- **SPECTRALRADAR\_API** void **setColoredDataRange** (**ColoredDataHandle** Data, double range1, double range2, double range3)  
*Sets the range in mm in the 3 axes represented in the data object buffer.*
- **SPECTRALRADAR\_API** **DataOrientation** **getColoredDataOrientation** (**ColoredDataHandle** Data)  
*Returns the data orientation of the colored data object.*
- **SPECTRALRADAR\_API** void **setColoredDataOrientation** (**ColoredDataHandle** Data, **DataOrientation**)  
*Sets the data orintation of the colored data object to the given orientation.*
- **SPECTRALRADAR\_API** void **getRawDataSize** (**RawDataHandle** Raw, int \*SizeX, int \*SizeY, int \*SizeZ)  
*Returns the size of the specified raw data (**RawDataHandle**).*
- **SPECTRALRADAR\_API** void **copyRawDataContent** (**RawDataHandle** RawDataSource, void \*DataContent)  
*Copies the content of the raw data into the specified buffer. The suer needs to assure that enough memory is allocated.*
- **SPECTRALRADAR\_API** void \* **getRawDataPtr** (**RawDataHandle** RawDataSource)  
*Returns the pointer to the raw data content. The pointer might no longer after additional actions using the **RawDataHandle**.*
- **SPECTRALRADAR\_API** int **getRawDataPropertyInt** (**RawDataHandle** RawData, **RawDataPropertyInt** Property)  
*Returns a raw data property.*
- **SPECTRALRADAR\_API** void **setRawDataBytesPerPixel** (**RawDataHandle** Raw, int BytesPerPixel)  
*Sets the bytes per pixel for raw data.*
- **SPECTRALRADAR\_API** void **resizeRawData** (**RawDataHandle** Raw, int Size1, int Size2, int Size3)  
*Resizes the specified raw data buffer accordingly.*
- **SPECTRALRADAR\_API** void **reserveRawData** (**RawDataHandle** Raw, int Size1, int Size2, int Size3)
- **SPECTRALRADAR\_API** void **setRawDataContent** (**RawDataHandle** RawDataSource, void \*NewContent)  
*Sets the content of the raw data buffer. The size of the **RawDataHandle** needs to be adjusted first, as otherwise not all data might be copied.*
- **SPECTRALRADAR\_API** void **setScanSpectra** (**RawDataHandle** RawData, int NumberOfScanRegions, int \*ScanRegions)  
*Sets the number of the spectra in the raw data that are used for creating A-scan/B-scan data.*
- **SPECTRALRADAR\_API** void **setApodizationSpectra** (**RawDataHandle** RawData, int NumberOfScanRegions, int \*ApodizationRegions)  
*Sets the number of the spectra in the raw data that contain data useful as apodization spectra.*
- **SPECTRALRADAR\_API** int **getNumberOfScanRegions** (**RawDataHandle** Raw)  
*Returns the number of regions that have been acquired that contain scan data, i. e. spectra that are used to compute A-scans.*
- **SPECTRALRADAR\_API** int **getNumberOfApodizationRegions** (**RawDataHandle** Raw)  
*Returns the number of regions in the raw data containing spectra that are supposed to be used for apodization.*
- **SPECTRALRADAR\_API** void **getScanSpectra** (**RawDataHandle** Raw, int \*SpectralIndex)  
*Returns the indices of spectra that contain scan data, i. e. spectra that are supposed to be used to compute A-scans.*
- **SPECTRALRADAR\_API** void **getApodizationSpectra** (**RawDataHandle** Raw, int \*SpectralIndex)  
*Returns the indices of spectra that contain apodization data, i. e. spectra that are supposed to be used as input for apodization.*
- **SPECTRALRADAR\_API** **RawDataHandle** **createRawData** (void)  
*Creates a raw data object (**RawDataHandle**).*
- **SPECTRALRADAR\_API** void **clearRawData** (**RawDataHandle** Raw)  
*Clears a raw data object (**RawDataHandle**).*
- **SPECTRALRADAR\_API** **DataHandle** **createData** (void)  
*Creates a 1-dimensional data object, containing floating point data.*
- **SPECTRALRADAR\_API** **DataHandle** **createGradientData** (int Size)

- Creates a 1-dimensional data object, containing floating point data with equidistant arranged values between [0, size-1] with distance 1/(size-1).*
- **SPECTRALRADAR\_API** void **clearData** ([DataHandle](#) Data)  
*Clears the specified [DataHandle](#), [DataHandle](#), [DataHandle](#) or [DataHandle](#) objects.*
  - **SPECTRALRADAR\_API** [ColoredDataHandle](#) **createColoredData** (void)  
*Creates a colored data object ([ColoredDataHandle](#)).*
  - **SPECTRALRADAR\_API** void **clearColoredData** ([ColoredDataHandle](#) Volume)  
*Clears a colored volume object.*
  - **SPECTRALRADAR\_API** [ComplexDataHandle](#) **createComplexData** (void)  
*Creates a data object holding complex data.*
  - **SPECTRALRADAR\_API** void **clearComplexData** ([ComplexDataHandle](#) Data)  
*Clears a data object holding complex data.*
  - **SPECTRALRADAR\_API** **BOOL** **isDeviceAvailable** ()
  - **SPECTRALRADAR\_API** void **expectedDevice** (char \*Buffer, int Size, int \*Rev)
  - **SPECTRALRADAR\_API** [OCTDeviceHandle](#) **initDevice** (void)  
*Initializes the installed device.*
  - **SPECTRALRADAR\_API** void **getDeviceType** ([OCTDeviceHandle](#) Dev, char DevName[], int BufferSize)  
*Gives the name of the device type that is given by the [OCTDeviceHandle](#).*
  - **SPECTRALRADAR\_API** int **getDeviceRevision** ([OCTDeviceHandle](#) Dev)  
*Returns the revision of the device given by the [OCTDeviceHandle](#).*
  - **SPECTRALRADAR\_API** void **getDeviceSerialNumber** ([OCTDeviceHandle](#) Dev, char DevName[], int BufferSize)  
*Returns the serial number of the device given by the [OCTDeviceHandle](#).*
  - **SPECTRALRADAR\_API** int **getDevicePropertyInt** ([OCTDeviceHandle](#) Dev, [DevicePropertyInt](#) Selection)  
*Returns properties of the device belonging to the specified [OCTDeviceHandle](#).*
  - **SPECTRALRADAR\_API** double **getDevicePropertyFloat** ([OCTDeviceHandle](#) Dev, [DevicePropertyFloat](#) Selection)  
*Returns properties of the device belonging to the specified [OCTDeviceHandle](#).*
  - **SPECTRALRADAR\_API** void **closeDevice** ([OCTDeviceHandle](#) Dev)  
*Closes the device opened previously with [initDevice](#).*
  - **SPECTRALRADAR\_API** **BOOL** **isDeviceOn** ([OCTDeviceHandle](#) Handle)  
*Returns if the device is switched on.*
  - **SPECTRALRADAR\_API** **BOOL** **isVideoCameraAvailable** ([OCTDeviceHandle](#) Dev)  
*Returns if the video camera is available.*
  - **SPECTRALRADAR\_API** **BOOL** **isSLDAvailable** ([OCTDeviceHandle](#) Dev)  
*Returns whether the SLD is available.*
  - **SPECTRALRADAR\_API** void **setSLD** ([OCTDeviceHandle](#) Dev, **BOOL** OnOff)  
*switches the SLD of the SpectralRadar device on and off.*
  - **SPECTRALRADAR\_API** void **moveScanner** ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe, [ScanAxis](#) Axis, double Position)  
*manually moves the scanner to a given position*
  - **SPECTRALRADAR\_API** void **setLaserDiode** ([OCTDeviceHandle](#) Dev, **BOOL** OnOff)  
*switches the LaserDiode of the SpectralRadar device on and off.*
  - **SPECTRALRADAR\_API** double **getWavelengthAtPixel** ([OCTDeviceHandle](#) Dev, int Pixel)  
*Returns the wavelength at a specified pixel of the spectrometer.*
  - **SPECTRALRADAR\_API** void **setCameraPreset** ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe, [ProcessingHandle](#) Proc, int Preset)
  - **SPECTRALRADAR\_API** int **getCameraPreset** ([OCTDeviceHandle](#) Dev)  
*Gets the currently used device preset.*
  - **SPECTRALRADAR\_API** void **getCameraPresetDescription** ([OCTDeviceHandle](#) Dev, int Preset, char \*Description, int BufferSize)
  - **SPECTRALRADAR\_API** int **getNumberOfCameraPresets** ([OCTDeviceHandle](#) Dev)

- **SPECTRALRADAR\_API** int **getNumberOfInternalValues** (OCTDeviceHandle Dev)  
*Returns the number of Analog-to-Digital Converter present in the device.*
- **SPECTRALRADAR\_API** void **getInternalValueName** (OCTDeviceHandle Dev, int Index, char \*Name, int NameStringSize, char \*Unit, int UnitStringSize)  
*Returns names and unit for the specified Analog-to-Digital Converter.*
- **SPECTRALRADAR\_API** double **getInternalValueByName** (OCTDeviceHandle Dev, const char \*Name)  
*Returns the value of the specified Analog-to-Digital Converter (ADC);.*
- **SPECTRALRADAR\_API** double **getInternalValueByIndex** (OCTDeviceHandle Dev, int Index)  
*Returns the value of the selected ADC.*
- **SPECTRALRADAR\_API** ProbeHandle **initProbe** (OCTDeviceHandle Dev, const char \*ProbeFile)  
*Initializes a probe specified by ProbeFile.*
- **SPECTRALRADAR\_API** ProbeHandle **initStandardProbe** (OCTDeviceHandle Dev)  
*Creates a standard probe using the Probe.ini file. If this configuration file is not found, standard parameters without valid calibration will be used.*
- **SPECTRALRADAR\_API** ProbeHandle **initProbeWithType** (OCTDeviceHandle Dev, ProbeType Type)  
*Creates a standard probe for the given probe type but without valid calibration data .*
- **SPECTRALRADAR\_API** void **saveProbe** (ProbeHandle Probe, const char \*ProbeFile)  
*Saves the current properties of the [ProbeHandle](#) to a specified INI file to be reloaded using the [initProbe\(\)](#) function.*
- **SPECTRALRADAR\_API** void **setProbeParameterInt** (ProbeHandle Probe, ProbeParameterInt Selection, int Value)  
*Sets.*
- **SPECTRALRADAR\_API** void **setProbeParameterFloat** (ProbeHandle Probe, ProbeParameterFloat Selection, double Value)  
*Sets floating point parameters of the specified probe.*
- **SPECTRALRADAR\_API** int **getProbeParameterInt** (ProbeHandle Probe, ProbeParameterInt Selection)  
*Gets integer parameters of the specified probe.*
- **SPECTRALRADAR\_API** double **getProbeParameterFloat** (ProbeHandle Probe, ProbeParameterFloat Selection)  
*Gets floating point parameters of the specified probe.*
- **SPECTRALRADAR\_API** BOOL **getProbeFlag** (ProbeHandle Probe, ProbeFlag Selection)  
*Returns the selected boolean value of the specified probe.*
- **SPECTRALRADAR\_API** void **getProbeName** (ProbeHandle Probe, char ProbeName[], int BufferSize)  
*Returns the name of the specified probe.*
- **SPECTRALRADAR\_API** void **setProbeName** (ProbeHandle Probe, const char \*ProbeName)  
*Sets the given name of the specified probe.*
- **SPECTRALRADAR\_API** void **getProbeSerialNo** (ProbeHandle Probe, char SerialNo[], int BufferSize)  
*Gets the serial number of the specified probe.*
- **SPECTRALRADAR\_API** void **setProbeSerialNo** (ProbeHandle Probe, const char \*SerialNo)  
*Gets the serial number of the specified probe.*
- **SPECTRALRADAR\_API** void **getProbeType** (ProbeHandle Probe, char Type[], int BufferSize)  
*Gets the type of the specified probe.*
- **SPECTRALRADAR\_API** void **setProbeType** (ProbeHandle Probe, const char \*Type)  
*Sets the type of the specified probe.*
- **SPECTRALRADAR\_API** void **getProbeObjective** (ProbeHandle Probe, char Objective[], int BufferSize)  
*Gets the objective of the specified probe.*
- **SPECTRALRADAR\_API** void **setProbeObjective** (ProbeHandle Probe, const char \*Objective)  
*Sets the given objective of the specified probe.*
- **SPECTRALRADAR\_API** void **closeProbe** (ProbeHandle Probe)  
*Closes the probe and frees all memory associated with it.*
- **SPECTRALRADAR\_API** void **blendEnFaceInCamera** (ProbeHandle Probe, ScanPatternHandle Pattern, ColoredDataHandle EnFace2D, ColoredDataHandle Image, float Ratio, BOOL DenseView)

*Blends the en-face image of a given volume acquisition on top of the video image. Can be used to calibrate the probe manually.*

- **SPECTRALRADAR\_API** void **CameraPixelToPosition** (**ProbeHandle** Probe, **ColoredDataHandle** Image, int PixelX, int PixelY, double \*PosX, double \*PosY)

*Computes the physical position of a camera pixel of the video camera in the probe. It needs to be assured that the device is properly calibrated.*

- **SPECTRALRADAR\_API** void **PositionToCameraPixel** (**ProbeHandle** Probe, **ColoredDataHandle** Image, double PosX, double PosY, int \*PixelX, int \*PixelY)

*Computes the pixel of the video camera corresponding to a physical position. It needs to be assured that the device is properly calibrated.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createPointScanPattern** (**ProbeHandle** Probe, int Size, double PosX, double PosY)

*Creates a scan pattern that consists of a single point (PosX, PosY). The galvo doesn't move from there. Use this pattern for point scans and/or non-scanning probes.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createNoScanPattern** (**ProbeHandle** Probe, int Scans, int NumberOfScans)

*Creates a simple scan pattern that does not move the galvo. Use this pattern for point scans and/or non-scanning probes.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createTriggerPattern** (**ProbeHandle** Probe, int Scans)

*Creates a pattern only consisting of a specified amount of trigger signals.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createBScanPattern** (**ProbeHandle** Probe, double Range, int AScans, **BOOL** apodization)

*Creates a simple B-scan pattern that moves the galvo over a specified range.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createBilateralBScanPattern** (**ProbeHandle** Probe, double Range, int AScans, double Shift)

*Creates a bilateral scan pattern. The contouring error can be influenced using the Shift parameter.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createBScanPatternManual** (**ProbeHandle** Probe, double StartX, double StartY, double StopX, double StopY, int AScans, **BOOL** apodization)

*Creates a B-scan pattern specified by start and end points.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createIdealBScanPattern** (**ProbeHandle** Probe, double Range, int AScans)

*Creates an ideal B-scan pattern assuming scanners with infinite speed. No correction factors are taken into account. This is only used for internal purposes and not as a scan pattern designed to be output to the galvo drivers.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createCirclePattern** (**ProbeHandle** Probe, double Radius, int AScans)

*Creates a circle scan pattern.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createVolumePattern** (**ProbeHandle** Probe, double RangeX, int SizeX, double RangeY, int SizeY)

*Creates a simple volume pattern.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createBScanStackPattern** (**ProbeHandle** Probe, double RangeX, int SizeX, double RangeY, int SizeY)

*Creates a simple stack pattern.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createFreeformScanPattern** (**ProbeHandle** Probe, float \*positions, int size\_x, int size\_y, **BOOL** apodization)

*Creates a freeform scan pattern based on an array of positions.*

- **SPECTRALRADAR\_API** **ScanPatternHandle** **createFragmentedScanPattern** (**ProbeHandle** Probe, int ChunkSize, int NumberOfChunks)

*Creates a scan pattern which can be used to acquire a dataset of <NumberOfChunks> times <ChunkSize> A-scans at position 0/0. The Fragmented scan pattern can be compared in structure to a B-scan stack pattern with x and y ranges of 0; however the fragmented scan pattern behaves like a volume pattern in that it shows no delay between the respective chunks.*

- **SPECTRALRADAR\_API** void **updateScanPattern** (**ScanPatternHandle** Pattern)

*Updates the specified pattern (ScanPatternHandle);.*

- **SPECTRALRADAR\_API** void **rotateScanPattern** (**ScanPatternHandle** Pattern, double Angle)



- Rotates the specified pattern (ScanPatternHandle);.*

  - SPECTRALRADAR\_API void **rotateScanPatternExt** (ScanPatternHandle Pattern, double Angle, int index)

*Rotates the scan #index (0-based) of the specified pattern (ScanPatternHandle).*
- SPECTRALRADAR\_API void **shiftScanPattern** (ScanPatternHandle Pattern, double ShiftX, double ShiftY)

*Shifts the specified pattern (ScanPatternHandle).*
- SPECTRALRADAR\_API void **shiftScanPatternExt** (ScanPatternHandle Pattern, double ShiftX, double ShiftY, BOOL ShiftApo, int Index)

*Shifts the scan #index (0-based) of the specified pattern (ScanPatternHandle).*
- SPECTRALRADAR\_API void **zoomScanPattern** (ScanPatternHandle Pattern, double Factor)

*Zooms the specified pattern (ScanPatternHandle).*
- SPECTRALRADAR\_API int **getScanPatternLUTSize** (ScanPatternHandle Pattern)

*Returns the number of data points the specified pattern (ScanPatternHandle) used.*
- SPECTRALRADAR\_API void **getScanPatternLUT** (ScanPatternHandle Pattern, double \*PosX, double \*PosY)

*Returns the actual positions to be scanned with the specified pattern (ScanPatternHandle).*
- SPECTRALRADAR\_API void **clearScanPattern** (ScanPatternHandle Pattern)

*Clears the specified scan pattern (ScanPatternHandle).*
- SPECTRALRADAR\_API size\_t **projectMemoryRequirement** (OCTDeviceHandle Handle, ScanPatternHandle Pattern, AcquisitionType type)
- SPECTRALRADAR\_API void **startMeasurement** (OCTDeviceHandle Dev, ScanPatternHandle Pattern, AcquisitionType type)

*starts a continuous measurement BScans.*
- SPECTRALRADAR\_API void **getRawData** (OCTDeviceHandle Dev, RawDataHandle RawData)

*Acquires data and stores the data unprocessed.*
- SPECTRALRADAR\_API void **getRawDataEx** (OCTDeviceHandle Dev, RawDataHandle RawData, int CameraIdx)

*Acquires data with the specific camera given with camera index and stores the data unprocessed.*
- SPECTRALRADAR\_API void **stopMeasurement** (OCTDeviceHandle Dev)

*stops the current measurement.*
- SPECTRALRADAR\_API void **measureSpectra** (OCTDeviceHandle Dev, int NumberOfSpectra, RawDataHandle Raw)

*Acquires N spectra of raw data without moving galvo scanners.*
- SPECTRALRADAR\_API void **measureSpectraEx** (OCTDeviceHandle Dev, int N, RawDataHandle Raw, int CameraIndex)

*Acquires N spectra of raw data without moving galvo scanners. Supports multiple cameras (e.g. PS-OCT).*
- SPECTRALRADAR\_API ProcessingHandle **createProcessing** (int SpectrumSize, int BytesPerRawPixel, BOOL Signed, float ScalingFactor, float MinElectrons, ProcessingType Type, float FFTOversampling)
- SPECTRALRADAR\_API ProcessingHandle **createProcessingForDevice** (OCTDeviceHandle Dev)

*Creates suitable standard processing routines for the specified device (OCTDeviceHandle).*
- SPECTRALRADAR\_API ProcessingHandle **createProcessingForDeviceEx** (OCTDeviceHandle Dev, int CameraIndex)

*Creates suitable standard processing routines for the specified device (OCTDeviceHandle) with camera index.*
- SPECTRALRADAR\_API int **getInputSize** (ProcessingHandle Proc)

*Returns the expected input size (pixels per spectrum); of the processing algorithms.*
- SPECTRALRADAR\_API int **getAScanSize** (ProcessingHandle Handle)

*gives the number of pixels in an A-Scan of the SpectralRadar device. This number is identical to the number of rows in a finished B-Scan.*
- SPECTRALRADAR\_API int **getApodizationWindow** (ProcessingHandle Proc)

*Gets the window function that is being used for apodization.*
- SPECTRALRADAR\_API void **setApodizationWindow** (ProcessingHandle Proc, ApodizationWindow Window)

*Sets the window function that is to be used for apodization. The selected function will be used in all subsequent processings.*

- [SPECTRALRADAR\\_API](#) void [setApodizationWindowParameter](#) ([ProcessingHandle](#) Proc, [ApodizationWindowParameter](#) Selection, double Value)  
*Sets the apodization window parameter, such as window width or ratio between constant and cosine part.*
- [SPECTRALRADAR\\_API](#) double [getApodizationWindowParameter](#) ([ProcessingHandle](#) Proc, [ApodizationWindowParameter](#) Selection)  
*Gets the apodization window parameter, such as window width or ratio between constant and cosine part.*
- [SPECTRALRADAR\\_API](#) void [setDechirpAlgorithm](#) ([ProcessingHandle](#) Proc, [ProcessingType](#) Type)  
*Sets the algorithm that is to be used for dechirping the input spectra.*
- [SPECTRALRADAR\\_API](#) void [setProcessingParameterInt](#) ([ProcessingHandle](#) Proc, [ProcessingParameterInt](#) Selection, int Value)  
*Sets the specified integer value processing parameter.*
- [SPECTRALRADAR\\_API](#) int [getProcessingParameterInt](#) ([ProcessingHandle](#) Proc, [ProcessingParameterInt](#) Selection)  
*Returns the specified integer value processing parameter.*
- [SPECTRALRADAR\\_API](#) void [setProcessingParameterFloat](#) ([ProcessingHandle](#) Proc, [ProcessingParameterFloat](#) Selection, double Value)  
*Sets the specified processing parameter.*
- [SPECTRALRADAR\\_API](#) double [getProcessingParameterFloat](#) ([ProcessingHandle](#) Proc, [ProcessingParameterFloat](#) Selection)  
*Gets the specified processing parameter.*
- [SPECTRALRADAR\\_API](#) void [setProcessingFlag](#) ([ProcessingHandle](#) Proc, [ProcessingFlag](#) Flag, [BOOL](#) Value)  
*Sets the specified processing flag.*
- [SPECTRALRADAR\\_API](#) [BOOL](#) [getProcessingFlag](#) ([ProcessingHandle](#) Proc, [ProcessingFlag](#) Flag)  
*Returns TRUE if the specified processing flag is set, FALSE otherwise.*
- [SPECTRALRADAR\\_API](#) void [setProcessingAveragingAlgorithm](#) ([ProcessingHandle](#) Proc, [ProcessingAveragingAlgorithm](#) Algorithm)  
*Sets the algorithm that is used for averaging by the processing.*
- [SPECTRALRADAR\\_API](#) void [setCalibration](#) ([ProcessingHandle](#) Proc, [CalibrationData](#) Selection, [DataHandle](#) Data)  
*Sets the current active calibration data.*
- [SPECTRALRADAR\\_API](#) void [getCalibration](#) ([ProcessingHandle](#) Proc, [CalibrationData](#) Selection, [DataHandle](#) Data)  
*Returns the currently active calibration parameter.*
- [SPECTRALRADAR\\_API](#) void [measureCalibration](#) ([OCTDeviceHandle](#) Dev, [ProcessingHandle](#) Proc, [CalibrationData](#) Selection)  
*Measures the specified calibration parameters and uses them in subsequent processing.*
- [SPECTRALRADAR\\_API](#) void [measureCalibrationEx](#) ([OCTDeviceHandle](#) Dev, [ProcessingHandle](#) Proc, [CalibrationData](#) Selection, int CameraIndex)  
*Measures the specified calibration parameters and uses them in subsequent processing with specified camera index.*
- [SPECTRALRADAR\\_API](#) void [measureSpectrum](#) ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe, [ProcessingHandle](#) Proc, [BOOL](#) moveToApoPos)  
*Measures apodization spectrum and uses them in subsequent processing.*
- [SPECTRALRADAR\\_API](#) void [saveCalibrationAuto](#) ([ProcessingHandle](#) Proc, [CalibrationData](#) Selection)  
*Saves the selected calibration in its default path.*
- [SPECTRALRADAR\\_API](#) void [saveCalibration](#) ([ProcessingHandle](#) Proc, [CalibrationData](#) Selection, const char Path[])  
*Saves the selected calibration in the specified path.*
- [SPECTRALRADAR\\_API](#) void [loadCalibration](#) ([ProcessingHandle](#) Proc, [CalibrationData](#) Selection, const char Path[])  
*Will load a specified calibration file and use for subsequent processing.*
- [SPECTRALRADAR\\_API](#) void [setSpectrumOutput](#) ([ProcessingHandle](#) Proc, [DataHandle](#) Spectrum)  
*Sets the location for the resulting spectral data.*

- **SPECTRALRADAR\_API** void **setOffsetCorrectedSpectrumOutput** (**ProcessingHandle** Proc, **DataHandle** OffsetCorrectedSpectrum)  
*Sets the location for the resulting offset corrected spectral data.*
- **SPECTRALRADAR\_API** void **setDCCorrectedSpectrumOutput** (**ProcessingHandle** Proc, **DataHandle** ProcessingCorrectedSpectrum)  
*Sets the location for the resulting DC removed spectral data.*
- **SPECTRALRADAR\_API** void **setApodizedSpectrumOutput** (**ProcessingHandle** Proc, **DataHandle** Apodized-Spectrum)  
*Sets the location for the resulting apodized spectral data.*
- **SPECTRALRADAR\_API** void **setProcessedDataOutput** (**ProcessingHandle** Proc, **DataHandle** Scan)  
*Sets the pointer the resulting B-Scan of the next processing is written to.*
- **SPECTRALRADAR\_API** void **setHorMirroredDataOutput** (**ProcessingHandle** Proc, **DataHandle** Scan)  
*Sets the pointer the resulting B-Scan of the next processing is written to. The result will be written mirrored at the horizontal axis.*
- **SPECTRALRADAR\_API** void **setColoredDataOutput** (**ProcessingHandle** Proc, **ColoredDataHandle** BScan, **Coloring32BitHandle** Color)  
*Sets the pointer the resulting colored B-Scan of the next processing is written to.*
- **SPECTRALRADAR\_API** void **setTransposedColoredDataOutput** (**ProcessingHandle** Proc, **ColoredDataHandle** BScan, **Coloring32BitHandle** Color)  
*Sets the pointer the resulting colored B-Scan of the next processing is written to. The data will be transposed so that the first axis is the x-axis.*
- **SPECTRALRADAR\_API** void **setComplexDataOutput** (**ProcessingHandle** Proc, **ComplexDataHandle** ComplexBScan)  
*Sets the pointer the resulting complex B-Scan of the next processing is written to.*
- **SPECTRALRADAR\_API** void **executeProcessing** (**ProcessingHandle** Proc, **RawDataHandle** RawData)  
*Execute the processing.*
- **SPECTRALRADAR\_API** void **closeProcessing** (**ProcessingHandle** Proc)  
*Closes the processing and frees all temporary memory that was associated with it. Processing threads will be stopped.*
- **SPECTRALRADAR\_API** void **computeDispersion** (**DataHandle** Spectrum1In, **DataHandle** Spectrum2In, **DataHandle** ChirpOut, **DataHandle** DispOut)  
*Computes the dispersion and chirp of the two provided spectra, where both spectra need to have been subjected to same dispersion mismatch. Both spectra need to have been acquired for different path length differences.*
- **SPECTRALRADAR\_API** void **computeDispersionByCoeff** (double QuadraticIn, **DataHandle** ChirpIn, **DataHandle** DispOut)  
*Computes dispersion by a quadratic approximation specified by the quadratic factor.*
- **SPECTRALRADAR\_API** void **computeDispersionByImage** (**DataHandle** LinearKSpectralIn, **DataHandle** ChirpIn, **DataHandle** DispOut)  
*Guesses the dispersion based on the raw data specified. The raw data needs to be linearized in k before applying to this function.*
- **SPECTRALRADAR\_API** int **getNumberOfDispersionPresets** (**ProcessingHandle** Proc)  
*Gets the number of dispersion presets.*
- **SPECTRALRADAR\_API** const char \* **getDispersionPresetName** (**ProcessingHandle** Proc, int Index)  
*Gets the name of the dispersion preset specified with index.*
- **SPECTRALRADAR\_API** void **setDispersionPresetByName** (**ProcessingHandle** Proc, const char \*Name)  
*Sets the dispersion preset specified with name.*
- **SPECTRALRADAR\_API** void **setDispersionPresetByIndex** (**ProcessingHandle** Proc, int Index)  
*Sets the dispersion preset specified with index.*
- **SPECTRALRADAR\_API** void **setDispersionPresets** (**ProcessingHandle** Proc, **ProbeHandle** Probe)  
*Sets the dispersion presets for the probe.*
- **SPECTRALRADAR\_API** void **exportData1D** (**DataHandle** Data, **Data1DExportFormat** Format, const char \*Path)  
*Exports 1-dimensional data (**DataHandle**).*



- **SPECTRALRADAR\_API** void **exportData2D** ([DataHandle](#) Data, [Data2DExportFormat](#) Format, const char \*Path)  
*Exports 2-dimensional data ([DataHandle](#)).*
- **SPECTRALRADAR\_API** void **exportData3D** ([DataHandle](#) Volume, [Data3DExportFormat](#) Format, const char \*Path)  
*Exports 3-dimensional data ([DataHandle](#)).*
- **SPECTRALRADAR\_API** void **exportData2DAsImage** ([DataHandle](#) Data, [Coloring32BitHandle](#) Color, [ColoredDataExportFormat](#) format, const char \*fileName, [BOOL](#) drawScale, [BOOL](#) drawMarkers, [BOOL](#) physicalAspectRatio)
- **SPECTRALRADAR\_API** void **exportData3DAsImage** ([DataHandle](#) Data, [Coloring32BitHandle](#) Color, [ColoredDataExportFormat](#) format, [Direction](#) SliceNormalDirection, const char \*fileName, [BOOL](#) drawScale, [BOOL](#) drawMarkers, [BOOL](#) physicalAspectRatio)
- **SPECTRALRADAR\_API** void **importData** ([DataHandle](#) Data, [DataImportFormat](#) Format, const char \*Path)  
*Imports data with the specified format and copies it into a dat data object ([DataHandle](#)).*
- **SPECTRALRADAR\_API** void **exportComplexData** ([ComplexDataHandle](#), [ComplexDataExportFormat](#), const char \*)  
*Exports 1-, 2- and 3-dimensional complex data ([ComplexDataHandle](#)).*
- **SPECTRALRADAR\_API** void **exportColoredData** ([ColoredDataHandle](#) Image, [ColoredDataExportFormat](#) Format, const char \*fileName)  
*Exports colored data ([ColoredDataHandle](#)).*
- **SPECTRALRADAR\_API** void **importColoredData** ([ColoredDataHandle](#) ColoredData, [DataImportFormat](#) Format, const char \*Path)  
*Imports colored data ([ColoredDataHandle](#)) with the specified format and copied it into a data object ([ColoredDataHandle](#)).*
- **SPECTRALRADAR\_API** void **exportRawData** ([RawDataHandle](#) Raw, [RawDataExportFormat](#) Format, const char \*Path)  
*Exports the specified data to disk.*
- **SPECTRALRADAR\_API** void **importRawData** ([RawDataHandle](#) Raw, [RawDataImportFormat](#) Format, const char \*Path)  
*Imports the specified data from disk.*
- **SPECTRALRADAR\_API** void **appendRawData** ([RawDataHandle](#) Data, [RawDataHandle](#) DataToAppend, [Direction](#) direction)  
*Appends the new raw data to the old raw data in the specified direction.*
- **SPECTRALRADAR\_API** void **getRawDataSliceIndex** ([RawDataHandle](#) Data, [RawDataHandle](#) Slice, [Direction](#) SliceNormalDirection, int Index)  
*Returns a slice of raw data in the specified direction at the specified index.*
- **SPECTRALRADAR\_API** double **analyzeData** ([DataHandle](#) Data, [DataAnalyzation](#) Selection)  
*Performs the selected analyzation of the specified data and returns the resulting value.*
- **SPECTRALRADAR\_API** double **analyzeAScan** ([DataHandle](#) Data, [AScanAnalyzation](#) Selection)  
*Performs the selected analyzation of the specified A-scan and returns the resulting value.*
- **SPECTRALRADAR\_API** void **determineDynamicRange** ([DataHandle](#) Data, float \*MinRange\_dB, float \*MaxRange\_dB)  
*Gives a rough estimation of the dynamic range of the specified data object.*
- **SPECTRALRADAR\_API** void **transpose** ([DataHandle](#) DataIn, [DataHandle](#) DataOut)  
*Transposes the given data and writes the result to DataOut.*
- **SPECTRALRADAR\_API** void **transposenplace** ([DataHandle](#) Data)  
*Transposes the given Data.*
- **SPECTRALRADAR\_API** void **transposeAndScaleData** ([DataHandle](#) DataIn, [DataHandle](#) DataOut, float Min, float Max)  
*Transposes the given data and scales it to the range [Min, Max].*
- **SPECTRALRADAR\_API** void **normalizeData** ([DataHandle](#) Data, float Min, float Max)  
*Scales the given data to the range [Min, Max].*
- **SPECTRALRADAR\_API** void **lockData** ([DataHandle](#) Data)

- Locks the given data.*

  - **SPECTRALRADAR\_API** void **unlockData** (**DataHandle** Data)

*Unlocks the given data.*
- **SPECTRALRADAR\_API** void **getDataSlicePos** (**DataHandle** Data, **DataHandle** Slice, **Direction** SliceNormalDirection, double Pos)

*Returns a slice of data in the specified direction at the specified position.*
- **SPECTRALRADAR\_API** void **getDataSliceIndex** (**DataHandle** Data, **DataHandle** Slice, **Direction** SliceNormalDirection, int Index)

*Returns a slice of data in the specified direction at the specified index.*
- **SPECTRALRADAR\_API** void **getDataSliceAnalyzed** (**DataHandle** Data, **DataHandle** Slice, **Direction** SliceNormalDirection, **DataAnalysis** Selection)

*Returns a slice of data that has been computed of all slice using the specified analyzation method.*
- **SPECTRALRADAR\_API** void **appendData** (**DataHandle** Data, **DataHandle** NewData, **Direction** direction)

*Appends the new data to the old data in the specified direction.*
- **SPECTRALRADAR\_API** void **cropData** (**DataHandle** Data, **Direction** direction, int Index)

*Crops the data at the specific direction at the given index. The result will contain the data with range [0, index] at the cropping direction.*
- **SPECTRALRADAR\_API** void **cropDataEx** (**DataHandle** Data, **Direction** direction, int IndexMax, int IndexMin)

*Crops the data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping direction.*
- **SPECTRALRADAR\_API** void **separateData** (**DataHandle** Data1, **DataHandle** Data2, int SeparationIndex, **Direction** Dir)

*Separates the data at the given index at specific separation direction. The first part of the separated data will be in Data1, the second separated in Data2.*
- **SPECTRALRADAR\_API** void **flipData** (**DataHandle** Data, **Direction** FlippingDirection)

*Flips the data around the specific direction.*
- **SPECTRALRADAR\_API** void **appendComplexData** (**ComplexDataHandle** Data, **ComplexDataHandle** DataToAppend, **Direction** direction)

*Appends the new data to the old data in the specified direction.*
- **SPECTRALRADAR\_API** void **getComplexDataSlicePos** (**ComplexDataHandle** Data, **ComplexDataHandle** Slice, **Direction** SliceNormalDirection, double Pos)

*Returns a slice of data in the specified direction at the specified position.*
- **SPECTRALRADAR\_API** void **getComplexDataSliceIndex** (**ComplexDataHandle** Data, **ComplexDataHandle** Slice, **Direction** SliceNormalDirection, int Index)

*Returns a slice of data in the specified direction at the specified index.*
- **SPECTRALRADAR\_API** void **cropComplexData** (**ComplexDataHandle** Data, **Direction** CroppingDirection, int IndexMax, int IndexMin)

*Crops the complex data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping direction.*
- **SPECTRALRADAR\_API** void **appendColoredData** (**ColoredDataHandle** Data, **ColoredDataHandle** DataToAppend, **Direction** AppendingDirection)

*Appends the new colored data to the old colored data in the specified direction.*
- **SPECTRALRADAR\_API** void **cropColoredData** (**ColoredDataHandle** Data, **Direction** CroppingDirection, int IndexMax, int IndexMin)

*Crops the colored data at the specific direction at the given indeces. The result will contain the data with range [IndexMin, IndexMax] at the cropping direction.*
- **SPECTRALRADAR\_API** void **getColoredDataSlicePos** (**ColoredDataHandle** Data, **ColoredDataHandle** Slice, **Direction** SliceNormalDirection, double Pos)

*Get a slice of the colored data with specific slicing direction at given index position.*
- **SPECTRALRADAR\_API** void **getColoredDataSliceIndex** (**ColoredDataHandle** Data, **ColoredDataHandle** Slice, **Direction** SliceNormalDirection, int Index)

*Get a slice of the colored data with specific slicing direction at given index.*
- **SPECTRALRADAR\_API** **ImageFieldHandle** **createImageField** (void)

- Creates an object holding image field data.*
- **SPECTRALRADAR\_API** void **clearImageField** (**ImageFieldHandle** ImageField)
- Frees an object holding image field data.*
- **SPECTRALRADAR\_API** void **saveImageField** (**ImageFieldHandle** ImageField, const char \*Path)
- Saves data containing image field data.*
- **SPECTRALRADAR\_API** void **loadImageField** (**ImageFieldHandle** ImageField, const char \*Path)
- Loads data containing image field data.*
- **SPECTRALRADAR\_API** void **determinelImageField** (**ImageFieldHandle** ImageField, **DataHandle** Surface)
- Determines the image field correction of the surface.*
- **SPECTRALRADAR\_API** void **determinelImageFieldForProbe** (**ProbeHandle** Probe, **DataHandle** Surface)
- Determines the image field correction of the surface for the specified probe handle.*
- **SPECTRALRADAR\_API** void **determinelImageFieldForProbeWithMap** (**ProbeHandle** Probe, **DataHandle** Surface, **DataHandle** Map)
- Determines the image field correction of the surface for the specified probe handle using the given map. Values != 0 in the map specifies to use the data in the surface handle otherwise thex will be interpolated.*
- **SPECTRALRADAR\_API** void **correctImageField** (**ImageFieldHandle** ImageField, **ScanPatternHandle** Pattern, **DataHandle** Data)
- Applies the image field correction to the B-Scan or volume data .*
- **SPECTRALRADAR\_API** void **correctImageFieldFromProbe** (**ProbeHandle** Probe, **ScanPatternHandle** Pattern, **DataHandle** Data)
- Applies the image field correction saved in the probe handle to the B-Scan or volume data .*
- **SPECTRALRADAR\_API** void **correctSurface** (**ImageFieldHandle** ImageField, **DataHandle** Surface)
- Applies the image field correction to the given Surface.*
- **SPECTRALRADAR\_API**  
**VisualCalibrationHandle** **createVisualCalibration** (**OCTDeviceHandle** Device, double TargetCornerLength\_  
mm, **BOOL** CheckAngle, **BOOL** SaveData)
- **SPECTRALRADAR\_API** void **clearVisualCalibration** (**VisualCalibrationHandle** Handle)
- **SPECTRALRADAR\_API** **BOOL** **visualCalibrate\_1st\_CameraScaling** (**VisualCalibrationHandle** Handle,  
**ProbeHandle** Probe, **ColoredDataHandle** Image)
- **SPECTRALRADAR\_API** **BOOL** **visualCalibrate\_2nd\_Galvo** (**VisualCalibrationHandle** Handle, **Probe-**  
**Handle** Probe, **ColoredDataHandle** Image)
- **SPECTRALRADAR\_API** **BOOL** **visualCalibrate\_previewImage** (**VisualCalibrationHandle** Handle, **Colored-**  
**DataHandle** Image)
- **SPECTRALRADAR\_API** void **visualCalibration\_getHoles** (**VisualCalibrationHandle** Handle, int \*x0, int \*y0,  
int \*x1, int \*y1, int \*x2, int \*y2)
- **SPECTRALRADAR\_API** const char \* **visualCalibrate\_Status** (**VisualCalibrationHandle** Handle)
- **SPECTRALRADAR\_API** **BOOL** **visualCalibrate\_CameraCenter** (**VisualCalibrationHandle** Handle, **OCT-**  
**DeviceHandle** Device, **ColoredDataHandle** Image)
- **SPECTRALRADAR\_API** **BOOL** **visualCalibrate\_getVideoCameraCenterImage** (**VisualCalibrationHandle**  
Handle, **OCTDeviceHandle** Device, **ColoredDataHandle** Image)
- **SPECTRALRADAR\_API**  
**DopplerProcessingHandle** **createDopplerProcessing** (void)
- **SPECTRALRADAR\_API** void **createDopplerProcessingForProcessing** (**DopplerProcessingHandle** \*Doppler,  
**ProcessingHandle** Proc)
- **SPECTRALRADAR\_API** void **closeDopplerProcessing** (**DopplerProcessingHandle** Doppler)
- Closes the Doppler processing routines and frees the memory that has been allocated for these to work properly.*
- **SPECTRALRADAR\_API** void **setDopplerPropertyInt** (**DopplerProcessingHandle** Doppler, **DopplerPropertyInt**  
Property, int Value)
- Sets Doppler processing properties.*
- **SPECTRALRADAR\_API** void **setDopplerPropertyFloat** (**DopplerProcessingHandle** Doppler, **Doppler-**  
**PropertyFloat** Property, float Value)
- Sets Doppler processing properties.*
- **SPECTRALRADAR\_API** void **setDopplerFlag** (**DopplerProcessingHandle** Doppler, **DopplerFlag** Flag, **BOOL**  
OnOff)

- Sets the Doppler processing flags.*

  - **SPECTRALRADAR\_API** void **setDopplerAmplitudeOutput** (**DopplerProcessingHandle** Doppler, **DataHandle** AmpOut)

*Sets the location of the resulting doppler amplitude output.*
- **SPECTRALRADAR\_API** void **setDopplerPhaseOutput** (**DopplerProcessingHandle** Doppler, **DataHandle** PhasesOut)

*Sets the location of the resulting doppler phase output.*
- **SPECTRALRADAR\_API** void **getDopplerOutputSize** (**DopplerProcessingHandle** Doppler, int Size1In, int Size2In, int \*Size1Out, int \*Size2Out)

*Returns the final size of the Doppler output if executeDopplerProcessing is executed using data of the specified input size.*
- **SPECTRALRADAR\_API** void **executeDopplerProcessing** (**DopplerProcessingHandle** Doppler, **ComplexDataHandle** Input)

*Executes the Doppler processing of the input data and returns phases and amplitudes.*
- **SPECTRALRADAR\_API** void **calcContrast** (**DataHandle** ApodizedSpectrum, **DataHandle** Contrast)

*Computes the contrast for the specified (apodized); spectrum.*
- **SPECTRALRADAR\_API** **SettingsHandle** **loadSettingsFile** (const char \*Path)

*Loads a settings file (usually \*.ini); and prepares its properties to be read.*
- **SPECTRALRADAR\_API** int **getSettingsEntryInt** (**SettingsHandle** SettingsFile, const char \*Node, int DefaultValue)

*Gets an integer number from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** double **getSettingsEntryDouble** (**SettingsHandle** SettingsFile, const char \*Node, double DefaultValue)

*Gets an floating point number from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** void **getSettingsEntryFloatArray** (**SettingsHandle** SettingsFile, const char \*Node, const double \*DefaultValues, double \*Values, int \*Size)
- **SPECTRALRADAR\_API** void **getSettingsEntryString** (**SettingsHandle** SettingsFile, const char \*Node, const char \*Default, char \*Data, int MaxDataSize)

*Gets a string from the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** void **setSettingsEntryInt** (**SettingsHandle** SettingsFile, const char \*Node, int Value)

*Sets an integer entry in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** void **setSettingsEntryFloat** (**SettingsHandle** SettingsFile, const char \*Node, double Value)

*Sets a floating point entry in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** void **saveSettings** (**SettingsHandle** SettingsFile)

*Saves the changes to the specified Settings file.*
- **SPECTRALRADAR\_API** void **closeSettingsFile** (**SettingsHandle** SettingsFile)

*Closes the specified ini file and stores the set entries (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** void **setSettingsEntryString** (**SettingsHandle** SettingsFile, const char \*Node, const char \*Value)

*Sets a string in the specified ini file (see [SettingsHandle](#) and [loadSettingsFile](#));.*
- **SPECTRALRADAR\_API** **Coloring32BitHandle** **createColoring32Bit** (**ColorScheme** Color, **ColoringByteOrder** ByteOrder)

*Creates processing that can be used to color given floating point B-scans to 32 bit colored images.*
- **SPECTRALRADAR\_API** **Coloring32BitHandle** **createCustomColoring32Bit** (int LUTSize, unsigned long LUT[])

*Create custom coloring using the specified color look-up-table.*
- **SPECTRALRADAR\_API** void **setColoringBoundaries** (**Coloring32BitHandle** Colorng, float Min\_dB, float Max\_dB)

*Sets the boundaries in dB which are used by the coloring algorithm to map colors to floating point values in dB.*
- **SPECTRALRADAR\_API** void **setColoringEnhancement** (**Coloring32BitHandle** Colorng, **ColorEnhancement** Enhancement)

*Selects a function for non-linear coloring to enhance (subjective) image impression.*

- **SPECTRALRADAR\_API** void **colorizeData** (**Coloring32BitHandle** Coloring, **DataHandle** Data, **ColoredDataHandle** ColoredData, **BOOL** Transpose)  
*Colors a given data object (**DataHandle**) into a given colored object (**ColoredDataHandle**).*
- **SPECTRALRADAR\_API** void **colorizeDopplerData** (**Coloring32BitHandle** AmpColoring, **Coloring32BitHandle** PhaseColoring, **DataHandle** AmpData, **DataHandle** PhaseData, **ColoredDataHandle** Output, double MinSignal\_dB, **BOOL** Transpose)  
*Colors a two given data object (**DataHandle**) using overlay and intensity to represent phase and amplitude data. Used for Doppler imaging.*
- **SPECTRALRADAR\_API** void **clearColoring32Bit** (**Coloring32BitHandle** Coloring)  
*Clears the coloring previously created by **createColoring32Bit()**.*
- **SPECTRALRADAR\_API** void **getMaxCameraImageSize** (**OCTDeviceHandle** Dev, int \*SizeX, int \*SizeY)  
*Returns the maximum possible camera image size for the current device.*
- **SPECTRALRADAR\_API** void **getCameraImage** (**OCTDeviceHandle** Dev, int SizeX, int SizeY, **ColoredDataHandle** Image)  
*Gets a camera image.*
- **SPECTRALRADAR\_API** void **getFlippedCameraImage** (**OCTDeviceHandle** Dev, int SizeX, int SizeY, **ColoredDataHandle** Image, **BOOL** InvertX, **BOOL** InvertY)
- **SPECTRALRADAR\_API** void **getMirroredCameraImage** (**OCTDeviceHandle** Dev, int SizeX, int SizeY, **ColoredDataHandle** Image)  
*Gets a camera image.*
- **SPECTRALRADAR\_API** void **setCameraPropertyFloat** (**OCTDeviceHandle** Dev, **CameraPropertyFloat** Section, double Value)  
*Sets saturation, brightness and contrast for the camera images if this option is available for the current device.*
- **SPECTRALRADAR\_API** void **setCameraShowScanPattern** (**OCTDeviceHandle** Dev, **BOOL** Value)  
*Enables to turn on/off the scan pattern overlay in the video camera image.*
- **SPECTRALRADAR\_API** void **visualizeScanPattern** (**OCTDeviceHandle** Dev, **ProbeHandle** Probe, **ScanPatternHandle** Pattern, **BOOL** showRawPattern)  
*Visualizes the scan pattern in top of the camera image.*
- **SPECTRALRADAR\_API** double **getReferenceIntensity** (**ProcessingHandle** Proc)  
*Returns an absolute value that indicates the reference intensity that was present when the currently used apodization was determined.*
- **SPECTRALRADAR\_API** double **getRelativeReferenceIntensity** (**OCTDeviceHandle** Dev, **ProcessingHandle** Proc)  
*Returns a value larger than 0.0 and smaller than 1.0 that indicates the reference intensity that was present when the currently used apodization was determined.*
- **SPECTRALRADAR\_API** void **getConfigPath** (char \*Path, int StrSize)  
*Returns the path that hold the config files.*
- **SPECTRALRADAR\_API** void **getPluginPath** (char \*Path, int StrSize)  
*Returns the path that hold the plugins.*
- **SPECTRALRADAR\_API** void **getInstallationPath** (char \*Path, int StrSize)  
*Returns the installation path.*
- **SPECTRALRADAR\_API** **BufferHandle** **createMemoryBuffer** (void)  
*Creates a buffer holding data and colored data.*
- **SPECTRALRADAR\_API** void **appendToBuffer** (**BufferHandle**, **DataHandle**, **ColoredDataHandle**)  
*Appends specified data and colored data to the requested buffer.*
- **SPECTRALRADAR\_API** int **getBufferSize** (**BufferHandle**)  
*Returns the currently available data sets in the buffer.*
- **SPECTRALRADAR\_API** **DataHandle** **getBufferData** (**BufferHandle**, int Index)  
*Returns the data in the buffer.*
- **SPECTRALRADAR\_API** **ColoredDataHandle** **getColoredBufferData** (**BufferHandle**, int Index)  
*Returns the colored data in the buffer.*
- **SPECTRALRADAR\_API** void **clearBuffer** (**BufferHandle**)

- Clears the buffer and frees all data and colored data objects in it.*
- **SPECTRALRADAR\_API** int **getNumberOfOutputValues** (OCTDeviceHandle Dev)
  - Returns the number of output values.*
  - **SPECTRALRADAR\_API** void **getOutputValueName** (OCTDeviceHandle Dev, int Index, char \*Name, int NameStringSize, char \*Unit, int UnitStringSize)
  - Returns names and units of the requested output values.*
  - **SPECTRALRADAR\_API** **BOOL** **doesOutputExist** (OCTDeviceHandle Dev, const char \*Name)
  - **SPECTRALRADAR\_API** void **setOutputValueByName** (OCTDeviceHandle Dev, const char \*Name, double value)
  - Sets the specified output value.*
  - **SPECTRALRADAR\_API** void **setOutputValueByIndex** (OCTDeviceHandle Dev, int Index, double Value)
  - Sets the specified output value.*
  - **SPECTRALRADAR\_API** void **getOutputValueRangeByName** (OCTDeviceHandle Dev, const char \*Name, double \*Min, double \*Max)
  - Gives the range of the specified output value.*
  - **SPECTRALRADAR\_API** void **getOutputValueRangeByIndex** (OCTDeviceHandle Dev, int Index, double \*Min, double \*Max)
  - Gives the range of the specified output value.*
  - **SPECTRALRADAR\_API** void **computeLinearKRawData** (ComplexDataHandle ComplexDataAfterFFT, DataHandle LinearKData)
  - Computes the linear k raw data of the complex data after FFT by an inverse Fourier transform.*
  - **SPECTRALRADAR\_API** void **linearizeSpectralData** (DataHandle SpectralIn, DataHandle SpectraOut, DataHandle Chirp)
  - Linearizes the spectral data using the given chirp vector.*
  - **SPECTRALRADAR\_API** void **TestFileEngine** ()
  - **SPECTRALRADAR\_API** OCTFileHandle **createOCTFile** (OCTFileFormat format)
  - **SPECTRALRADAR\_API** void **clearOCTFile** (OCTFileHandle handle)
  - **SPECTRALRADAR\_API** int **getFileDataFileCount** (OCTFileHandle handle)
  - **SPECTRALRADAR\_API** void **loadFile** (OCTFileHandle handle, const char \*fileName)
  - **SPECTRALRADAR\_API** void **saveFile** (OCTFileHandle handle, const char \*fileName)
  - **SPECTRALRADAR\_API** void **copyFileMetadata** (OCTFileHandle src, OCTFileHandle dest)
  - **SPECTRALRADAR\_API** double **getFileMetadataFloat** (OCTFileHandle handle, FileMetadataFloatField floatfield)
  - **SPECTRALRADAR\_API** void **setFileMetadataFloat** (OCTFileHandle handle, FileMetadataFloatField floatfield, double value)
  - **SPECTRALRADAR\_API** int **getFileMetadataInt** (OCTFileHandle handle, FileMetadataIntField intfield)
  - **SPECTRALRADAR\_API** void **setFileMetadataInt** (OCTFileHandle handle, FileMetadataIntField intfield, int value)
  - **SPECTRALRADAR\_API** size\_t **getFileMetadataString** (OCTFileHandle handle, FileMetadataStringField stringfield, char \*content, int length)
  - **SPECTRALRADAR\_API** void **setFileMetadataString** (OCTFileHandle handle, FileMetadataStringField stringfield, const char \*content)
  - **SPECTRALRADAR\_API** **BOOL** **getFileMetadataBool** (OCTFileHandle handle, FileMetadataBoolField boolfield)
  - **SPECTRALRADAR\_API** void **setFileMetadataBool** (OCTFileHandle handle, FileMetadataBoolField boolfield, **BOOL** value)
  - **SPECTRALRADAR\_API** DataHandle **getFileRealData** (OCTFileHandle handle, size\_t index)
  - **SPECTRALRADAR\_API** ColoredDataHandle **getFileColoredData** (OCTFileHandle handle, size\_t index)
  - **SPECTRALRADAR\_API** ComplexDataHandle **getFileComplexData** (OCTFileHandle handle, size\_t index)
  - **SPECTRALRADAR\_API** RawDataHandle **getFileRawData** (OCTFileHandle handle, size\_t index)
  - **SPECTRALRADAR\_API** ColoredDataHandle **getFileBinaryData** (OCTFileHandle handle, size\_t index)
  - **SPECTRALRADAR\_API** void **getFile** (OCTFileHandle handle, size\_t index, const char \*filenameOnDisk)
  - **SPECTRALRADAR\_API** void **addFileRealData** (OCTFileHandle handle, DataHandle data, const char \*filename)



- [SPECTRALRADAR\\_API](#) void **addFileColoredData** ([OCTFileHandle](#) handle, [ColoredDataHandle](#) data, const char \*filename)
- [SPECTRALRADAR\\_API](#) void **addFileComplexData** ([OCTFileHandle](#) handle, [ComplexDataHandle](#) data, const char \*filename)
- [SPECTRALRADAR\\_API](#) void **addFileRawData** ([OCTFileHandle](#) handle, [RawDataHandle](#) data, const char \*filename)
- [SPECTRALRADAR\\_API](#) void **addFileBinaryData** ([OCTFileHandle](#) handle, [ColoredDataHandle](#) data, const char \*filename)
- [SPECTRALRADAR\\_API](#) void **addFileBinary** ([OCTFileHandle](#) handle, const char \*filenameOnDisk, const char \*filename)
- [SPECTRALRADAR\\_API](#) void **addFileText** ([OCTFileHandle](#) handle, const char \*filenameOnDisk, const char \*filename)
- [SPECTRALRADAR\\_API](#) [DataKind](#) **getFileDataFileType** ([OCTFileHandle](#) handle, int index)
- [SPECTRALRADAR\\_API](#) void **getFileDataFileName** ([OCTFileHandle](#) handle, int index, char \*filename, int length)
- [SPECTRALRADAR\\_API](#) int **getFileDataSizeX** ([OCTFileHandle](#) handle, [size\\_t](#) index)
- [SPECTRALRADAR\\_API](#) int **getFileDataSizeY** ([OCTFileHandle](#) handle, [size\\_t](#) index)
- [SPECTRALRADAR\\_API](#) int **getFileDataSizeZ** ([OCTFileHandle](#) handle, [size\\_t](#) index)
- [SPECTRALRADAR\\_API](#) float **getFileDataRangeX** ([OCTFileHandle](#) handle, [size\\_t](#) index)
- [SPECTRALRADAR\\_API](#) float **getFileDataRangeY** ([OCTFileHandle](#) handle, [size\\_t](#) index)
- [SPECTRALRADAR\\_API](#) float **getFileDataRangeZ** ([OCTFileHandle](#) handle, [size\\_t](#) index)
- [SPECTRALRADAR\\_API](#) void **setMarkerListFromRealData** ([OCTFileHandle](#) handle, [DataHandle](#) data)
- [SPECTRALRADAR\\_API](#) void **setMarkerListInRealData** ([OCTFileHandle](#) handle, [DataHandle](#) data)
- [SPECTRALRADAR\\_API](#)  
[SpeckleVarianceHandle](#) **initSpeckleVariance** (void)
- [SPECTRALRADAR\\_API](#) void **closeSpeckleVariance** ([SpeckleVarianceHandle](#) [SpeckleVar](#))
- [SPECTRALRADAR\\_API](#) void **setAveraging** ([SpeckleVarianceHandle](#) [SpeckleVar](#), int Av1, int Av2, int Av3)
- [SPECTRALRADAR\\_API](#) void **setSpeckleVarianceType** ([SpeckleVarianceHandle](#) [SpeckleVar](#), [SpeckleVarianceType](#) Type)
- [SPECTRALRADAR\\_API](#) void **computeSpeckleVariance** ([SpeckleVarianceHandle](#) [SpeckleVar](#), [ComplexDataHandle](#) [CompDataIn](#), [DataHandle](#) [DataOutMean](#), [DataHandle](#) [DataOutVar](#))
- [SPECTRALRADAR\\_API](#) void **setSpeckleVarianceThreshold** ([SpeckleVarianceHandle](#) [SpeckleVar](#), double Threshold)
- [SPECTRALRADAR\\_API](#) void **setTriggerMode** ([OCTDeviceHandle](#) Dev, [Device\\_TriggerType](#) TriggerMode)
- [SPECTRALRADAR\\_API](#)  
[Device\\_TriggerType](#) **getTriggerMode** ([OCTDeviceHandle](#) Dev)  
*Returns the trigger mode used for acquisition.*
- [SPECTRALRADAR\\_API](#) bool **isTriggerModeAvailable** ([OCTDeviceHandle](#) Dev, [Device\\_TriggerType](#) TriggerMode)  
*Returns whether the specified trigger mode is possible or not for the used device.*
- [SPECTRALRADAR\\_API](#) void **setTriggerTimeoutSec** ([OCTDeviceHandle](#) Dev, int Timeout)
- [SPECTRALRADAR\\_API](#) int **getTriggerTimeoutSec** ([OCTDeviceHandle](#) Dev)  
*Returns the timeout of the camera (not used in trigger mode [Trigger\\_FreeRunning](#)).*
- [SPECTRALRADAR\\_API](#) int **getScanPatternPropertyInt** ([ScanPatternHandle](#) [ScanPattern](#), [ScanPatternPropertyInt](#) Property)
- [SPECTRALRADAR\\_API](#) double **expectedAcquisitionTimeSec** ([ScanPatternHandle](#) [ScanPattern](#), [OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) void **interpolateScanPattern** ([DataHandle](#) [DataIn](#), [DataHandle](#) [DataOut](#), [InterpolationMethod](#) method, [BoundaryCondition](#) condition)
- [SPECTRALRADAR\\_API](#) double **meanFreeformScanPatternPathLength** ([DataHandle](#) [ScanPattern](#))
- [SPECTRALRADAR\\_API](#) [BOOL](#) **checkAvailableMemoryForScanPattern** ([OCTDeviceHandle](#) Dev, [ScanPatternHandle](#) [Pattern](#), [ptrdiff\\_t](#) [AdditionalMemory](#))  
*Checks whether sufficient memory is available for acquiring the specified scan pattern.*
- [SPECTRALRADAR\\_API](#) void **startMeasurePulseResponse** ([OCTDeviceHandle](#) Dev, double time\_ms, double voltage, int \*samplesPerChannel)

- [SPECTRALRADAR\\_API](#) void **getPulseResponseInput** ([OCTDeviceHandle](#) Dev, double \*Raw)
- [SPECTRALRADAR\\_API](#) void **getPulseResponse** ([OCTDeviceHandle](#) Dev, double \*RawX, double \*RawY)
- [SPECTRALRADAR\\_API](#) void **stopMeasurePulseResponse** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) void **startScanBenchmark** ([OCTDeviceHandle](#) Dev, double \*RawX, double \*RawY, int Size, double RateHz)
- [SPECTRALRADAR\\_API](#) void **stopScanBenchmark** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) void **getScanFeedback** ([OCTDeviceHandle](#) Dev, double \*RawX, double \*RawY)
- [SPECTRALRADAR\\_API](#) int **getScanFeedbackSize** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) double **QuantumEfficiency** ([OCTDeviceHandle](#) Dev, double CenterWavelength\_nm, double PowerIntoSpectrometer\_W, [DataHandle](#) Spectrum\_e)
- [SPECTRALRADAR\\_API](#) FullRangeHandle **initFullRange** ()
- [SPECTRALRADAR\\_API](#) void **executeFullRange** (FullRangeHandle FullRange, [DataHandle](#) ApodizedDataIn, [ComplexDataHandle](#) ApodizedDataOut)
- [SPECTRALRADAR\\_API](#) void **closeFullRange** (FullRangeHandle FullRange)
- [SPECTRALRADAR\\_API](#) void **setFullRangeSensitivity** (FullRangeHandle FullRangeObject, float CutOff, float Smoothness)
- [SPECTRALRADAR\\_API](#) void **executeComplexProcessing** ([ProcessingHandle](#) Proc, [ComplexDataHandle](#) ApodizedSpectralData)
- [SPECTRALRADAR\\_API](#) void **determineThickness** ([DataHandle](#) Data, float \*front, float \*back)
- [SPECTRALRADAR\\_API](#) void **flattenImage** ([DataHandle](#) ImageData)
- [SPECTRALRADAR\\_API](#) void **determineSurface** ([DataHandle](#) Volume, [DataHandle](#) Surface)
- [SPECTRALRADAR\\_API](#) unsigned long long **getFreeMemory** ()
- [SPECTRALRADAR\\_API](#) void **absComplexData** ([ComplexDataHandle](#) ComplexData, [DataHandle](#) Abs)
- [SPECTRALRADAR\\_API](#) void **logAbsComplexData** ([ComplexDataHandle](#) ComplexData, [DataHandle](#) dB)
- [SPECTRALRADAR\\_API](#) void **argComplexData** ([ComplexDataHandle](#) ComplexData, [DataHandle](#) Arg)
- [SPECTRALRADAR\\_API](#) void **realComplexData** ([ComplexDataHandle](#) ComplexData, [DataHandle](#) Real)
- [SPECTRALRADAR\\_API](#) void **imagComplexData** ([ComplexDataHandle](#) ComplexData, [DataHandle](#) Imag)
- [SPECTRALRADAR\\_API](#) void **equalizeColoredDataHistogram** ([ColoredDataHandle](#) ColoredData)
- [SPECTRALRADAR\\_API](#) void **equalizeDataHistogram** ([DataHandle](#) Data, double Min, double Max)
- [SPECTRALRADAR\\_API](#) void **medianFilter** ([DataHandle](#) Data, int Rank)  
*Computes a median filter on the specified 2D data.*
- [SPECTRALRADAR\\_API](#) void **pepperFilter** ([DataHandle](#) Data, [PepperFilterType](#) Type, float Threshold)  
*Removes pepper-noise (very low values, i. e. dark spots in the data). This enhances the visual (colored) representation of the data.*
- [SPECTRALRADAR\\_API](#) void **polynomialFilter** ([DataHandle](#) Data, int SizeX, int SizeY)
- [SPECTRALRADAR\\_API](#) void **gaussianFilter** ([DataHandle](#) Data, [GaussianFilterType](#) Type)
- [SPECTRALRADAR\\_API](#) void **prewittFilter** ([DataHandle](#) Data, [PrewittFilterType](#) Type)
- [SPECTRALRADAR\\_API](#) void **sobelFilter** ([DataHandle](#) Data, [SobelFilterType](#) Type)
- [SPECTRALRADAR\\_API](#) void **laplacianFilter** ([DataHandle](#) Data, [LaplacianFilterType](#) Type)
- [SPECTRALRADAR\\_API](#) void **applyNonlinearSobelFilter2D** ([DataHandle](#) Data)
- [SPECTRALRADAR\\_API](#) void **applyNonlinearPrewittFilter2D** ([DataHandle](#) Data)
- [SPECTRALRADAR\\_API](#) void **applyConvolutionFilter2D** ([DataHandle](#) Data, int \*Filter)
- [SPECTRALRADAR\\_API](#) void **applyMedianFilter2D\_1x** ([DataHandle](#) Data, int Rank)
- [SPECTRALRADAR\\_API](#) void **applyFilter1D** ([DataHandle](#) Data, int \*Size, float \*FilterKernel, [Direction](#) FilterDirection=[Direction\\_1](#), bool Normalization=false)
- [SPECTRALRADAR\\_API](#) void **applyFilter2D** ([DataHandle](#) Data, int \*Size, float \*FilterKernel, [Direction](#) FilterNormalDirection=[Direction\\_3](#), bool Normalization=false)
- [SPECTRALRADAR\\_API](#) void **applyFilter3D** ([DataHandle](#) Data, int \*Size, float \*FilterKernel, bool Normalization=false)
- [SPECTRALRADAR\\_API](#) void **applyFilter** ([ComplexDataHandle](#) ComplexData, [FilterType](#) Type, double FilterParameter)
- [SPECTRALRADAR\\_API](#) void **normalizeFilter** (float \*FilterKernel, int FilterSize)
- [SPECTRALRADAR\\_API](#) void **smoothCurve1D** (int Size, float \*Curve, int DegreePolynom, int NumberPolynoms=1)



- [SPECTRALRADAR\\_API](#) float \* **polynomialFitAndEval1D** (int Size, float \*OrigPosX, float \*OrigY, int DegreePolynom, int EvalSize, float \*EvalPosX)
- [SPECTRALRADAR\\_API](#) float **calcParableMaximum** (float x0, float y0, float yp1, float ym1, float \*peak\_height=nullptr)
- [SPECTRALRADAR\\_API](#) void **thresholdData** ([DataHandle](#) Data, double Threshold)
- [SPECTRALRADAR\\_API](#) void **levelData** ([DataHandle](#) Data)
- Levels the specified data and removes tilt.*
- [SPECTRALRADAR\\_API](#) void **importRealBinaryData** ([DataHandle](#) RealData, int Size1, int Size2, int Size3, const char \*Path)
- [SPECTRALRADAR\\_API](#) void **filterDC** ([DataHandle](#) Data)
- [SPECTRALRADAR\\_API](#) void **crossCorrelatedProjection** ([DataHandle](#) DataIn, [DataHandle](#) Res)
- [SPECTRALRADAR\\_API](#) void **thresholdDataPtByPt** ([DataHandle](#) Phase, [DataHandle](#) Intensity, float threshold, float targetValue)
- [SPECTRALRADAR\\_API](#) void **analyzeScatteringSample** ([DataHandle](#) Data)
- [SPECTRALRADAR\\_API](#) void **getCurrentIntensityStatistics** ([OCTDeviceHandle](#) Dev, [ProcessingHandle](#) Proc, float \*relToRefIntensity, float \*relToProjAbsIntensity)
- [SPECTRALRADAR\\_API](#) void **getCurrentApodizationEdgeChannels** ([ProcessingHandle](#) Proc, int \*LeftPix, int \*RightPix)
- [SPECTRALRADAR\\_API](#) int **getNumberOfProbeConfigs** ()
- [SPECTRALRADAR\\_API](#) void **getProbeConfigName** (int index, char \*ProbeName, int StringSize)
- [SPECTRALRADAR\\_API](#) [PolarizationProcessingHandle](#) **createPolarizationProcessing** (void)
- [SPECTRALRADAR\\_API](#) void **closePolarizationProcessing** ([PolarizationProcessingHandle](#) Polarization)
- [SPECTRALRADAR\\_API](#) void **setPolarizationPropertyInt** ([PolarizationProcessingHandle](#) Polarization, [PolarizationPropertyInt](#) Property, int Value)
- [SPECTRALRADAR\\_API](#) void **setPolarizationPropertyFloat** ([PolarizationProcessingHandle](#) Polarization, [PolarizationPropertyFloat](#) Property, float Value)
- [SPECTRALRADAR\\_API](#) void **setPolarizationFlag** ([PolarizationProcessingHandle](#) Polarization, [PolarizationFlag](#) Flag, [BOOL](#) OnOff)
- [SPECTRALRADAR\\_API](#) void **setPolarizationIntensityOutput** ([PolarizationProcessingHandle](#) Polarization, [DataHandle](#) AmpOut)
- [SPECTRALRADAR\\_API](#) void **setPolarizationRetardationOutput** ([PolarizationProcessingHandle](#) Polarization, [DataHandle](#) PhasesOut)
- [SPECTRALRADAR\\_API](#) void **executePolarizationProcessing** ([PolarizationProcessingHandle](#) Polarization, [ComplexDataHandle](#) SData, [ComplexDataHandle](#) PData)
- [SPECTRALRADAR\\_API](#) [BOOL](#) **initUSBProbeCtrl** ([OCTDeviceHandle](#) dev)
- [SPECTRALRADAR\\_API](#) [BOOL](#) **configureUSBProbeCtrlButton** ([OCTDeviceHandle](#) dev, [USBProbeButtonID](#) btn, [USBProbeCommand](#) cmd)
- [SPECTRALRADAR\\_API](#) [BOOL](#) **getLastUSBProbeMessage** ([OCTDeviceHandle](#) dev, char \*msg, size\_t size)
- [SPECTRALRADAR\\_API](#) [BOOL](#) **toggleUSBProbeLED** ([OCTDeviceHandle](#) Dev, int LED, [BOOL](#) OnOff)
- [SPECTRALRADAR\\_API](#) [BOOL](#) **refstageAvailable** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) [RefstageStatus](#) **refstageGetStatus** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) double **refstageGetLength\_mm** ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe)
- [SPECTRALRADAR\\_API](#) double **refstageGetPosition\_mm** ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe)
- [SPECTRALRADAR\\_API](#) void **refstageHome** ([OCTDeviceHandle](#) Dev, [BOOL](#) wait)
- [SPECTRALRADAR\\_API](#) void **refstageMoveLonger** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) void **refstageMoveShorter** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) void **refstageMoveAbsolute** ([OCTDeviceHandle](#) Dev, [ProbeHandle](#) Probe, double pos\_mm)
- [SPECTRALRADAR\\_API](#) void **refstageStop** ([OCTDeviceHandle](#) Dev)
- [SPECTRALRADAR\\_API](#) void **refstageSetSpeed** ([OCTDeviceHandle](#) Dev, [RefstageSpeed](#) speed)
- [SPECTRALRADAR\\_API](#) void **refstageSetStatusCallback** ([OCTDeviceHandle](#) Dev, [cbRefstageStatusChanged](#) Callback)
- [SPECTRALRADAR\\_API](#) void **refstageSetPosChangeCallback** ([OCTDeviceHandle](#) Dev, [cbRefstagePositionChanged](#) Callback)

### 7.1.1 Detailed Description

Header containing all functions of the Spectral Radar SDK. This SDK can be used for Callisto, Ganymede, Hyperion and Telesto devices.

### 7.1.2 Macro Definition Documentation

#### 7.1.2.1 `#define FALSE 0`

FALSE for use with data type `BOOL`.

#### 7.1.2.2 `#define SPECTRALRADAR_API __declspec(dllexport)`

Export/Import of define of DLL members.

#### 7.1.2.3 `#define TRUE 1`

TRUE for use with data type `BOOL`.

### 7.1.3 Typedef Documentation

#### 7.1.3.1 `BOOL`

A standard boolean data type used in the API.

#### 7.1.3.2 `MarkerListHandle`

Handle to the marker list class.

#### 7.1.3.3 `VisualCalibrationHandle`

Handle to the visual galvo calibration class.

### 7.1.4 Enumeration Type Documentation

#### 7.1.4.1 `enum Device_TriggerType`

Enumerator

***Trigger\_FreeRunning*** Standard mode.

***Trigger\_TrigBoard\_ExternalStart*** Used to trigger the start of an acquisition. Additional hardware is needed.

***Trigger\_External\_AScan*** Mode to trigger the acquisition of each A-scan. An external trigger signal is needed. Please see the software manual for detailed information.

#### 7.1.4.2 `enum ScanPattern_AcquisitionOrder`

Enumerator

***ScanPattern\_AcqOrderFrameByFrame*** The scan pattern will be acquired slice by slice which means that the function `GetRawData()` needs to be called more than once to get the data for the whole scan pattern.

***ScanPattern\_AcqOrderAll*** The scan patten will be acquired in one piece.

#### 7.1.4.3 `enum ScanPatternPropertyInt`

Enumerator

***ScanPattern\_SizeTotal*** Total count of trigger pulses needed for acquisition of the scan pattern once. The acquisition will start again after finishing for continuous acquisition mode.

**ScanPattern\_Cycles** Count of cycles for the scan pattern.

**ScanPattern\_SizeCycle** Count of trigger pulses needed to acquire one cycle, e.g. one B-scan in a volume scan.

**ScanPattern\_SizePreparationCycle** Count of trigger pulses needed before the scanning of the sample starts. The OCT beam needs to be positioned and the apodization scans used for processing need to be acquired. The flyback time is the time used to reach the position of apodization and start of scan pattern.

**ScanPattern\_SizeImagingCycle** Count of trigger pulses to acquire the sample depending on averaging and size-x of the scan pattern.

### 7.1.5 Function Documentation

#### 7.1.5.1 **BOOL** checkAvailableMemoryForScanPattern ( **OCTDeviceHandle** *Dev*, **ScanPatternHandle** *Pattern*, **ptrdiff\_t** *AdditionalMemory* )

Checks whether sufficient memory is available for acquiring the specified scan pattern.

**AdditionalMemory** The parameter specifies additional memory that will be required during the measurement (from `startMeasurement()` to `stopMeasurement()`) unknown to the SDK and/or memory that will be freed/available prior to the call of `startMeasurement()`.

#### 7.1.5.2 **void** getConfigPath ( **char \*** *Path*, **int** *StrSize* )

Returns the path that hold the config files.

#### 7.1.5.3 **void** getDopplerOutputSize ( **DopplerProcessingHandle** *Doppler*, **int** *Size1In*, **int** *Size2In*, **int \*** *Size1Out*, **int \*** *Size2Out* )

Returns the final size of the Doppler output if `executeDopplerProcessing` is executed using data of the specified input size.

Doppler

#### 7.1.5.4 **void** getInstallationPath ( **char \*** *Path*, **int** *StrSize* )

Returns the installation path.

#### 7.1.5.5 **int** getNumberOfOutputValues ( **OCTDeviceHandle** *Dev* )

Returns the number of output values.

#### 7.1.5.6 **void** getOutputValueName ( **OCTDeviceHandle** *Dev*, **int** *Index*, **char \*** *Name*, **int** *NameStringSize*, **char \*** *Unit*, **int** *UnitStringSize* )

Returns names and units of the requested output values.

#### 7.1.5.7 **void** getOutputValueRangeByIndex ( **OCTDeviceHandle** *Dev*, **int** *Index*, **double \*** *Min*, **double \*** *Max* )

Gives the range of the specified output value.

#### 7.1.5.8 **void** getOutputValueRangeByName ( **OCTDeviceHandle** *Dev*, **const char \*** *Name*, **double \*** *Min*, **double \*** *Max* )

Gives the range of the specified output value.

#### 7.1.5.9 **void** getPluginPath ( **char \*** *Path*, **int** *StrSize* )

Returns the path that hold the plugins.

#### 7.1.5.10 double getReferenceIntensity ( ProcessingHandle Proc )

Returns an absolute value that indicates the reference intensity that was present when the currently used apodization was determined.

#### 7.1.5.11 double double getRelativeReferenceIntensity ( OCTDeviceHandle Dev, ProcessingHandle Proc )

Returns a value larger than 0.0 and smaller than 1.0 that indicates the reference intensity that was present when the currently used apodization was determined.

#### 7.1.5.12 Device\_TriggerType getTriggerMode ( OCTDeviceHandle Dev )

Returns the trigger mode used for acquisition.

#### 7.1.5.13 double getTriggerTimeoutSec ( OCTDeviceHandle Dev )

Returns the timeout of the camera (not used in trigger mode Trigger\_FreeRunning).

#### 7.1.5.14 unsigned long InterpretReferenceIntensity ( float intensity )

interprets the reference intensity and gives a color code that reflects its state.

Possible colors include:

- red = 0x00FF0000 (bad intensity);
- orange = 0x00FF7700 (okay intensity);
- green = 0x0000FF00 (good intensity);

##### Parameters

|                  |                                                                |
|------------------|----------------------------------------------------------------|
| <i>intensity</i> | the current reference intensity as a value between 0.0 and 1.0 |
|------------------|----------------------------------------------------------------|

##### Returns

the color code reflecting the state of the reference intensity

#### 7.1.5.15 bool isTriggerModeAvailable ( OCTDeviceHandle Dev, Device\_TriggerType TriggerMode )

Returns whether the specified trigger mode is possible or not for the used device.

#### 7.1.5.16 void setOutputValueByIndex ( OCTDeviceHandle Dev, int Index, double Value )

Sets the specified output value.

#### 7.1.5.17 void setOutputValueByName ( OCTDeviceHandle Dev, const char \* Name, double value )

Sets the specified output value.